

High Performance Routine Characterization of Proteins and Peptides by LC/MS

John C Gebler, Weibin Chen, and Asish Chakraborty

Waters Corporation, Biopharmaceutical Sciences, 34 Maple Street, Milford MA 01757, USA

ABSTRACT

All too often high performance combined with routine analysis methods result in a compromise in system performance. This work will illustrate how one can obtain high-performance in a routine manner without compromises. This LC/MS system and methods will be demonstrated with intact proteins and peptide maps from antibodies and recombinant proteins. The value to a work-flow of an optimized LC, columns, MS, and software for data analysis will also be presented.

INTRODUCTION

Comprehensive characterization of biopolymers as drug candidates is a requirement for safety and regulatory agencies. LC/MS of proteins and peptides is a powerful method for characterization. The ability to reproducibly generate and interpret high quality LC/MS data in a timely manner is challenging. This is often hindered with the lack of tools addressing a complete work flow. This poster outlines a complete set of work flow tools for time efficient high performance LC/MS characterization of proteins and peptides.



Figure 1. Waters biopharmaceutical LC/MS system: ACQUITY UPLC[®] LC and Xevo[™] QToF MS

METHODS

Samples:

- Humanized IgG (intact and reduced)
- RapiGest SF[™] aided Trypsin digested IgG (Peptide mapping)

Instrumentation:

- Liquid chromatography, ACQUITY UPLC[®]
- Mass spectrometer, Xevo[™] QTof


Columns:

- Intact IgG, MassPREP[™] Micro desalting column polymeric 20μm particle size 1000Å pore size, 2.1×5 mm
- Separation of IgG heavy/light chains, Protein Separation Technology, ACQUITY UPLC[®] C4 BEH300, 1.7μm, 2.1×50 mm
- Peptide mapping, Peptide Separation Technology, ACQUITY UPLC[®] C18 BEH300, 1.7μm, 2.1×150 mm


Informatics:

- BiopharmaLynx[™] ver 1.2 application manager


EXPERIMENTAL



INTACT LC/MS of IgG



PEPTIDE MAP



CONCLUSION

- High performance routine characterization of intact proteins and peptide maps
- No compromise LC/MS performance
- Featuring
 - Waters ACQUITY UPLC[®] LC System
 - ACQUITY UPLC Columns (1.7μm particles)
 - Xevo[™] QTof Mass Spectrometer
 - BiopharmaLynx[™] informatics
- Ideal for characterization of:
 - Intact Proteins
 - Peptide Mapping
 - Oligonucleotides
- Advance characterization of peptides with MS^E enabled with BiopharmaLynx
- Complete solution (work flow) tightly integrated

Protein	Peptide	Fragment Number	Modifiers	Control b/y Found	Control b/y List
humanized mab	DIQMTQSPSSLSASV	1:T001		28 1/b2;1/b3;1/b4;1/b5;1/b6;1/b7;1/b9;1/b10;1/b11;1/b12;1/b13;1/b14;1/b15	
humanized mab	DIQMTQSPSSLSASV	1:T001*		9 1/b2;1/b7;1/b10;1/y3;1/y4;1/y5;1/y6;1/y7;1/y8	
humanized mab	DIQMTQSPSSLSASV	1:T001-002*			
humanized mab	VITTCR	1:T002*	Carbamidomethyl C(1)	4 1/b2;1/y1;1/y2";1/y3"	
humanized mab	VITTCRASQDVNTAV	1:T002-003*	Carbamidomethyl C(1)		
humanized mab	ASQDVNNTAVAVYQQ	1:T003		22 1/b2;1/b3;1/b4;1/b5;1/b8;1/y1;1/y2;1/y3;1/y4;1/y5;1/y6;1/y7;1/y8;1/y9	
humanized mab	ASQDVNNTAVAVYQQ	1:T003*	Deamidation N(1)	14 1/b2;1/b4;1/b8";1/y2;1/y3;1/y4;1/y6;1/y7;1/y8;1/y13;1/y14;1/y15	
humanized mab	ASQDVNNTAVAVYQQ	1:T003-004			
humanized mab	API	1:T004			
humanized mab	APKLIIYASAFLYSGV	1:T004-005			
humanized mab	LLYSASFYSGVPSI	1:T005		14 1/b2;1/b3;1/b9;1/b13;1/y3;1/y5;1/y6;1/y7;1/y8;1/y10;1/y11;1/y12	
humanized mab	LLYSASFYSGVPSI	1:T005*	carbamidomethyl Y	0	
humanized mab	LLYSASFYSGVPSI	1:T005-006			
humanized mab	FSGSR	1:T006			
humanized mab	FSGSRSGTDFLTISS	1:T006-007*	Carbamidomethyl C(1)	4 1/b5;1/y1;1/y3;1/y4	
humanized mab	SGSR	1:T006/y4			
humanized mab	SGTDFITLTISSLQPE	1:T007*	Carbamidomethyl C	2 1/y1;1/y3	
humanized mab	SGTDFITLTISSLQPE	1:T007-008*	Carbamidomethyl C(1)	24 1/b4;1/b5;1/b6;1/b7;1/b8;1/b9;1/b10;1/b11;1/b12;1/b13;1/b28";1/b37";1/b38	
humanized mab	VEIK	1:T008			
humanized mab	VEIKR	1:T008-009			
humanized mab	R	1:T009			
humanized mab	RTVAAPS/FIFFFFSDI	1:T009-010			
humanized mab	TVAPS/FIFFFFSDI	1:T010		16 1/b2;1/b3;1/b4;1/b6;1/b8;1/b9;1/b10;1/y7;1/y8;1/y9;1/y10;1/y11;1/y12	
humanized mab	TVAPS/FIFFFFSDI	1:T010-011*	Carbamidomethyl C(1)		
humanized mab	TVAPS/FIFFFFSDI	1:T010/010			
humanized mab	PPSDEQLK	1:T010/y8		6 1/b2;1/b3;1/b4;1/b6;1/b9	
humanized mab	PPSDEQLK	1:T010/y8	1/y7		
humanized mab	SGTASIVCLNNFYP	1:T011*	Carbamidomethyl C	25 1/b2;1/b3;1/b4;1/b5;1/b6;1/b7;1/b8";1/b9";1/b10";1/b11;1/b12;1/b13;1/b14	
humanized mab	SGTASIVCLNNFYP	1:T011*	Deamidation N(1);C	10 1/b14";1/y2;1/y3;1/y4;1/y5;1/y6";1/y7;1/y8";1/y9";1/y10	
humanized mab	SGTASIVCLNNFYP	1:T011/b16*	Carbamidomethyl C	13 1/b2;1/b3;1/b4;1/b5;1/b6;1/b7;1/b8";1/b9";1/b10";1/b11;1/b12;1/b13	

Table 1
Results of BiopharmaLynx processed Xevo QTof data posing MS and MS/MS (MS^E) data. Fragment ions (b/y ions) are identified and listed in the Peak Match Data Table.

