

# AN AUTOMATED UPLC/MS DATA ANALYSIS USING BIOPHARMALYNX ; A CASE STUDY FOR CHARACTERIZATION OF THERAPEUTIC INTERFERON PROTEIN

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Joomi Ahn, Asish B. Chakraborty, Ying Qing Yu, Beth L. Gillece-Castro and John Gebler  
Waters Corporation; Milford, MA

## OVERVIEW

- Therapeutic Interferon proteins were analyzed to demonstrate the utility of BiopharmaLynx™ software.
- Fast data process, accurate protein MW determination and peptide mapping to confirm protein sequence were demonstrated.
- Efficient data interpretation was achieved using automated features for peptide peak annotation and data comparison.

## INTRODUCTION

Interferon alpha (IFN α) proteins are a group of genetically related cytokine with a wide spectrum of biological activities including antiviral, antiproliferative, and immunomodulatory properties.<sup>1,2</sup> The goal of the study is to chracterize and compare the recombinant Interferon α-2b proteins harvested from different cell growth conditions.

Intact Interferon protein and the enzymatically digested protein were analyzed using Liquid Chromatography-Time of Flight Mass Spectrometry (LC/TOF-MS). LC/TOF-MS has been extensively used for biotherapeutic drug development. Conventionally, LC/TOF-MS data generated for these studies is interpreted manually or being processed with a limited automation to confirm protein masses and peptide maps. They are often time-consuming processes and offer limited information due to the high complexity nature of the data.

BiopharmaLynx, a new application manager, significantly reduces data analysis time for UPLC/TOF-MS peptide mapping and intact protein analyses. It automatically annotates the peaks from LC/TOF-MS data. This software also automates batch data processing of multiple UPLC/MS experiments. The display tools for data interpretation in user friendly interface are some of the features designed to reduce the burden of data analysis.



Waters ACQUITY UPLC® system with SYNAPT™ MS System and BiopharmaLynx™ Application Manager for MassLynx™ software

## METHODS

### Intact Protein LC and MS conditions

Waters ACQUITY UPLC System  
Waters MassPREP™ Micro Desalting column 2.1 x 5 mm  
Mobile phase A : 0.1% Formic acid in water  
Mobile phase B : 0.1% Formic acid in acetonitrile  
Column Temperature : 80 °C  
Gradient <sup>3</sup> : Total run time in 10 min  
5—10 %B in 0.25 min (0.2 mL/min, waste)  
10—20 %B in 0.01 min ( 0.2 mL/min, MS)  
20—45 %B in 7.1 min ( 0.2 mL/min, MS)  
3-sawtooth 90—5 %B in 0.5 min (0.5 mL/min, wash)

Waters LCT Premier™ ESI-TOF Mass Spectrometry  
Capillary Voltage : 3.2 kV  
Cone Voltage : 40 V  
Desolvation Temperature : 350 °C  
Desolvation Gas : 800 L/hr  
Source Temperature : 150 °C  
Acquisition : 600-3000 m/z

### Peptide Mapping LC and MS Conditions

Waters ACQUITY UPLC System  
Waters ACQUITY UPLC BEH 300 C18, 1.7 µm 300 Å  
Peptide Separation Technology, 2.1 x 150 mm column  
Mobile phase A : 0.1% formic acid in water  
Mobile phase B : 0.1% formic acid in acetonitrile  
Flow : 200 µL/min.  
Gradient : 0-50 %B over 120min.  
Column Temperature : 60 °C  
Weak wash : 95 % buffer A / 5 % buffer B  
Strong wash : 20 % buffer A / 80 % buffer B

Waters SYNAPT ESI-TOF Mass Spectrometry  
Capillary Voltage : 3.0 kV  
Cone Voltage : 35 V  
Desolvation Temperature : 250 °C  
Desolvation Gas : 350 L/hr  
Source Temperature : 120 °C  
Acquisition : 50-1700 m/z for peptide,  
Lockmass : 100 fmol/µL Glu-Fibrinopeptide B directly infused

### Sample preparation

Production batches of recombinant human interferon α-2b protein were used. A control IFN and an analyte IFN were expressed in two different cell conditions in X and Y. For peptide mapping, a control and one of analytes were dissolved in 50 mM ammonium bicarbonate and denatured with 0.025 % RapiGest™ SF (Waters), then reduced and alkylated with DTT and iodoacetamide (Sigma) respectively. Finally the protein was digested with trypsin (Promega). For intact protein analysis, the samples were dissolved in 50 mM ammonium bicarbonate and directly injected into a column.

### BiopharmaLynx 1.1 Application Manager

LC/MS data were processed in BiopharmaLynx 1.1. The screen captures from BiopharmaLynx browser were shown in the results section.

## RESULTS

### INTACT PROTEIN ANALYSIS

Figure 1. Raw spectra of intact IFN control and analyte




Figure 2. Processed spectra in mirror view; The processed spectrum is a MaxEnt1 spectrum




Figure 3. Peak match data table

Chromatogram	Spectrum	Coverage Map	Protein Digests	Peak Match Data
Protein	Modifiers	Calculated Protein Mass (Da)		
IFN-α2b		19265.1543		
IFN-α2b	Oxidation M(1)	19281.1484		
IFN-α2b	Acetyl N-TERM(1)	19307.1641		
IFN-α2b	Acetyl N-TERM(1), Oxidation M(1)	19323.1602		

The intact protein analysis revealed that the control is the expected IFN and its oxidized forms. The analyte was confirmed as acetylated IFN and oxidized form of acetylated IFN.

## PEPTIDE MAPPING ANALYSIS

Figure 4. Automated peak annotation in raw TIC chromatograms in Mirror view

