Beth L. Gillece-Castro, Thomas E. Wheat, Kim Van Tran, Joomi Ahn, Jonathan Turner, Diane Diehl Waters Corp., Milford, MA 01757

INTRODUCTION

Glycosylation is a post-translational modification of proteins. The sugar chains on glycoproteins can mediate biological activity and are, therefore, associated with safety and efficacy attributes of many biopharmaceuticals. The relative amounts of the individual glycan structures must be monitored at all stages of biopharmaceutical research and development.

In the typical HPLC workflow, *N*-linked glycans are released from glycoproteins and are labeled with 2-aminobenzamide (2-AB). The derivatized oligosaccharides have been separated using HPLC in HILIC mode on amidephase chromatography columns. The HPLC columns separate many glycan structures, but with the advent of sub-2 micron particles, and appropriate instrumentation, increased resolution and sensitivity are achieved. The isomeric structures that are present on most glycoprotein drugs can be quantitated separately.

METHODS

Instrument

Pump:	Waters [®] ACQUITY [®] Binary			
-	Solvent Manager			
Injector:	ACQUITY Sample Manager			
Detector:	ACOUITY Fluorescence Detecto			

Typical Chromatography Conditions

49

50

63

Sample: Sample Concentra Sample Volume: Injection Mode:	ition:	2-AB 10 pr 1.5 µ Partia	labeled <i>N</i> -lir nol/µL L al Loop	nked glycan	IS	
Column:		ACQUITY UPLC [®] BEH Glycan, $1.7\mu m$ 2.1 x 150 mm column				
Eluent A: Eluent B: Weak Needle Wash: Strong Needle Wash: Seal Wash:		100 mM Ammonium Formate, pH 4.5 Acetonitrile Acetonitrile/100 mM Ammonium Formate, pH 4.5 (90:10 v/v) Acetonitrile/100 mM Ammonium Formate, pH 4.5 (10:90 v/v) Methanol/Water (10:90 v/v)				
Temperature:		60°C				
Fluorescence:	λex =	330	nm, λem = 4	420 nm		
Gradient:	<u>Time</u> (minu Init 46.5 48	ites)	Flow rate (mL/min) 0.5 0.5 0.25	Buffer A % 25 40 100	Buffer B % 75 60 0	

0.25

0.5

0.5

100

25

25



Typical Sample Preparation

<u>Hydrophilic Interaction Chromatography</u> (HILIC)



Figure 1. The primary mechanism of hydrophilic interaction liquid chromatography (HILIC) is the partitioning of a polar analyte between bulk mobile phase and the immobilized water layer. Analytes may also hydrogen bond to the amide stationary phase directly or interact by an ion-exchange mechanism.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

75

75

A NOVEL COLUMN FOR ANALYSES OF *N*-LINKED GLYCANS OF GLYCOPROTEINS

RESULTS







Figure 3. The 2-AB labeled glycans released from pooled human IgGs are typical of the analytical challenge associated with biopharmaceuticals. The N-linked glycans in this mixture include the high mannose, neutral, and sialylated complex structures shown above.



Figure 4. 2-AB labeled IgG glycans are typically analyzed with ammonium formate buffer. Changing the mobile phase modifier from ammonium formate to formic acid may be used for LC/MS applications as long as high ionic strength is maintained for retention sialylated (acidic) glycans.



Figure 5. Decreasing the ionic strength of ammonium formate results in decreased retention of glycans with particularly significant loss of retention and yield of sialylated (acidic) glycan structures.



Figure 6. The 2-AB labeled high mannose glycans released from bovine ribonuclease B were separated under similar conditions as the IgG glycans. The Man7 region (zoom) shows clearly that there are three isomers present.





Figure 7B. Changing the temperature of the glycan column separation results in significant changes in backpressure.







Figure 7A. 2-AB labeled IgG glycans were analyzed at increasing temperatures. Increasing the column temperature to 60°C narrowed peak widths and decreased retention without significant effects on selectivity.

Figure 8. If highest resolution is not required, shorter columns and much shorter run times are possible. In all cases the Man5 glycan is fully separated from GOF, and the G1F isomers can be quantitated.



Figure 9. In this example, sialylated 2-AB fetuin oligosaccharides were separated using the same mobile phase and gradient conditions that were used for the IgG glycans. Note that the isomeric triantennary structures (*) were particularly well-resolved.



Figure 10. Neutral 2-AB labeled oligosaccharides, represented by this series of homologous glucose oligomers, also separate on the ACQUITY UPLC BEH Glycan column. The separation gives good peak shape and identification from 1 to more than 22 Glucose Units.

CONCLUSION

- Waters has developed a glycan analysis solution using the new ACQUITY UPLC BEH Glycan column containing 1.7µm particles.
- The Waters ACQUITY UPLC Glycan Separation Technology columns when operated with the ACQUITY UPLC instrument and fluorescence detector provide a high resolution and rapid method for separating and profiling glycans.
- The 2-AB labeled glycans are separated in HILIC mode on an amide stationary phase.
- Positional and linkage isomers can be resolved.
- The same column and mobile phase can be used for 2-AB labeled neutral and charged oligosaccharides such as highly sialylated glycans.