

## DETERMINATION OF MELAMINE RESIDUE IN WATER SAMPLES BY UPLC/MS/MS

Wang Jing, Liu Zhengzheng (Environmental Monitoring Center of Zhejiang Province)

The melamine contamination of milk that occurred in 2008 was undoubtedly a serious food safety situation in China as safe drinking water is an important part of food safety. Since both the drinking water and the food chain may cause a direct or indirect impact on human health, it is urgent to carry out the study on the analysis method of melamine residue in water samples.

### EXPERIMENTAL

#### Equipment and Reagents

An Ultra Performance Liquid Chromatography/tandem mass spectrometry (ACQUITY®/Quattro Premier™); MassLynx™ 4.1 workstation; solid-phase extraction (Zymark); Waters Oasis® MCX SPE column (200 mg, 6 mL); pressure blowing concentrator (Organomation).

Methanol and acetonitrile are both in liquid chromatography pure (Fisher); water is obtained from a Millipore water purification system; ammonium acetate is in guaranteed grade (Merck), hydrochloric acid and ammonia are in analytical pure (purchased from East China Pharmaceutical Company); and the standard reference is purchased from J&K.

#### Sample of the Pretreatment Methods

**Solid-phase extraction:** Add 5 mL of methanol and 5 mL water, at the speed of 5 mL/min to the column to make it active.  
**Acidification of water sample:** Add 0.25 mL of concentrated hydrochloric acid to 500 mL of water. Add the water sample with the speed of 4 mL/min and the volume of 500 mL. Rinse the sample with 5 mL of water and methanol respectively. Elute with 5 mL of methanol (containing 5% of ammonia) twice and concentrate the eluent with a pressure blowing concentrator until it is almost dried up. Dilute to 1 mL with mobile phase and filter it. Inject.

### CHROMATOGRAPHY CONDITIONS

UPLC®/MS method: column: UPLC column (ACQUITY UPLC® BEH HILIC, 1.7  $\mu$ m, 2.1 50 mm); column temperature: 40 °C; mobile phase: a mixture of acetonitrile and water (containing 10 mM of ammonium acetate) (90:10); flow speed: 0.4 mL/min.

The target contaminator is quantitatively analyzed with the MRM mode. Detect the target with ESI<sup>+</sup>. Capillary voltage: 3.0 KV, source temperature: 120 °C, the temperature of desolvent gas: 400 °C, flow speed: 800 L/h, cone gas flow: 50 L/h; argon flow speed: 0.38 mL/min at the MRM mode; cone current: 40 V, CID: 17; characteristic ions pair: 127>85, as in the following Figure 1.




Figure 1. Mass spectra of melamine.

### PERFORMANCE

#### Accuracy

Add 6.0  $\mu$ g of target substance to 500 mL of testing water, and the recovery of a series of pre-sampling is  $101.8 \pm 10.7\%$  ( $n=6$ ), as showed in Table 1.

#### Precision

Inject 1.0  $\mu$ g/mL of standard solution for 6 times repeatedly and the RSD of peak area is 3.2%.

## Detective limit

The signal/noise ratio is 15 for 1.0 ng/mL of standard solution and 3 for 0.2 ng/mL of standard solution, as shown in Figure 2. Taking the pre-sampling procedure into consideration, the quantitative detective limit of this method is 0.4 ng/L.

## Linearity

Prepare 4 standard solutions with different concentrations ranging from 1.2 ng/mL to 2400 ng/mL, the equation for the standard curve is  $y=1350.3x+161822$ ,  $R=0.9925$ .

## Blank

The blank is lower than the detective limit.

## Determination of sample

The water sample is obtained from a reservoir in Zhejiang province with a concentration of melamine of <0.4 ng/L. Add some standard solution to this water sample and analyze the content of melamine. The concentration of melamine is 0.03 µg/L and the total ion flow is showed in Figure 3.




Figure 3. TIC of testing sample.

## CONCLUSIONS

The UPLC/MS/MS method used for the determination of melamine residue in water sample has been proved to be practicable. The analyses can be fulfilled in 1.0 min with a recovery of  $101.8 \pm 10.7\%$  and a precision of 3.2%. With the range covering three magnitudes, the linear coefficient is greater than 0.99 and the detective limit is up to 0.4 ng/L. All of these indicated that this method meets the requirement of practical analyses.

Analytes	Precision % (n=6)	Recover% (n=6)	Detective limit (ng/L)	Linearity
Melamine	3.2	101.8±10.7	0.4	$y = 1350.3x + 161822$ , $r=0.9925$

Figure 1. Performance of this method.




Figure 2. The S/N of melamine in different concentration.

## Sales Offices

**Austria and European Export**  
(Central South Eastern Europe, CIS  
and Middle East) 43 1 877 18 07

**Australia** 61 2 9933 1777

**Belgium** 32 2 726 1000

**Brazil** 55 11 5094-3788

**Canada** 1 800 252 4752 x2205

**China** 86 21 6879 5888

**CIS/Russia** +497 727 4490/290 9737

**Czech Republic** 420 2 617 1 1384

**Denmark** 45 46 59 8080

**Finland** 09 5659 6288

**France** 33 1 30 48 72 00

**Germany** 49 6196 400600

**Hong Kong** 852 2964 1800

**Hungary** 36 1 350 5086

**India and India Subcontinent**  
91 80 2837 1900

**Ireland** 353 1 448 1500

**Italy** 39 02 265 0983

**Japan** 81 3 3471 7191

**Korea** 82 2 6300 4800

**Mexico** 52 55 5524 7636

**The Netherlands** 31 76 508 7200

**Norway** 47 6 384 60 50

**Poland** 48 22 6393000

**Puerto Rico** 1 787 747 8445

**Singapore** 65 6593 7100

**Spain** 34 936 009 300

**Sweden** 46 8 555 11 500

**Switzerland** 41 56 676 70 00

**Taiwan** 886 2 2543 1898

**United Kingdom** 44 208 238 6100

**All other countries:**

**Waters Corporation U.S.A.**

1 508 478 2000

1 800 252 4752

[www.waters.com](http://www.waters.com)

# Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



© 2010 Waters Corporation. Waters, The Science of What's Possible, ACQUITY, ACQUITY UPLC, Quattro Premire, MassLynx, Oasis, and UPLC are trademarks of Waters Corporation

©2010 Waters Corporation. Produced in the U.S.A.  
April 2010 720003440EN VW-PDF

**Waters Corporation**  
34 Maple Street  
Milford, MA 01757 U.S.A.  
T: 1 508 478 2000  
F: 1 508 872 1990  
[www.waters.com](http://www.waters.com)