# MAPPING HOST CELL PROTEIN VARIATIONS IN MONOCLONAL ANTIBODY SAMPLES DERIVED FROM DIFFERENT PURIFICATION PROTOCOLS

Catalin E. Doneanu<sup>1</sup>, <u>Weibin Chen<sup>1</sup></u>, Alex Xenopoulos<sup>2</sup>, Holly Prentice<sup>2</sup>, StJohn Skilton<sup>1</sup> and Jeff Mazzeo<sup>1</sup> <sup>1</sup> Waters Corporation, Milford, MA 01757; <sup>2</sup> Millipore Corporation, Bedford, MA, USA

# **OVERVIEW**

A generic, efficient UPLC/MS<sup>E</sup> assay for identification and quantification of HCPs over 4 orders of magnitude in concentration is presented.

# INTRODUCTION

- Residual host cell proteins (HCPs) from recombinant production are a major component of biopharmaceutical process-related impurities. HCPs can elicit an unpredictable immune response in patients
- The composition of HCPs is extremely heterogeneous and changes with production and purification procedures.
- Current analytical methods for measuring HCPs, (typically ELISA, blots), are expensive, subjective, time-consuming to develop, and require prior knowledge about the identities of the contaminant proteins.

# WORKFLOW

- Spike a mixture of 5 proteins in the concentration range of 10-1,000 ppm (ng/mg) in the protein biopharmaceutical
- Enzymatic digestion of the protein sample to peptides
- 2D-LC/MS<sup>E</sup> with IDENTITY<sup>E</sup> to discover HCPs
- Use the top 3 peptides from each protein for absolute label-free quantitation based on the known concentration of the spiked proteins
- Use VERIFY<sup>E</sup> for mining the MS<sup>E</sup> dataset to extract appropriate MRMs for monitoring HCPs by Xevo TQ MS
- Analyze the sample digests by tandem quadrupole MS for targeted, high-throughput quantitation (e.g. using isotopically labeled peptides)

# **METHODS**

#### Sample Preparation

A chimeric anti-phosphotyrosine IgG1 mAb (PTG1) was expressed in two CHO cell lines (CHO-S and DG44) and purified by Protein A chromatography using two different protocols. Same cell lines, expressing no PTG1 protein, were also grown under the same conditions and purified on a Protein A column. Five protein standards were spiked in PTG1 and the resulting protein mixture was denatured with 0.1% RapiGest (15 min at 60 °C), reduced with 10 mM dithiothreitol (DTT) for 30 min at 60°C, alkylated with 20 mM iodoacetamide (IAM) for 30 min (at RT) and enzymatically digested with Promega trypsin (1:20 w/w ratio) at 37°C overnight. Following digestion, the protein amounts loaded on-column on the 2D-LC system were: 4,000 fmoles LA (bovine alpha-lactoglobulin), 800 fmoles PHO (rabbit glycogen phosphorylase b), 320 fmoles ADH (yeast alcohol dehydrogenase), 80 fmoles BSA (bovine serum albumin), and 16 fmoles ENL (yeast enolase).

#### LC Conditions

A Waters nanoACQUITY UPLC<sup>®</sup> system with 2D technology was used for peptide separations. A reversed-phase/reversed-phase (RP/RP) method was developed that uses the pH of the mobile phases to change the selectivity of a peptide separation in two separate dimensions [1,2]

- First Dimension (1D) pH=10: 1.0 x 50 mm XBridge<sup>™</sup> C<sub>18</sub> column (5 µm particles); flow 10 µL/min. Mobile phase 20 mM ammonium formate in water (Solvent A) and ACN (Solvent B).
- Online dilution (1:10) of the eluent from 1D before analyte trapping onto the 2D column.
- Trap column: A 5-μm Symmetry C<sub>18</sub> trap (2.7 μl volume) was used to trap peptides between the two LC dimensions.
- Second Dimension (2D) pH=2.4: 0.3 mm x 150 mm analytical column BEH C<sub>18</sub> 1.7 µm, kept at 65 °C and operated at 12 µL/ min
- Fractions were eluted in ten steps (1: 10.8% Eluent B, 2: 12.4%,
- 3: 14.0%, 4: 15.4%, 5: 16.7%, 6: 18.6%, 7: 20.4%, 8: 25.0%,

9: 30.0%, 10: 50.0% B. Each step was mixed in a 1:10 ratio with 0.1% TFA in water (pH=2.1) before trapping. Low pH separations in the second chromatographic dimension used a 30 min gradient from 7 to 35% acetonitrile (0.1% FA-formic acid). See Figure 1 for the diagram of the 2D-LC system.

#### MS conditions

A) Discovery and Label-free Quantitation: Data independent, alternate scanning LC/MS<sup>E</sup> experiments were performed with a

Waters SYNAPT<sup>®</sup> HDMS<sup>™</sup> mass spectrometer.

- Acquisition time was 0.5 sec
- m/z range: 50-1990
- Fixed CE at 5 eV for MS; CE ramp 15 35 eV for high-energy fragmentation
- Database search—ProteinLynx Global Server (PLGS) 2.4 with Identity<sup>E</sup> Informatics [3].

B) Targeted Quantitation: MRM analyses were performed on a Waters Xevo<sup>®</sup> TQ MS tandem guadrupole mass spectrometer coupled with the ACQUITY UPLC<sup>®</sup> system.

- 2.1 x 150 mm BEH130  $C_{18}$  column (1.7  $\mu$ m particles) at 35 °C.
- Flow rate of 300 µL/min; 3-40% B over 10 min (B: acetonitrile with 0.1% FA; A: water with 0.1% FA)
- Other parameters: Potential 3.5 kV, CV 37 V, source 90 °C, MS1/ MS2 isolation window 0.75 Da (FWHM), adjustable collision energy depending of precursors/fragment (15-30 eV) and 10 ms dwell time.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS







Figure 1: Fluidic configuration for 2-dimensiona chromatography with on-line dilution: (A) Sample loading; (B) Step elution of peptides from 1D (in 10 steps) followed by trapping; (C) Peptide separation in the second dimension.



Figure 2. SYNAPT HDMS mass spectrometer coupled to the 2D nanoACQUITY UPLC system.



Figure 3. Reproducibility of 2D chromatography: extracted mass chromatograms of the T43 ENL peptide (VNQIGTLSESIK,  $[M+2H]^{2+}= 644.86$ ) from four consecutive 2D-LC experiments. All second dimension chromatography runs were performed at 4 µL/min using a 90 min gradient (3-40% ACN, 0.1% FA). T43 peptide eluted only in Fraction 3/5 (16.7% B) in the 1st dimension chromatography. The amount of ENL digest loaded on column in each 2D experiment was 24 fmoles.



Figure 4. Chromatographic performance (e.a. reproducibility and peak width) is maintained during 1D fractionation: mass chromatograms of T43 peptide obtained under four fractionation conditions:

(A) "simulated" 1D run using a single elution step (50% ACN) (B) fraction 3 out of 5 (2.7% B step elution) (C) fraction 5 out of 10 (1.3% B step elution) (D) fraction 7 out of 20 (0.6% B step elution)



*Figure 5. Typical second dimension chromatography separation* demonstrating a peak capacity of 160. When coupled with extensive 1D fractionation (eg. 20 fractions), the 2D-LC setup can achieve very high peak capacities (3,200 in this case)







### **HCP Identification Method**

• LC-MS<sup>E</sup> Data: a combined MS<sup>E</sup> dataset collected from 10step 2D RP/RP separations performed in triplicate.

 Protein Identification Informatics: ProteinLynx Globlal SERVER<sup>™</sup> (PLGS) 2.4 featuring Identity<sup>E</sup> Bioinformatics.

• **Database**: protein database was compiled from 12,943 mouse Swiss Prot entries, 654 Golden hamster proteins the sequences of LA, ADH, PHO, BSA, ENL, porcine trypsin, protein A (*S aureus*) and the heavy and light chain sequence of PTG1. The final randomized database contained 27,212 entries.

Table I. HCP concentrations (ng/mL) measured in PTG1 mAb expressed in CHO (samples I-IV) and DG-44 cells (samples A-D). Two protein A purification protocols were compared across ten samples: I vs II - CHO cells grown under regular conditions; III vs IV - CHO cells non-expressing the PTG1 target; A1/2 vs B1/2 biological replicates for DG-44 cells grown under regular conditions; C vs D - non-expressing DG-44 cells. High abundance proteins are highlighted in RED (>1,000 ng/mL), medium abundance are either in YELLOW (500-1000 ng/mL) or GREEN (100-500 ng/mL) and low-abundance HCPs are GREY (< 100 ng/mL).

Cell	PTG1	No of HCPs	Common	HCP concentration (ng/mL)	
line	expressed?	Identified	HCPs	Mass Spec	ELISA
CHO	Y	5	4	4,870	7,013
CHO	Y	14	4	9,752	12,228
CHO	Ν	0	0	0	30
СНО	Ν	31	0	9,093	4,138
DG-44	Y	34	21	61,138	3,684
DG-44	Y	50	21	109,437	8,042
DG-44	Y	22	15	35,031	4,282
DG-44	Y	32	15	105,123	8,041
DG-44	Ν	0	0	0	6
DG-44	N	28	0	1 454	1 318

Table II. Comparison of total HCP concentrations (ng/mL) measured by the mass spectrometry based assay and ELISA.

The MS<sup>E</sup> data was automatically mined using VERIFY<sup>E</sup> for the best proteotypic peptides for each MIX-5 spiked protein. MRM assays comprising 2 peptides/protein and two transitions per peptide (for a total of 20 MRMs) were then performed on a Xevo TQ mass spectrometer with a 10-min LC separation.



Figure 6. Reproducibility of the MRM assay: three replicates showing the MRM chromatograms for an MRM transition of the T4 peptide (TGNPTVEVELTTEK) from ENL (200 fmoles ENL digest loaded on a 2.1 x 150 mm BEH130 C18 column (1.7 µm particles) operated at a flow rate of 300 µL/min using a short 10 min gradient (from 3% to 40% ACN, 0.1% FA)

# CONCLUSIONS

- The LC/MS<sup>L</sup>-based assay allows the identification and quantification of low-abundance HCP contaminants in biopharmaceuticals over four orders of magnitude in concentration
- Protein A purification of mAbs using different purification protocols produces different HCP patterns. Cell lines used for mAb purification can affect the HCPs identity and concentration significantly
- ELISA assay seems to underestimate total HCP concentration
- A high-throughput MRM assay for monitoring and quantification of HCPs using a tandem quadrupole mass spectrometer can be easily implemented.

#### References

- 1. Gilar M, Olivova P, Daly AE, Gebler JC J. Sep. Sci, 2005, 1694.
- 2. Stapels MD, Fadgen K Current Trends in Mass Spectrometry, Spectroscopy supplement, March 2009.
- 3. Li GZ, Vissers JP, Silva JC, Golick D, Gorenstein MV, Geromanos SJ Proteomics, 2009, 1696