# APPLYING LC-MS WORKFLOWS TOWARDS DEVELOPING A WELL-CHARACTERIZED RECOMBINANT SUBUNIT VACCINE CANDIDATE

Hongwei Xie<sup>1</sup>; Catalin Doneanu<sup>1</sup>; Weibin Chen<sup>1</sup>; St John Skilton<sup>1</sup>; Joseph Rininger<sup>2</sup>; Martin Gilar<sup>1</sup>; Jeff Mazzeo<sup>1</sup>. 1: Waters Corporation, 34 Maple Street, Milford, MA 01757. 2: Protein Sciences, Meriden, CT

# CHARACTERIZATION BY PEPTIDE MAPPING WITH RP LC-MS<sup>E</sup>

- Label-free quantification of rHA antigens and processrelated impurity proteins with ProteinLynx Global SERVER™ (PLGS) Software for automated data processing<sup>1</sup> (Table 1);
- Confirmation of sequences of rHA/B proteins with BiopharmaLynx<sup>™</sup> for automated data processing<sup>2</sup> (Figure 1);
- Confirmation of site-specific glycosylation and glycoforms on each rHA antigen (Diagram 1, Figures 2-4 and 7, Table 2);
- Identification of degradations on each rHA antigen (Fig. 5).



Diagram 1: Algorithm for the proposal of glycosylation in PLGS 2.4. *The algorithm uses accurate masses to propose glycan masses* matching identified peptides to within a reconstructed spectrum.



# LC-MS<sup>E</sup> METHODOLOGY

**Sample:** 4-hour tryptic digest of a research vaccine formulated with purified recombinant hemaglutinins rH1, rH3 and rB.

**LC Conditions:** Waters ACQUITY UPLC<sup>®</sup> System. 120-min gradient (0 to 40%) ACN in 0.1% FA) on a 2.1 x 150 mm BEH C18 1.7 µm column at a flow rate of 0.2 mL/min and a column temperature of 60 °C.

**MS Conditions**: Waters SYNAPT<sup>®</sup> HDMS<sup>™</sup> System. MS<sup>E</sup> data acquired in ESI positive ion mode, collision cell energy alternating between low energy (5 V, for peptide precursor MS data) and elevated energy (rampping from 20 to 40 V, for fragments MS<sup>E</sup> data).

Data processing: Protein identification and absolute quantification by searching against a combined sequence database with PLGS 2.4. Sequence verification, mining of site-specific PTMs (N-linked glycosylation and degradations) by BiopharmaLynx 1.2. Strict tryptic cleavage rule was applied.

Protein	Organism	M.W.	PLGS Score	Unique Peptides	Normalised Conc (%)*	R (%		
H1-H1N1A/Brisbane/59/2007	rHA	61229	2166	45	64.8	14.9		
H3-H3N2A/Brisbane/16/2007	rHA	61814	675	29	35.8	6.2		
B-B/Florida/4/2006	rHA	61554	2606	44	100	9.9		
P17501-Major envelope glycoprotein	A. Californica	58566	1410	16	6.7	12.1		
P32651-Structural glycoprotein gp41	A. Californica	45381	204	7	2.7	16.9		
P41678-Capsid protein p24	A. Californica	22110	223	3	2.3	3.9		
Proteins Identified by Homology:					2.5			
Q05825-ATP synthase subunit beta	Drosophila	54108	555	8	2.5	14.8		
(fragment)	Drosophila	34850	463		2.2	15.6		
BI1502-GM10439	Drosophila	36921	161	3	1.9	4.0		
Q5KMQ8-Expressed protein (proposed uncharacterized protein)	Yeast	30298	94	4	6.3	4.5		
* Normalized to HA protein B and averaged from 4 repeated injections; Calculated using sum of the intensities of 3 most abundant identified unique peptides (ref 1)								

Table 1: Proteins identified from the tryptic digest of the recombinant vaccine candidate by PLGS 2.4. Protein concentrations were calculated using the intensity sum of the 3 most abundant unique peptides of each protein, and normalized to the most abundant protein (rB).



Figure 2: Glycan-related fragment ion spectrum of a glycosylated tryptic peptide rH3\_T21 proposed by PLGS2.4. The 'Y1' (3071 Da) and precursor ions (4732 Da) are at different relative abundances and scaled separately. 8 of the 11 theoretical possibilities are automatically proposed as 'gd' [glycan difference], providing an almost complete fragment ladder with no pre-knowledge of the glycosylations, nor of the protein sequence.



system<sup>3</sup>. A) XIC chromatograms of tryptic peptide rB\_T8 using Figure 5: Separation and quantification of unmodified and modified accurate masses. B)  $MS^{\epsilon}$  spectrum of unmodified rB\_T8. C)  $MS^{\epsilon}$ tryptic peptide rH1\_T22 (shown as XIC chromatogram, but originally spectrum of glycosylated rB\_T8 by BiopharmaLynx 1.2 with identified and tabulated by BiopharmaLynx 1.2.) proposed glycan assignments based on accurate masses.

# TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

Figure 1: Sequence coverage of rB determined from the LC- $MS^{E}$  data (98.8%) and displayed by BiopharmaLynx 1.2. Similarly high sequence coverage for rH1 and rH3 was obtained.

# INTRODUCTION

- Recombinant vaccines have emerged as an alternative to eggbased production systems.
- The need for annual re-evaluation of influenza vaccines for efficacy calls for appropriate methods for characterization of hemaglutinin (HA) antigens and final vaccine products.
- Better characterization may help speed up vaccine development and reduce safety concerns.
- This study describes advanced LC-MS methods for characterization (by LC-MS<sup>E</sup>) and monitoring (by LC-MRM) of a trivalent influenza vaccine candidate formulated from purified recombinant HA proteins (rH1, rH3, and rB corresponding to influenza A subtypes H1N1 and H3N2 and influenza B viruses, respectively).



Diagram 2: An illustration of the different analyses possible from the same sample preparation starting with the formulated vaccine. Paths towards complete characterization are possible on tof-based instruments, and targeted analyses are performed on tandem guadrupole instruments

# CONCLUSIONS

- Peptide mapping by RP LC-MS<sup>E</sup> was used to:
  - Verify rHA primary sequences with high sequence coverage,
  - Discover site-specific glycosylation for each rHA antigen;
  - Examine potential degradations on each rHA antigen;
  - Quantify rHA antigens;
  - Identify and quantify impurity proteins.
- LC-MRM was used to confirm certain critical attributes of vaccines:
  - ◆ Reverse phase (RP) LC-MRM was used to quantify rHA antigens
  - HILIC MRM was used to confirm site-specific glycoforms that were poorly separated by RP, and low abundance glycoforms without supporting MS<sup>E</sup> spectra.
- LC-MS<sup>E</sup> and LC-MRM are complementary analytical tools and contribute to assays appropriate for a well-characterized vaccine.

# THE SCIENCE OF WHAT'S POSSIBLE.™

# RAPID QUANTIFICATION AND MONITORING BY LC-MRM

Quantification of rHA antigens and ratios of the antigens by RP LC-MRM (Figure 6);



*Figure 6: Calibration curve for quantification of rH3. Two* additional MRM transitions were used to confirm and quantify the peptide (both native and isotopically labeled forms). The concentration of rH3 was calculated to be 7.9  $\pm$  0.8 pmole/ml. Similar methods were used to quantify rH1 and rB. The concentration ratio obtained between the antigens is consistent with the result quantified by  $RP LC-MS^{E}$ .

(R)DNVSCSICL

ntaining -GlcNAC[Fuc]-GlcNAC-(Man)x; X-X or XF-XF - peptide with two glycans on two differe



Figure 7: HILIC LC-MRM separation and confirmation of rH3\_T21 glycoforms (C). The glycoforms were detected by RP LC-MS<sup>E</sup> and BiopharmaLynx 1.2 (B), but have no separation (A) and except for the most abundant glycoform (rH3\_T21-GlcNAc2Man9), no  $MS^{E}$  spectra were available for the confirmation of less-abundant glycoforms.

Protein	Peptide	Residue	Peptide Sequence containing -NXS/T- motif)	Peptide	Modification Site	Proposed	Motif
		Numbers		IVI VV		Giycolomis	
rH1	TI	1-22	(-)DTICIGYHANNSTDTVDTVLEK(N)	2408.12	NII	3,3F,5,5F,7,8	NNSTD
	T2	23-40	(K)NVTVTHSVNLLENSHNGK(L)	1961.99	N23	3F,6,8,9	(K)NVTV
			(K)GIAPLQLGNCSVAGWILGNPECELLISK(E)	2894.50		0,3,3F, 5,8	GNCSV
			(K)ESWSYIVEKPNPENGTCYPGHFADYEELR(E)			2,3,5,6,7,8,9	ENGTC
	Т8	120-145	(K)ESSWPNHTVTGVSASCSHNGESSFYR(N)	2825.21	N125	1,1F,2,2F,3,3F	PNHTV
	T10	154-162	(K)NGLYPNLSK(S)	1004.53	N159	0,3,3F,5,6,7,8,9	PNLSK
	T24	278-304	(K)CQTPQGAINSSLPFQNVHPVTIGECPK(Y)	2864.39	N286	3,3F,5,8	INSSL
	T47	480-487	(K)NGTYDYPK(Y)	956.42	N480	3,3F,5,7,8,9	(K)NGTY
	T52	504-545	(K)LESMGVYQILAIYSTVASSLVLLVSLGAISFWMCSNGSLQCR(I)	4509.28		nd	SNGSL
rH3	T2	3-27	(K)LPGNDNSTATLCLGHHAVPNGTIVK(T)	2528.28	N8-N22	2F-2F	DNSTA/PNGTI
	T3	28-82	(K)TITNDQIEVTNATELVQSSSTGEICDSPHQILDGENCTLIDALLG DPQCDGFQNK(K)		N38-N63		TNATE/ENCTL
	Т8	110-141	(R)SLVASSGTLEFNNESFNWTGVTQNGTSSACIR(R)	3376.56	N122/N126/N133	6-6, 7-7, 7-8	NNESF/FNWTG/QNGTS
	T10	143-150	(R)SNNSFFSR(L)	957.43	N144	0,3,3F,5,7,8,9	SNNSF
	T13	161-173	(K)YPALNVTMPNNEK(F)	1489.72		3,3F,5,6,9	LNVTM
	T21	230-255	(R)ISIYWTIVKPGDILLINSTGNLIAPR(G)	2866.63	N246	5,6,7,8,9	INSTG
	T27	277-299	(K)CNSECITPNGSIPNDKPFQNVNR(I)	2546.16	N285	3F,6,7,8,9	PNGSI
	T52	483-492	(R)NGTYDHDVYR(D)	1238.53	N483	3F,5,5F,6,8	(R)NGTY
rB	T3	18-38	(K)TATQGEVNVTGVIPLTTTPTK(S)	2127.14	N25	3,3F,5,6,7,8,9	VNVTG
	T8	53-80	(K)LCPDCLNCTDLDVALGRPMCVGTTPSAK(A)	2892.33	N59	0,3F	LNCTD
	T16	137-149	(R)LGTSGSCPNATSK(S)	1221.57	N145	3,3F,5,6,7,8,9	PNATS
	T19	167-197	(K)NATNPLTVEVPYICTEGEDQITVWGFHSDDK(T)	3477.60	N167	3,7,8,9	(K)NATN
	T30	299-304	(K)YGGLNK(S)	650.34	N303	3,3F,5,6,7,8,9	LNK(SK)
	T34	330-335	(K)LANGTK(Y)	602.34	N332	0,3,5,6,7,8,9	ANGTK
	T51	490-497	(K)CNQTCLDR(I)	951.39		nd	CNQTC
	T52	498-560	(R)IAAGTFNAGEFSLPTFDSLNITAASLNDDGLDNHTILLYYSTAA	6699.38		nd	LNITA/DNHTI

Table 2: Site-specific glycoforms detected by LC-MS<sup> $\pm$ </sup> and BiopharmaLynx 1.2, and verified by HILIC LC-MRM (except for those glycopeptides with multiple modification sites) for the 3 rHA antigens.

# LC-MRM METHODOLOGY

**ACQUITY UPLC Conditions:** For quantification of rHA proteins, a 30-min gradient (0 to 35% ACN in 0.1% FA) on a 2.1 x 150 mm BEH C18 1.7 µm column at a flow rate of 0.3 mL/min and a column temperature of 35 °C; For monitoring and verification of site-specific glycoforms, a 60-min gradient (90 to 70% B in 10 min, then to 50% B in 50 min) on a 2.1 x 150 mm BEH glycan 1.7 μm column at a flow rate of 0.2 mL/min and a column temperature of 60 °C. B—10 mM NH4-formate in 90/10 CAN/water and A—10 mM NH4-formate in water. MS Conditions: Waters<sup>®</sup> Xevo TQ MS. MRM acquisition was implemented in an ESI+ mode, Capillary voltage 3.5 V, Cone voltage 35 V, Dwell time 5 ms, collision energy was optimized for each transition

**Data process:** QuanLynx for chromatographic peak integration and obtained peak areas were exported to MS Excel for subsequent analysis.

## REFERENCE

I.Xie HW, Gilar M, Gebler JC. "Characterization of protein impurities and site-specific modifications using peptide mapping with liquid chromatography and data independent acquisition mass spectrometry", Anal. Chem., 2009, 81:5699-5708

2.Xie HW, Chakraborty A, Ahn J, Yu YQ, Dakshinamoorthy DP, Gilar M, Chen WB, Skilton SJ, Mazzeo JR. "Rapid comparison of a candidate biosimilar to an novator monoclonal antibody with advanced liquid chromatography and mass spectrometry technologies", mAbs, 2010, 2:379-394

3. Tomiya NT, Narang S, Lee YC, Betenhaugh MJ. "Comparing N-glycan processing in mammalian cell lines and engineered lepidopteran insect cell lines Glycoconjugate J. 2004; 21:343-360

©2010 Waters

