

# Eliminating Phospholipids in Drug Discovery Extractions Using a Fast, Generic Sample Clean-up Method

Jessalynn P. Wheaton, Erin E. Chambers, John Martin and Kenneth J. Fountain

## APPLICATION BENEFITS

- Simple, universal sample prep protocol
- Speeds up workflow through direct injection of eluates
- Eliminates the vast majority (>99%) of plasma phospholipids
- Reduces sample variability, eliminating a major source of suppression
- Facilitates use of shorter runtimes, improving throughput

## WATERS SOLUTIONS

- ACQUITY UltraPerformance LC® technology
- Xevo® TQ MS system
- Ostro™ sample preparation products
- IntelliStart™ software

## KEY WORDS

Phospholipid removal, LC/MS/MS, Ostro, matrix effects, drug discovery, screening

## INTRODUCTION

Drug discovery is a vital segment of pharmaceutical research where vast numbers of compounds are screened to determine therapeutic efficacy, activity, and ADME properties. This process helps identify the handful of drug candidates that will progress further. Many closely related drug compounds must be rapidly analyzed and quick decisions must be made as to which drugs will continue into development and eventually clinical trials. During the drug discovery stage, speed, time, ease of use, and high throughput are key aspects of everyday work. There is little time for method development, making simple and universal sample prep methods such as protein precipitation (PPT) an attractive choice. Crude techniques such as PPT are often quite efficient in terms of generating high analyte recovery but result in relatively dirty samples. In particular, PPT does little to eliminate phospholipids (PLs), a major source of concern in bioanalysis. PLs build up in LC/MS/MS systems and are one of the major sources of matrix effects in plasma-based assays. Amongst other problems, matrix effects also alter mass spectrometry response in an unpredictable manner, decrease method robustness, and add to method variability. In this publication, Ostro™ 96-well sample preparation plates are used to eliminate both proteins and the vast majority of PLs while maintaining high analyte recovery, all with a simple single step method. A screening method for a group of 26 structural analogs and metabolites (see Figure 1 for representative structures) in plasma was developed using this technique.




Figure 1. Representative structures of 6 of the 26 analogs utilized in this study.

## EXPERIMENTAL

## ACQUITY UPLC Conditions

Column:	ACQUITY UPLC® BEH C <sub>18</sub> , 2.1 x 50 mm, 1.7 μm			
Mobile Phase A:	0.1% HCOOH in H <sub>2</sub> O			
Mobile Phase B:	Methanol			
Flow Rate:	0.6 mL/min			
Gradient:	Time (min)	Profile %A %B	Curve	
	0.0	98 2	6	
	2.0	1 99	6	
	2.5	1 99	6	
	2.6	98 2	6	
	3.0	98 2	6	
Injection Volume:	18.0 μL			
Injection Mode:	Partial Loop			
Column Temperature:	35 °C			
Sample Temperature:	15 °C			
Strong Needle Wash:	60:40 ACN:IPA + 0.2% conc. HCOOH (600 μL)			
Weak Needle Wash:	80/20 H <sub>2</sub> O/MeOH (200 μL)			

Waters Xevo™ TQ MS Conditions, ESI<sup>+</sup>

Capillary Voltage:	1.0 V
Desolvation Temp:	400 °C
Cone Gas Flow:	Not used
Desolvation Gas Flow:	1000 L/Hr
Collision Cell Pressure:	2.6 x 10 <sup>(-3)</sup> mbar
MRM transition monitored, ESI <sup>+</sup> :	See Table 1

Table 1. MRM transitions and MS conditions for the 26 structural analogs in plasma, generated using IntelliStart<sup>1</sup>

Compound Name	Precursor Ion	Product Ion	Cone Voltage (V)	Collision Energy (eV)
Triazolam	343.1	238.9	42	40
Alpha-hydroxymidazolam	342.0	203.0	38	28
2-hydroxyethylflurazepam	333.1	108.9	36	26
Clozapine	327.2	269.9	50	24
Midazolam	326.1	290.9	42	26
Prazepam	325.2	270.9	32	22
Alpha-hydroxyalprazolam	325.1	296.9	40	24
Bromazepam	318.1	181.9	38	32
Clonazepam	316.1	269.9	38	24
Flunitrazepam	314.1	268.0	42	24
Alprazolam	309.1	280.9	42	26
Temazepam	301.1	254.9	26	24
Clobazam	301.1	223.9	32	30
n-Desmethyflunitrazepam	300.2	253.9	38	22
Chlordiazepoxide	300.2	226.9	28	22
Estazolam	295.1	266.9	34	22
Desalkylflurazepam	289.1	139.9	42	26
Oxazepam	287.1	241.0	28	32
7-aminoflunitrazepam	284.2	135.0	46	26
Nitrazepam	282.2	207.0	38	34
Nordiazepam	271.1	139.9	32	30
7-aminonitrazepam	252.1	121.0	40	26
7-aminoclonazepam	286.1	120.9	38	28



The Ostro 96-well plate was used to remove residual PLs prior to LC/MS/MS analysis. Utilizing the generic, simple protocol provided by the manufacturer (Figure 2), a group of analogous compounds were extracted. The resulting analysis demonstrated an average recovery of 86% for the group of structural analogs in plasma (Figure 4), which is equal to or better than traditional PPT. To assess PL removal, eight individual PLs were summed. Results showed that Ostro plates remove >99% of the 8 PLs relative to traditional PPT (Figure 5). In addition, the MRM transition, 184>184, was monitored to visually demonstrate the significant decrease in residual PLs using Ostro plates compared to traditional PPT (Figure 6). One of the primary reasons to eliminate PLs is to improve method robustness. Overnight runs of both Ostro plates and traditional PPT samples were carried out using the generic gradient and PLs were monitored continuously. Figure 7 shows the LC/MS trace for a representative PL at the beginning and end of the runs. When Ostro plates are used, the PL levels are negligible and no build-up occurs. When PPT is used, a significant amount of PLs are present and accumulate throughout the run. The result of this undesirable build-up is a continuous downward trend in area counts throughout the duration of the run (Figure 8). This in turn results in high signal variability relative to samples prepared with Ostro plates, 33% using PPT vs. 9% for Ostro samples. In addition, area counts decrease by 57% from the first injection to the final injection when PPT is used. In discovery bioanalysis, high throughput is of utmost importance. If one tries to increase throughput by shortening gradient time, the impact of residual PLs is further magnified. To demonstrate the negative effect PLs have on analytical throughput, the gradient time was decreased by half. Flow rate was increased and organic content was ramped from 50 to 98% in 0.5 minutes. 200 Ostro samples and 200 PPT samples were injected using the shorter gradient. The MRM transition 184>184 was monitored to reveal overall PL build-up and elution in the shortened gradient window (Figure 9). Using the 2 minute gradient, PLs elute within 0.2 minutes of the end of the gradient. Under the truncated gradient conditions PLs continue to elute for more than 1 minute after the end of the gradient and well into the re-equilibration phase and beginning of the next




Figure 4. Average analyte recovery data for 26 structural analogs and metabolites in plasma using an Ostro 96-well plate. Overall average recovery was 86%.




Figure 5. Comparison of phospholipids remaining after traditional PPT and PPT using Ostro 96-well plate.

injection. These resultant chromatograms demonstrate the inability to shorten gradient time with PPT due to PLs which continue to elute significantly after the gradient ends at 0.5 minutes. Overall, the Ostro plate allows for increased method robustness and reduced variability. Additionally both improved instrument uptime and the ability to significantly shorten run times are realized through the elimination of PLs, all of which are highly desirable in a discovery setting. Calibration curves from 1-500 ng/mL for each of the 26 structural analogs had a resulting average  $r^2$  value of 0.925, sufficient for discovery screens.




Figure 6. Representative chromatograms of MRM transition 184>184 to demonstrate total remaining residual phospholipids from traditional PPT and Ostro eluates.




Figure 7. Representative chromatograms of MRM transition 758>184 to demonstrate build-up of an individual PL over subsequent injections using traditional PPT. A gradient from 2-98% B in 2 minutes was used followed by a 0.5 minute hold at high organic.






Figure 8. Variability of area counts over subsequent injections using Ostro 96-well plate vs. traditional PPT.




Figure 9. Representative chromatograms of the MRM transition 184 > 184 to demonstrate the overall build-up of PLs over subsequent injections using traditional PPT vs. Ostro 96-well plate. A fast gradient from 50-98% in 0.5 minutes was used. When PPT is used, PLs continue to elute well after the end of the gradient and when the Ostro plate is used, the short gradient can readily be implemented without concern over PL build-up.

---

## CONCLUSIONS

- Simple, universal sample preparation protocol
- Direct injection of the eluate streamlines workflow
- Removes >99% of phospholipids from plasma resulting in improved instrument uptime and more robust methods
- High recovery for analogs and metabolites with no method development
- Reduces sample variability
- Facilitates use of shorter runtimes, improving throughput

## REFERENCES

- <sup>1</sup> Rainville, P.D., Mather J, Waters Application Note; 2008, 720002569EN

## Sales Offices

**Austria and European Export**  
(Central South Eastern Europe, CIS  
and Middle East) 43 1 877 18 07

**Australia** 61 2 9933 1777

**Belgium** 32 2 726 1000

**Brazil** 55 11 4134 3788

**Canada** 1 800 252 4752

**China** 86 10 5293 6688

**Czech Republic** 420 2 617 11384

**Denmark** 45 46 59 8080

**Finland** 358 9 5659 6288

**France** 33 1 30 48 72 00

**Germany** 49 6196 400 600

**Hong Kong** 852 2964 1800

**Hungary** 36 1 350 5086

**India** 91 80 2837 1900

**Ireland** 353 1 448 1500

**Italy** 39 02 265 0983

**Japan** 81 3 3471 7191

**Korea** 82 2 6300 4800

**Mexico** 52 55 52 00 1860

**The Netherlands** 31 76 508 7200

**Norway** 47 6 384 6050

**Poland** 48 22 833 4400

**Puerto Rico** 1 787 747 8445

**Russia/CIS** 7 495 727 4490/  
290 9737

**Singapore** 65 6593 7100

**Spain** 34 93 600 9300

**Sweden** 46 8 555 115 00

**Switzerland** 41 56 676 7000

**Taiwan** 886 2 2501 9928

**United Kingdom** 44 208 238  
6100

**All other countries:**

**Waters Corporation U.S.A.**

1 508 478 2000

1 800 252 4752

[www.waters.com](http://www.waters.com)

# Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

© 2011 Waters Corporation. Waters, The Science of What's Possible, ACQUITY UltraPerformance LC, ACQUITY UPLC, Xevo, IntelliStart, and Ostro are trademarks of Waters Corporation.

©2011 Waters Corporation. Produced in the U.S.A.  
July 2011 720004046EN KK-PDF

**Waters Corporation**  
34 Maple Street  
Milford, MA 01757 U.S.A.  
T: 1 508 478 2000  
F: 1 508 872 1990  
[www.waters.com](http://www.waters.com)

