

Oasis PRIME HLB Cartridge for Cleanup of Infant Formula Extracts Prior to UPLC-MS/MS Multiresidue Veterinary Drugs Analysis

Michael S Young and Kim Van Tran Waters Corporation, Milford, MA, USA

APPLICATION BENEFITS

- Efficient, timesaving multiclass/ multiresidue methodology
- Simple, rapid, and effective sample cleanup suitable for a diverse range of analytes
- Fast, sensitive UPLC-MS/MS analysis

WATERS SOLUTIONS

ACQUITY UPLC® I-Class

Xevo® TQ-S Mass Spectrometer

Oasis[®] PRIME HLB Cartridge for SPE Clean-up

KEYWORDS

UPLC-MS/MS, Oasis PRIME HLB Cartridges, veterinary drugs, infant formula

OVERVIEW

In order to insure public health and safety, reliable analytical methods are necessary to determine veterinary drug residue levels in foods. Of particular importance is such residue analysis in foods for infants. The compounds of interest range from highly polar water-soluble compounds to very non-polar fat-soluble compounds. In order to maximize throughput and minimize costs it is desirable to determine the widest possible range of veterinary drug residues in such samples with a single analytical method. Powdered infant formula typically contains significant amounts of proteins, fats, and lecithin (phospholipids). These components can be detrimental to good instrumental performance and should be reduced or eliminated prior to LC-MS analysis.

INTRODUCTION

Infant formula powder contains significant amounts of protein, about 20% fat and 1-3 % phospholipids. During the sample pre-treatment, the protein is removed from the extract by precipitation and centrifugation. However, significant amounts of fat and phospholipids are co-extracted along with the target veterinary drugs. The presence of these co-extracted substances can lead to interference in the LC-MS analysis, contamination of the analytical column and other components of the UPLC® System, and contamination of the mass spectrometer itself. Fats have traditionally been removed from tissue extracts using cumbersome hexane defatting steps or by the use of reversed-phase sorbents such as C₁₈-silica in pass-through or dispersive cleanup. Although these techniques may be effective for fat removal, neither of these procedures removes phosholipids. In a prior study, sample preparation, cleanup, and analysis protocols were developed for tandem LC-MS determination of a wide variety of veterinary drug residues in seafood tissue samples. This cleanup protocol was effective for removal of both fats and phospholipids. In this study similar extraction and cleanup protocols were applied to the analysis of infant formula powder. Representative compounds were chosen from major classes of veterinary drugs including tetracyclines, fluoroquinolones, sulfonamides, macrolides, beta-lactams, NSAIDS, steroids, and beta-andrenergics. These compounds were spiked into the infant formula samples prior to extraction and cleanup.

EXPERIMENTAL

UPLC conditions

LC system:			ACQUITY UPLC I-Class				
			ACQUITY UPLC CSH [™] C ₁₈ , 1.7μm, 100 mm x 2.1 mm ID				
Mobile	phase A:		0.19	6 formi	c in water		
Mobile phase B:			0.1% formic acid in acetonitrile				
Injectio	n volume:		5 µl	_			
Injectio	n mode:		partial loop injection				
Column temp.:			30 °C				
Weak needle wash:			10:90				
			ace	tonitril	e:water		
			(60	0 µL)			
Strong needle wash:			50:30:40				
			water:acetonitrile:				
			IPA (200 µL)				
Seal wash:		10:90					
			ace	tonitril	e: water		
<u>Time</u> (min)	<u>Flow</u> (mL/min)	0/	<u>%A</u>	<u>%B</u>	<u>Curve</u>		
Initial	0.4	;	85	15	Initial		
2.5	0.4	(60	40	6		
3.9	0.4		5	95	6		
4.9	0.4		5	95	6		
5.0	0.4	1	85	15	6		

MS conditions

0.4

85

15

6

7.0

Mass spectrometer:	Xevo TQ-S
Positive ion	
electrospray	
source temp.:	150 °C
Desolvation temp.:	500 °C
Desolvation gas	
flow:	1000 L/Hr
Cone gas flow:	30 L/Hr
Collision gas flow:	0.15 mL/Min
Data management:	MassLynx® v4.1

Sample preparation

1. Initial Extraction/Precipitation:

Place a 0.5 g sample of infant formula into a 50 mL centrifuge tube. For standards or QC samples spike with appropriate amounts of desired analytes. Add 3 mL extraction solvent (0.2% formic acid in 70:30 acetonitrile/water). Vortex for 30 seconds and place on mechanical shaker for 30 minutes. Centrifuge at 3220 rcf for 5 minutes.

Note: The extraction/precipitation step gives good recovery of most compounds of interest but also extracts significant amounts of fat and phospholipid.

2. SPE Cleanup:

Mount an Oasis PRIME HLB Cartridge (3 cc, 60 mg) on a pre-cleaned vacuum manifold. Cartridge conditioning is NOT required and is not performed. The vacuum is set to 1–2 psi. Approximately 0.5 mL of the supernatant is passed-through the Oasis PRIME Cartridge and collected. A 0.3 mL aliquot of the pass-thru cleanup sample is taken and diluted three-fold with aqueous 10 mM ammonium formate buffer (pH 4.5) prior to UPLC-MS/MS analysis.

[APPLICATION NOTE]

Compounds	MRM	Cone (V)	Collision (eV)	Spike level (low, high) µg/kg	Calibration range µg/kg	Corr (R²)	RT
Carbadox	263.0>231.0 263.0>145.0	25 25	15 20	25, 100	12.5-200	0.9990	1.37
Ceftiofur	524.3>241.1 524.3>285.0	30 30	16 16	250, 1000	125-2000	0.9980	2.79
Chlortetracycline	479.3>444.2 479.3>462.2	30 30	21 18	25, 100	12.5–200	0.9977	1.53
Ciprofloxacin	332.1>288.1 332.1>231.1	30 30	18 40	25, 100	12.5–200	0.9929	0.84
Cortisol	363.2>121.0 363.2>91.03	42 30	52 22	50, 200	25-400	0.9961	2.93
Dexamethasone	393.2>373.2 393.2>355.3	30 30	10 15	25, 100	12.5-200	0.9917	3.39
Enrofloxacin	360.4>245.0 360.4>316.1	50 50	25 25	50, 200	25-400	0.9991	0.98
Erythromycin	734.4>158.1 734.4>576.5	30 30	32 20	2.5, 10	1.25-20	0.9978	2.16
Lincomycin	407.2>126.1 407.2>359.3	36 36	34 20	12.5, 50	6.25-100	0.9962	0.57
Lomefoxacin	352.1>265.1 352.1>308.1	31 31	22 16	50, 200	25-400	0.9991	0.90
Oxacillin	402.2>160.0 402.2>243.1	30 30	12 15	25, 100	12.5-200	0.9990	3.75
Oxytetracyline	461.2>426.2 461.2>408.11	30 30	21 13	25, 100	12.5-200	0.9971	0.96
Pennicillin	335.16>160.1 335.15>176.1	20 20	30 30	12.5, 50	6.25-100	0.9956	3.41
Phenylbutazone	309.4>160.0 309.4>103.9	37 37	20 20	25, 100	12.5-200	0.9962	4.22
Ractopamine	302.2.164.1 302.2>107.0	30 30	15 27	75, 300	37.5-600	0.9990	0.90
Salbutamol	240.2>148.1 240.2>222.1	30 30	20 12	25, 100	12.5-200	0.9941	0.55
Sulfamerazine	265>92.0 265>156.0	30 30	28 15	25, 100	12.5-200	0.9998	1.50
Sulfamethazine	279.1>186.0 279.1>92.0	30 30	16 28	25, 100	12.5-200	0.9996	1.67
Sulfanilamide	156>92.0 156>65.0	30 30	15 25	25, 100	12.5-200	0.9964	0.86
Tetracycline	445.3>154.0 445.3.410.2	30 30	26 21	25, 100	12.5-200	0.9964	1.05
Tylosin	916.5>174.1 916.5>101.1	57 57	40 45	5, 20	2.5-40	0.9920	2.38

Table 1. Matrix matched calibration data, MRM transitions (primary transition first), instrument parameters, and retention times (RT) used for this study.

RESULTS

Table 2 shows the recovery data obtained from replicate analysis of spiked tissue samples. Matrix effects averaged about 40% for infant formula. The chromatograms shown in Figure 1 show the effectiveness of the Oasis PRIME HLB Cartridge for removal of \geq 95% of phopholipids from the infant formula extracts.

	Low le	vel	High Level		
Compounds	Recovery (%)	RSD (%)	Recovery (%)	RSD (%)	
		n=6		n=6	
Carbadox	110	5	115	3	
Ceftiofur	84	4	71	10	
Chlortetracycline	47	9	46	9	
Ciprofloxacin	107	8	102	6	
Cortisol	111	8	117	5	
Dexamethasone	113	14	121	5	
Enrofloxacin	113	4	110	5	
Erythromycin	126	6	125	6	
Lincomycin	50	7	52	5	
Lomefloxacin	120	2	111	2	
Oxacillin	117	6	114	4	
Oxytetracycline	29	13	27	13	
Penicillin	120	10	116	7	
Phenylbutazone	105	13	94	6	
Ractopamine	115	2	117	5	
Salbutamol	82	7	83	4	
Sulfamerazine	130	4	126	4	
Sulfamethazine	128	2	129	3	
Sulfanilamide	105	12	120	5	
Tetracyline	41	13	40	17	
Tylosin	110	13	116	7	
Carbadox	110	5	115	3	
Ceftiofur	84	4	71	10	

Table 2. Recovery data obtained from replicate analysis of spiked infant formula samples (n = 6).

4

DISCUSSION

The procedure utilized in this study was developed from methods presented previously.^{1,2} The overall method recoveries are generally above 70% but lower recovery was observed for some of the more polar compound classes, such as tetracyclines. However, the Oasis PRiME HLB Cartridge cleanup contributes very little to any method recovery losses. As shown in Figure 2, the measured recovery for the SPE cleanup step is better than 80%.




Figure 1. LC-MS/MS chromatograms showing effective removal of \geq 95% of phospholipids from shrimp extract.




Figure 2. Recovery of veterinary compounds from prepared extracts subjected to Oasis PRiME HLB pass-through clean-up.



- A simple and effective extraction/protein precipitation procedure was applied to the analysis of infant formula
- A simple one-step pass-thru cleanup protocol using Oasis PRIME HLB Cartridges was employed to remove greater than 90% of fats and phospholipids from the initial extracts
- The sample preparation methodology produced an extract that was free of particulates and required no subsequent filtration prior to LC-MS analysis
- High and consistent recoveries were observed for a wide range of veterinary drugs using the simple one-step pass-through clean-up protocol with Oasis PRIME HLB Cartridges

References

- M. Young and K. Tran, Oasis PRIME HLB Cartridge for Effective Cleanup of Meat Extracts Prior to Multi-Residue Veterinary Drug UPLC-MS Analysis, Waters Technology Brief <u>720005411EN</u>, 2015.
- M. Young and K. Tran, Rapid, Simple, and Effective Cleanup of Seafood Extracts Prior to UPLC-MS/MS Multiresidue Veterinary Drugs Analysis, Waters Application Note <u>720005488EN</u>, 2015.



Waters, The Science of What's Possible, UPLC, Oasis, ACQUITY UPLC, Xevo, and MassLynx are registered trademarks of Waters Corporation. CSH is a trademark of Waters Corporation. All other trademarks are the property of their respective owners.

©2016 Waters Corporation. Produced in the U.S.A. June 2016 720005730EN AG-PDF

Waters Corporation

34 Maple Street Milford, MA 01757 U.S.A. T: 1 508 478 2000 F: 1 508 872 1990 www.waters.com