Essentials in biorese

930654

nstunateran., Chanisties, and Copositios (c) Biosecentions

Poster Presentation

Prep '91 - May, 1991

Simplified Methods Development and Scale-Up for Organic Synthesis Using Segmented Column TechnologyTM

K. O'Connor and C. Phoebe, Jr. Millipore Corporation, Waters Chromatography Division 34 Maple Street, Milford, MA. 01757 U.S.A.

Waters Chromatography Division Millipore Corporation 34 Maple St. Milford, MA 01757 (508) 478-2000



Introduction

In order to increase the yield of a final product, many synthetic organic chemists need to purify reaction intermediates. The large scale purification of these synthetic intermediates is traditionally performed using flash chromatography. This methodology offers low resolution and can use excessive amounts of solvent. Presented here is an alternative methodology for the separation of reaction intermediates by HPLC using Segmented Column Technology (SCT). SCT is a novel approach for preparative chromatography that permits column length and/or diameter to be varied in order to achieve increased resolution and/or capacity. Up to three 8mm, 25mm or 40mm I.D. cartridge columns or segments in 100mm lengths can be connected to optimize a separation and can conserve solvent. To demonstrate this simple and direct approach to the preparative isolation of reaction intermediates, the cis- and trans- products of a well characterized Wadsworth-Emmons reaction employed in the synthesis of many pharmaceutically interesting retenoids were purified using normal phase chromatography and SCT (Figure 1).

Figure 1: *cis-* and *trans-* Isomers formed by Wadsworth-Emmons Reaction



Analytical Separation

A separation was first developed on a high resolution Nova-Pak® Silica, 4μ m, 3.9mm x 150mm analytical column, so that fractions ultimately collected from a preparative separation could be analyzed (Figure 2). The desired cis-isomer (peak 1) was well resolved from the trans- isomer (peak 2) and other reaction by-products (peaks 3-5). The analytical separation was then transferred to a preparative cartridge column. This transfer required the determination of which preparative silica packings should be used and the length and loading of the cartridge(s) required for this separation. In order to determine these parameters, 8mm I.D. scaling cartridges were utilized. Once the optimum conditions were determined, sample load and solvent consumption were decided by using several equations. The 8mm cartridge conditions could then be scaled to either 25mm or 40mm I.D. cartridges.

Figure 2: Analytical Separation of Crude Material



.

.

Preparative Silica Evaluation

Four different preparative silica packings, (6µm, 10µm, 15-20µm and 37-55µm), packed in 8mm x 100mm cartridges were evaluated to determine their ability to resolve the *cis*- and *trans*- isomers (Figure 3). Both the 6µm and the 10µm silica packings provided good resolution, whereas the 15-20µm and 37-55µm silica packings poorly resolved the pair for this length of cartridge column. The 10µm particle size silica was chosen to determine further scaleup parameters because it resolved the *cis/trans*-isomers and would provide a higher sample capacity than the 6µm silica.

Figure 3: Comparison of Four Silica Packings



Figure 4: Comparison of 100mm, 200mm and 300mm Cartridge Column Lengths



Loading Study

A loading study was performed on the 8mm x 200mm cartridge column packed with 10µm silica to determine the maximum amount of crude sample the column could tolerate without compromising resolution of the cis- and trans- isomers. Once established, the flow rate and the maximum load on the scaling column could now be used to determine the conditions to be employed for any larger volume column packed with the same material. The sample load and flow rate were increased proportionally to the volume of the packed bed (Table 1). This results in identical chromatographic performance (retention times, peak widths and resolution) at the preparative scale as on the scaling cartridge column. The maximum sample load, on the scaling cartridge column, was 2.12 mg (Figure 5). Using equations 1 & 2, the flow rate and load were calculated for a 200mm length (two cartridge column configuration) with a 40mm diameter cartridge column.

Figure 5: Separation of Crude Material on a Scaling Column



Table 1: Flow Rate and Sample Load IncreaseProportionally to Column Volume



Figure 6: Preparative Separation of *cis* and *trans*-isomers



Large Scale Purification

As predicted from Table 2, 52.88mg (3.75mls) could be loaded on the 40mm x 200mm cartridge column, resulting in the same chromatography (Figure 6). Fractions across the peaks of interest were collected, and aliquots were injected on the Nova-Pak Silica, 4µm, 3.9mm x 150mm analytical column to assess purity. Fractions 1A, 1B and 1C proved to be the *cis*-isomer while the remaining three were the *trans*-isomer (Figure 7).

Table 2: Choice of Cartridge Columns Which Could beUsed for the Separation of cis- and trans- isomers

Number of Cartridges	Column Volume (mls)	Flow Rate (ml/min)	Solvent Consumption (mls)	Load (mg)	Time of Analysis
1(8 x 100mm)	5.0	2.8	22.4	1.06	8.0
2 (8 × 100mm)	10.0	2.8	44.8	2.12	16.0
3 (8 × 100mm)	15.0	2.8	67.2	3.18	24.0
1 (25 x 100mm)	50.0	28.0	224.0	10.60	8.0
2(25 x 100mm)	100.0	28.0	448.0	21.20	16.0
3(25×100mm)	150.0	28.0	672.0	31.80	24.0
1 (40 x 100mm)	125.0	70.0	560.0	26.50	8.0
2 (40 x 100mm)	250.0	70.0	1120.0	52.88	16.0
3 (40 × 100mm)	375.0	70.0	1680.0	79.50	24.0

Shaded areas represent cartridge columns used in Figures 5&6. The sample load does not reflect the absolute capacity of the various cartridge columns.

Figure 7: Analysis of Fractions from Preparative Separation



Alternative Methodology

The pathway presented here represents only one way in which to perform this separation. The speed and capacity were optimized to offer a compromise between the two. Many other options exist when using Segmented Column Technology (SCT). For example, this separation could have been tailored to optimize speed or capacity. Using Equations 1 and 2 (Table 1) other options in column configuration can be employed (Table 2).

Determining Column Length: Optimizing the Separation to Suit Your Application

Injections were made on 100mm, 200mm and 300mm column lengths, using 8mm cartridge columns packed with 10µm silica to determine which would offer the largest sample capacity, resolution and fastest run time (Figure 4). The 100mm length offered a fast separation with low peak volumes but is limited in capacity. The 300mm length offered the greatest resolution and increased capacity but its long run time would use larger amounts of solvent and increase the peak volumes. The 200mm length was chosen for this application because it provided the best compromise between resolution of peaks of interest, separation time and solvent consumption.

Conclusion

The use of SCT offers the synthetic organic chemist a means to purify large amounts of material without the column defining the purification conditions. This separation was assessed by considering resolution, capacity, time and solvent consumption. The column length and packing material size were tailored to the separation and a capacity study was performed at the analytical scale. Having determined the optimum column length and capacity, the separation was scaled to a 40mm cartridge column with the same number of segments as in the scaling separation. Fractions collected from the preparative separation were analyzed by HPLC, confirming that pure products were obtained.

Acknowledgements

Sample courtesy of Hoffman LaRoche