

930786

Rx 029 12/90

# R Prescription for success

## **REDUCED ANALYSIS TIME FOR CYCLIC PEPTIDE**



#### CONDITIONS ON WATERS QUANTATM 4000 MODE: FZCE BUFFER: A: $100 \text{ mM H,PO}_4$ pH = 2B: $50 \text{ mM Na,HPO}_4$

#### MODIFIER: CAPILLARY: VOLTAGE: DETECTOR: INJECTION:

#### pH = 2 B: 50 mM Na,HPO pH = 9 with NaOH None 75 $\mu$ m x 60 cm + 17 KV 214 nm

10 sec Hystrostatic

### **PEAK IDENTIFICATION:**

- 1. Cyclic Peptide
- 2. 1.0% Impurity

#### REFERENCE: John Van Antwerp and Larry Mugavero, Application Chemists, Morristown, N.J. Laboratory

Author: Peter Rahn

SAMPLE MATRIX:

Millipore Corporation Waters Chromatography Division 34 Maple Street Milford, MA 01757 (50

## **INTERESTING FACTS**

1. In both CE separations, the Quanta 4000 easily detects impurities at very low concentrations. Only the large linear dynamic range of the UV detector, and the extremely quiet baseline of the Quanta 4000 provides good quantitation of these small impurities. Under both CE separation conditions, the impurity level was quantitated at 1 %.

2. The original separation was performed at pH = 2 where minimal electroosmotic flow occurs. This separation required 22 minutes to complete.

3. The same sample was also separated using a pH = 9 electrolyte. At this pH the electroosmotic flow was a powerful driving force for the electrolyte in the capillary and the analysis time was reduced to 8 minutes (a 60 % reduction).

4. Unlike HPLC, repeated injections on the CE can be made every 10 minutes since the CE separations do not require equilibration after a gradient separation.