


930787

**R** Prescription for success

Rx 027 10/90

### CAPILLARY ELECTROPHORESIS SEPARATION OF PENICILLIN RELATED ANTIBIOTICS



#### CONDITIONS ON WATERS QUANTA™ 4000

MODE: MECC  
BUFFER: 0.02M  $\text{NaH}_2\text{PO}_4$   
0.02M  $\text{Na}_2\text{B}_4\text{O}_7$   
pH = 9  
MODIFIER: 50 mM SDS  
CAPILLARY: 75  $\mu\text{m}$  i.d. x 60 cm  
VOLTAGE: + 18 KV  
DETECTOR: 214 nm  
INJECTION: 5 sec Hydrostatic

#### PEAK IDENTIFICATION:

1. Amoxicillin
2. Ampicillin
3. 6-Amino penicillanic acid
4. Oxacillin
5. Cloxacillin
6. Ticarcillin
7. Nafcillin
8. Dicloxacillin

SAMPLE MATRIX: Standard Solution  
0.12 mg/ml in  $\text{H}_2\text{O}$

REFERENCE:  
Michael Swartz, Pharmaceutical Laboratory, Waters  
Chromatography Division

Author: Peter Rahn

## INTERESTING FACTS

1. Capillary electrophoresis provides the efficiency and resolution to quickly separate related compounds. In this example peak 3, 6-amino penicillanic acid, is the building block for the seven antibiotics shown in this separation.
2. The Quanta 4000 provides excellent sensitivity for these antibiotics. Previous work has shown that the calibration curves are linear on the Quanta 4000 and therefore it is feasible to quantitate any impurities present.
3. Another important tool provided by CE is the identification and quantitation of related materials due to cross contamination. The selectivity of CE would provide a rapid tool to positively identify cross contamination.