

**Poster Presentation**

**The Protein Society  
Seventh Symposium  
July, 1993**

**Poster #301-T**

**Improved Cysteine Analysis Using  
Disulfide Interchange**

**K.M. De Antonis\*, S.A. Cohen\*\*,  
and P.R. Brown\***

**\* University of Rhode Island Department of Chemistry,  
Kingstown, RI 02813 U.S.A.**

**\*\*Millipore Corporation, Waters Chromatography Division,  
34 Maple Street, Milford, MA 01757 U.S.A.**

**MILLIPORE**

Waters Chromatography Division  
34 Maple Street  
Milford, MA 01757  
508 478-2000

## Abstract

Amino acid analysis has become a simple and effective tool for quantitation of protein and peptide amino acids. However, under standard hydrolysis conditions unprotected cysteine is highly susceptible to thermal degradation and is characteristically difficult to quantify. A recently reported method uses disulfide interchange reagents for cysteine protection (1). Cysteine reacts with each dithiocompound to form a thermally stable mixed disulfide. Recently, we have used this method in conjunction with the novel derivatizing agent, 6-aminoquinolyl-N-hydroxysuccinimidyl carbamate, for sensitive amino acid analysis of proteins (2, 3). This poster will compare the analyses of several protein samples hydrolyzed in the presence of three dithio-compounds dithiodiglycolic acid (DTDGA), dithiodipropionic acid (DTDPA), and dithiobutyric acid (DTDBA). As expected, the three cysteine disulfide adducts were detected. When necessary, modifications were made in the normal hydrolysate chromatographic conditions to improve resolution of the cysteine analog. In addition to the cysteine-derived peak, a smaller previously unidentified peak was detected, and has been tentatively identified for each of the three.

Protein fractions ( $\beta$  lactoglobulin and cytochrome C) were collected on a membrane fraction collector following separation by capillary electrophoresis. These proteins were then hydrolyzed on the membrane in the presence of DTDPA and analyzed using AQC derivatization.

- [1] Barkholt, V., and Jenson, A.E., *Anal. Biochem.* 177 (1989) 318.
- [2] Stridom, D.J., and Cohen, S.A., (1993) in *Techniques in Biochemistry IV*, (R.H. Angeletti, ed.), Academic Press, San Diego.
- [3] Cohen, S.A., and Michaud, D.P., *Anal. Biochem.*, 211 (1993) 279.

## Instrumentation

### Hydrolysis Instrumentation

Waters Pico-Tag Workstation  
Precision Scientific Model 16 hydrolysis oven (Chicago, Illinois)

### Chromatographic Instrumentation

Waters 715 autoinjector  
Waters 625 pump  
3.9 x 150 mm AccQ-Tag™ C18 reverse phase column  
Waters Temperature Control Module  
Column Temp. 37°C for DTDGA and DTDPA, 39°C for DTDBA  
Waters 470 Scanning Fluorescence Detector  
 $\lambda$  Excitation 250 nm  
 $\lambda$  Emission 395 nm  
Waters Expert Ease™ software

### CE Membrane Fraction Collection

Waters Quanta 4000E™  
Membrane Fraction Collector


### HPLC Analysis

Amino acid analysis was achieved using a ternary gradient.


Eluent A	140 mM NaOAc, 17 mM TEA, 3 μM EDTA, pH 5.00
Eluent B	for DTDGA or 5.05 for DTDPA and DTDBA
Eluent C	Acetonitrile
	Water
Flow Rate	1 mL/min

Different gradients were used for each of the three reagents in order to resolve each of the disulfide interchange products from the other amino acids.

## Derivatization Chemistry



## Proposed Disulfide Interchange Chemistry\*



X = 1 for DTDGA, 2 for DTDPA, and 3 for DTDBA

## Sample Preparation

### Hydrolysis

1. Deliver the protein containing section of the membrane (washed first with methanol and water) or an aliquot of amino acid or protein sample to a pyrolyzed 6 x 50 mm tube.
2. Put all tubes in a hydrolysis vial and dry under vacuum.
3. Add an aliquot of the disulfide interchange reagent to each tube.
4. Vortex each tube and dry under vacuum.
5. Add constant boiling HCl and phenol to bottom of hydrolysis vial.
6. Evacuate and purge with nitrogen several times.
7. Evacuate and seal under vacuum.
8. Hydrolyze at 110°C for 22 hours.

### Derivatization

#### Standards

1. Mix 1 volume of derivatizing reagent to 4 volumes of buffered (pH 8.8) amino acid standard.
2. Vortex.
3. Heat for 10 minutes at 50°C.
4. Inject.

#### Hydrolyzed Proteins


1. Dry contents of hydrolysis vial.
2. Dissolve hydrolyzed sample in dilute HCl.
3. Buffer solution to pH 8.8.
4. Mix 1 volume of derivatizing reagent to 4 volumes of sample.
5. Vortex.
6. Heat for 10 minutes at 50°C.
7. Inject.

#### Hydrolyzed Membrane Bound Proteins


1. Extract hydrolyzed amino acids from the membrane using 50:50 Acetonitrile:4% SDS.
2. Remove membrane.
3. Follow steps 1-7 for hydrolyzed proteins.

\* Hoogendoorn, J.G., and Campbell, C.M., *Anal. Biochem.*, 210 (1992) 146.

### DTDGA Analysis




### DTDBA Analysis



1. AMQ	10. Ala	Gradient:				
2. Asp	11. Pro	Time	%A	%B	%C	Curve
3. Ser	12. Tyr	Inrt.	100	0	0	L
4. Glu	13. Val	0.5	99	1	0	L
5. Gly	14. Met	18	95	5	0	L
6. His	15. Lys	19	91	9	0	L
7. NH3	16. Ile	24	90	10	0	L
8. Arg	17. Leu	28	83	17	0	L
9. Thr	18. Phe	30	83	17	0	L
• Reagent byproduct		34	0	60	40	S
• Hydrolysis byproduct						

Gradient:					
Time	%A	%B	%C	Curve	
Inrt.	100	0	0	L	
0.5	99	1	0	L	
18	95	5	0	L	
19	91	9	0	L	
22	86	14	0	L	
26	86	14	0	L	
34	83	17	0	L	
35	0	60	40	S	

### DTDPA Analysis



Gradient:					
Time	%A	%B	%C	Curve	
Inrt.	100	0	0	L	
0.5	99	1	0	L	
18	95	5	0	L	
19	91	9	0	L	
29.5	83	17	0	L	
33	0	60	40	S	

M	Comp.	B Lactoglobulin		Cytchrome C		Avg. Error 18 hyd.
		Comp.	Cytochrome C	Comp.	Cytochrome C	
Asp	7	15.5	15.6	8	8	7.8
Ser	6.6	6.6	6.5	9	9	7.7
Glu	24.9	24.8	24.6	12	12	12.0
Gly	3.1	3.0	3.1	11.7	11.7	11.9
His	2.0	2.0	2.0	10.8	10.8	11.1
Arg	3.9	3.9	3.1	3.0	3.0	3.1
Thr	7.8	7.8	8.0	10.4	10.4	10.3
Ala	13.3	13.3	13.4	5.8	5.7	5.9
Pro	7.9	7.9	7.9	3.9	3.9	4.0
Cys	5.2	5.2	5.3	2.3	2.3	2.4
Tyr	4.6	4.4	4.4	2.2	2.2	2.3
Val	9.7	9.7	9.6	3.2	3.2	3.3
Met	4.3	4.4	4.4	1.9	1.9	2.0
Ile	14.9	14.9	14.9	10.0	10.0	10.3
Leu	9.4	9.3	9.2	17.8	17.8	18.0
Phe	21.7	21.6	21.4	16.3	16.3	16.8
	4.0	4.1	4.0	1.5	1.5	1.5
	3.7	3.8	4.0	5.0	5.0	5.0
	4.0	4.1	4.0	7.1	7.1	7.1
	1.4	1.4	1.4	3.9	3.9	3.9
				4.0	4.0	4.0

### DTDGA - AAA

- Good average error values for both proteins.
- Cys value high due to slight coelution with Ala.
- Highest yield of intact Met for the three disulfide interchange reagents.
- Run not measured

## DTDPA - Membrane Fraction Collector

AA	Comp.	B Lactoglobulin				Cytochrome C			
		B Lactoglobulin		Cytochrome C		Comp.	Cytochrome C		
Asp	16	15.8	16.5	17.0	17.0	8	8.8	8.8	8.5
Ser	7	9.0	10.0	9.2	8.7	0	nm	nm	nm
Glu	25	22.8	23.8	24.1	24.9	12	11.9	11.6	11.9
Gly	3	4.7	6.9	6.4	4.4	12	12.9	12.5	12.8
His	2	1.3	1.7	1.5	1.5	3	2.4	2.2	2.4
Arg	3	4.2	6.0	4.1	4.0	2	1.7	3.0	3.2
Thr	8	8.4	7.8	7.9	7.9	10	9.5	9.4	9.5
Ala	14	15.0	14.3	14.6	15.1	6	6.4	6.1	6.3
Pro	8	7.1	7.2	7.4	7.6	4	3.9	3.7	4.1
Cys	5	11.5	7.3	10.3	7.5	2	4.4	3.4	3.7
Tyr	4	3.8	4.0	4.0	4.0	4	2.9	3.1	
Val	10	9.5	9.5	9.6	10.0	3	3.2	3.2	3.1
Met	4	4.0	3.7	3.6	3.9	2	1.9	1.7	1.8
Lys	15	12.8	12.4	12.2	12.9	19	16.4	14.6	15.0
Ile	10	8.7	8.8	8.8	9.2	6	6.0	5.7	5.7
Leu	22	20.7	20.7	21.0	21.8	6	6.1	6.2	6.0
Phe	4	3.9	4.0	4.1	4.1	4	4.1	3.8	3.8
Avg. Err <sup>a</sup>		22 (15)	24 (22)	23 (18)	14 (11)		15 (8)	17 (13)	16 (12)
μg Hyd.		0.19	0.31	0.23	0.38		0.63	0.68	0.71
Recovery		32%	52%	40%	64%		76%	83%	86%

- Average error of all amino acids except Cys in parenthesis.
- A hydrolysis or membrane related impurity coelutes with Cys-P and causes high Cys values.
- Good overall error for low level Cytochrome C without Cys-P.

nm not measured

## DTDPA - AAA

AA	Comp.	B Lactoglobulin				Cytochrome C			
		B Lactoglobulin		Cytochrome C		Comp.	Cytochrome C		
Asp	16	17.9	18.3	18.4	8	9.1	9.1	8.8	
Ser	7	7.0	7.0	7.1	0	nm	nm	nm	
Glu	25	25.7	26.1	26.0	12	12.6	12.0	12.2	
Gly	3	3.4	3.2	2.9	12	11.6	11.5	11.6	
His	2	2.0	1.9	1.9	3	3.0	3.0	2.0	
Arg	3	2.8	2.7	3.6*	2	1.9	2.3	2.3	
Thr	8	7.6	7.7	7.8	10	9.5	9.2	9.2	
Ala	14	13.6	14.0	14.0	6	6.1	5.8	5.7	
Pro	8	7.9	8.0	7.9	4	4.0	3.9	3.8	
Cys	5	5.0	5.1	5.0	2	2.0	2.0	2.0	
Tyr	4	4.0	4.0	4.0	4	3.7	3.9	4.0	
Val	10	9.4	9.2	9.6	3	2.9	2.9	2.9	
Met	4	4.1	4.1	4.0	2	2.0	2.0	2.0	
Lys	15	15.0	15.1	15.0	19	18.8	18.2	18.7	
Ile	10	8.9	9.0	9.0	6	5.8	5.7	5.8	
Leu	22	21.5	22.1	21.8	6	6.1	5.9	6.1	
Phe	4	4.1	4.1	4.1	4	4.0	4.2	4.0	
Avg. Error		2.2	2.3	2.3		2.3	2.5	2.4	
μg Hyd.		4.3	4.1	4.0		3.5	4.6	3.8	

- Arg results complicated by coelution with a hydrolysis related impurity.
- Good average error values for both proteins.

nm not measured

## 4 CHLORINE ACID OXIDATION

M	Comp.	B Lactoglobulin				Cytochrome C			
		B Lactoglobulin		Cytochrome C		Comp.	Cytochrome C		
Asp	16	17.1	17.5	17.3	8	8.4	8.3	8.2	
Ser	7	5.9	5.9	5.9	0	nm	nm	nm	
Glu	25	25.3	26.0	25.9	12	12.1	12.1	12.0	
Gly	3	3.5	3.2	3.3	12	12.3	12.3	13.1	
His	2	2.0	2.0	2.0	3	3.0	3.0	3.0	
Arg	3	3.0	2.9	3.0	2	2.0	1.9	2.1	
Thr	8	7.4	7.4	7.2	10	9.0	8.9	8.7	
Ala	14	13.4	13.6	13.4	6	5.8	5.7	5.6	
Pro	8	8.1	8.1	8.1	4	3.9	4.0	3.9	
Cys	5	4.2	4.3	4.3	2	5.8	5.9	5.9	
Val	10	9.5	9.6	9.6	3	2.9	3.0	3.0	
Met SO <sub>2</sub>	4	5.6	5.7	5.8	2	2.8	2.8	2.7	
Lys	15	21.6	21.8	21.7	19	20.3	21.5	21.1	
Ile	10	9.2	9.3	9.2	6	5.7	5.8	5.7	
Leu	22	21.5	21.7	21.6	6	5.8	5.9	5.9	
Phe	4	4.2	4.2	4.1	4	4.5	4.3	4.1	
Avg. Error		10.8	10.7	11.0		9.7	9.4	9.2	
μg Hyd.		3.0	3.0	3.0		4.8	4.5	4.7	

nm not measured

Comparison with Disulfide Interchange Data (reported as an average, n=3):

	# Cys	Cys	Cys-G	Cys-P	Cys-B
B Lac	5	6.0	5.3	5.0	5.7
Cyt C	2	1.1	2.2	2.0	2.7

## DTDBA - AAA


AA	Comp.	B Lactoglobulin				Cytochrome C			
		B Lactoglobulin		Cytochrome C		Comp.	Cytochrome C		
Asp	16	16.9	16.8	16.8	8	7.7	7.7	7.6	
Ser	7	6.0	7.0	7.0	0	nm	nm	nm	
Glu	25	26.0	25.3	25.8	12	14.3	13.8	13.8	
Gly	3	3.0	3.1	3.1	12	11.8	11.7	11.6	
His	2	1.9	1.9	2.0	3	2.9	2.9	2.9	
Arg	3	3.0	3.0	3.0	2	2.0	2.0	2.0	
Thr	8	8.0	8.2	8.0	10	10.2	9.9	9.7	
Ala	14	14.0	14.1	14.0	6	5.8	5.8	5.8	
Pro	8	8.0	8.0	8.0	4	3.8	3.9	3.8	
Cys	5	6.0	5.5	5.6	2	2.7	2.8	2.7	
Tyr	4	4.1	4.1	4.0	4	4.2	4.2	4.1	
Val	10	9.8	9.6	9.6	3	3.1	3.1	3.1	
Met	4	4.6	4.5	4.4	2	2.3	2.4	2.3	
Lys	15	13.2	13.0	13.1	19	15.5	15.4	15.3	
Ile	10	8.9	8.8	8.8	6	5.7	5.7	5.7	
Leu	22	22.1	21.7	21.5	6	5.8	5.8	5.8	
Phe	4	4.0	4.0	4.0	4	3.8	4.0	3.9	
Avg. Error		4.9	4.4	4.2		7.87	8.06	7.82	
μg Hyd.		1.5	1.6	1.6		1.48	1.56	1.61	

\* Cys values are high due to poor resolution between Cys and Val.


\* Low Lys values due to disulfide interchange related impurity which coeluted with Lys in the standard.

nm not measured

## Methionine Adduct



Proposed Methionine Reaction:



## Conclusions

- Modified chromatographic gradients were developed which resolved the disulfide interchange products from the other amino acids.

- Analysis using AQC and disulfide interchange gave good overall error results and good results for cysteine for both  $\beta$  lactoglobulin and cytochrome C at the 1-5  $\mu$ g level. The overall error results were in the 2 - 8% range.

- The disulfide interchange results for cysteine are as good as the results generated by the performic acid oxidation technique.

- A disulfide interchange reagent byproduct present in each of the standard mixtures appears to be a methionine adduct. When methionine was hydrolyzed in the presence of each of the disulfide interchange reagents, a product was generated which was coeluted with the reagent byproduct.

- Based on the relative size of the peaks generated by the reagent byproducts in the methionine samples, reaction with methionine appears to be the slowest with DTDGA. In order to minimize methionine loss, DTDGA might be the best disulfide interchange reagent to use.

- Membrane fraction collection was attempted. An impurity from either the membrane or the disulfide interchange reaction was coeluted with the Cys-P thus complicating the result. The overall error ignoring the Cys-P result was good for the low level of protein (100 - 700 ng) extracted from membrane. Chromatographic resolution of the impurity from the Cys-P peak or the use of DTDGA will enhance this technique for analysis of cysteine from membrane bound proteins.