


NOTES

PHARMACEUTICAL

94-0421



PUBLISHED BY WATERS CHROMATOGRAPHY DIVISION, MILLIPORE CORPORATION

VOL. 1, NO. 1

SUMMER 1989

New PowerLine Based Automated System For Dissolution Control: Waters 610/600E PowerLine System for Dissolution Testing

Introduction

In 1987 Waters first introduced an automated dissolution/HPLC (High Performance Liquid Chromatograph) interface based on the 840 Data and Chromatography Control Station. This system improves the efficiency of the dissolution lab by transferring samples directly from the dissolution test station to the HPLC, analyzing those samples and documenting the results.

In March of this year, at the 40th Annual Pittsburgh Conference and Exposition, Waters announced the incorporation of dissolution acquisition software in the PowerLine™ 610 isocratic and the 600E gradient solvent delivery and system control modules. The 610 isocratic PowerLine system module can be upgraded to a

four solvent gradient system (functionally equivalent to the 600E) with the addition of a gradient valve station. These PowerLine modules are a new way to integrate the latest technologies of today's HPLC system components into a combined centralized controller, providing complete system communication and system parameter documentation. This new approach from Waters provides power beyond that of stand-alone HPLC components. Additionally, PowerLine HPLC systems can be configured with either an integrator, or a chromatography workstation using Waters Maxima™ or Expert Ease™ software.

The PowerLine Automated Dissolution/HPLC system consists of a 600E or 610 PowerLine module (Rev. 2.0 Software), a 484 Tunable UV/VIS Absorbance detector, a WISP™ 712D

IN THIS ISSUE	
New PowerLine Based Automated HPLC System for Dissolution Control . . .	1-3
Dissolution Testing Requiring High Sensitivity Analysis . . .	4-5
The Waters 845 Chromatography Data Station . . .	5-6
Waters Extends its Family of HPLC Auto Injectors . . .	7

sample processor, and a Hanson Transfer Controller (TC)* which transfers dissolution samples from the dissolution test apparatus directly into the active 712D sample processor.


How the PowerLine System for Dissolution Testing Works

The following operational description is applicable to either the 610 or the 600E PowerLine solvent delivery and system control modules.

Dissolution control is accessed from the Autoinjector Setup screen of the 600E controller. When enabled, the 600E coordinates the transfer of dissolution samples into the active 712D autosampler carousel while performing chromatography on previously collected dissolution samples or calibration standards.

*Hanson Research Corporation
9810 Varol Avenue
Chertsworth, CA 91311


Waters PowerLine Dissolution System including Waters 600E Gradient PowerLine Module, Waters 484 detector, the WISP 712D, and the Hanson Transfer Controller. (Hanson SR-2 Dissolution Bath courtesy of Hanson Research Corporation.)



The Hanson Transfer Controller (TC), activated by the 600E, removes samples from the dissolution bath at times specified in the Dissolution Table editor. Autoprobes, equipped with 5–10 micron filters, are lowered simultaneously into each of the dissolution vessels at the specified sampling time. Excess filtered sample is withdrawn into six volumetric cylinders under vacuum. Each of the six cylinders have a liquid leveling rod to determine the sample volume that will be dispensed into the 712D autosampler carousel. Using positive pressure, the TC delivers the excess sample back into the dissolution vessels, simultaneously backflushing the autoprobe filters. The remaining measured and filtered sample is then dispensed into empty, capped one ml vials inside the active 712D sample processor.

Dissolution sampling times are entered in the Dissolution Table editor of the 600E Controller according to the protocol specified by the dissolution test. The injection sequence consisting of samples and standards, and chromatography run times are entered in the Program Methods editor.

The 600E dissolution software determines whether an HPLC injection can be made without interfering with the acquisition of dissolution samples. If it can, the HPLC injection is completed and the 712D sample carousel is rotated to the correct position to collect the dissolution samples. On the other hand, if the system determines that it cannot make the next scheduled injection before the next set of dissolution samples are acquired, it interrupts the injection sequence to rotate the 712D carousel to the correct dissolution sample acquisition position. Once dissolution samples are acquired, the carousel is rotated back to the position it was in before the sample transfer



Waters 600E PowerLine Module Autoinjector Setup Screen Editor allows you to activate dissolution software. In the pharmaceutical industry, it is common to periodically run calibration standards after a predetermined number of samples have been analyzed. A convenient feature of the 600E PowerLine control module is the Auto Standards routine, which allows the user to program standard vial injections between a set number of sample vials.



Waters 600E Dissolution Table screen editor allows entry of sampling times. The automated dissolution system uses the 96 vial carousel in the 712D autosampler to accommodate the many samples taken during a dissolution test. The Dissolution Table editor provides the user with documentation of sample times and vial location of dissolution samples. Typically, the first six positions of the carousel are reserved for multilevel calibration standards.

took place and the injection sequence continues. In all cases, the 600E software prioritizes dissolution sample acquisition over that of a sample or standard injection.

There are cases where it is possible that the samples scheduled to be injected into the HPLC have not yet been acquired. This may happen in a situation where the dissolution sampling intervals are hours apart and the chromatography is programmed at minutes per injection. In this case, the dissolution software recognizes that the samples scheduled for analysis have not been acquired and will delay injections from these vials until they have been filled with dissolution samples.

Through intelligent software routines and advanced sample handling techniques, the Waters Automated Dissolution/HPLC system provides greater confidence in sampling and analytical results. Should verification of specific sample data be necessary, identified samples can be re-injected, since all dissolution samples are retained until the analyst approves all the data.

Combined Dissolution Testing and Content Uniformity for Quality Control

In addition to dissolution analyses, content uniformity determinations are also routinely performed in the pharmaceutical laboratory. Content uniformity determinations are usually performed separately by HPLC techniques, while dissolution results are more often based on spectrophotometric techniques. Waters Automated Dissolution System lets the pharmaceutical chemist combine dissolution testing and content uniformity determinations using a single automated method.

In combining these two tests, three criteria must be met. First, two dissolution baths must be used. Most dissolution tests require that only six samples be analyzed; however content uniformity determinations are usually performed on a minimum of ten samples, necessitating the use of a second dissolution bath (six additional samples).

Second, content uniformity determinations must be performed when the active ingredient(s) are 100% dissolved.

Third, the same HPLC method must be capable of analyzing dissolution samples as well as performing content uniformity determinations. If all three of these criteria are met, content uniformity determinations can be made utilizing either the same samples analyzed for the dissolution test (if actives are known to be completely dissolved), or samples withdrawn from the dissolution test vessels at a later time point (if the last dissolution sampling time is known to occur before the actives are 100% dissolved).

To demonstrate the feasibility of combining these two test procedures, identical HPLC methods (See Figure 1) were used for both the dissolution test samples and content uniformity determinations for the antihypertensive product Timolide®. Timolide contains hydrochlorothiazide and timolol maleate at 25 mg and 10 mg per tablet respectively.


The Waters Automated Dissolution System was used to control the sampling of 12 dissolution test vessels, transfer the samples to the active 712D autosampler, and analyze the dissolution sample and standard solutions. In this case, the Timolide actives were completely dissolved at the final dissolution test time point. Content uniformity determinations were based on the HPLC results obtained from 10 of the 12 dissolution samples. Additional content uniformity samples prepared offline by accepted United States Pharmacopeia (USP) methodology, were manually loaded into the WISP carousel for comparison with the dissolution samples (see Table 1).

Conclusion

The Waters PowerLine System for Dissolution testing brings a new dimension to dissolution analyses. Its ease of use, versatility, precision and reliability surpass all other dissolution sample assay techniques. It is the only online dissolution HPLC system that performs a complete drug product release study. Now for the first time researchers can execute extended dissolution profiles, content uniformity studies and composite drug assays on a single automated system providing them with the specificity and sensitivity demanded by today's drug formulations.

The data in Table 1 indicates that the automated method (different sample preparation method, but same chromatographic method of analysis), gives a similar mean and variance indicating that this method does not significantly affect the outcome of the content uniformity test. From the analysis of the means and variances by T-tests and F-tests respectively, we conclude that the manual method of sample preparation for content uniformity could be substituted with the automated dissolution method without significant bias to the results.

Figure 1: Waters 860 Networking Workstation Chromatography Sample Report of Timolide Separation



HPLC separation of Timolide active ingredients, hydrochlorothiazide and timolol maleate on a Waters μ Bondapak C₁₈ column. The mobile phase consisted of water, acetonitrile, methanol and pH 3.0 phosphate buffer (38:8:2:2). The separation was monitored at 295 nm as specified in USP XXI.

Table 1: Comparison of Manually Prepared Content Uniformity Sample Data to the Combined 10 Vessel Dissolution Data

Hydrochlorothiazide				Timolol Maleate			
Tablet	Manual mg/ml	Tablet	Auto mg/ml	Tablet	Manual mg/ml	Tablet	Auto mg/ml
1	25.05	11	24.53	1	9.88	11	9.93
2	25.38	12	24.58	2	10.01	12	9.92
3	25.10	13	24.62	3	9.80	13	9.78
4	24.93	14	25.23	4	9.81	14	9.99
5	25.03	15	24.78	5	9.97	15	10.00
6	24.83	16	25.34	6	9.74	16	10.20
7	25.23	17	24.61	7	9.98	17	9.63
8	25.05	18	25.50	8	9.85	18	9.79
9	25.18	19	25.21	9	9.89	19	10.06
10	25.88	20	25.25	10	10.30	20	9.93
AVG	25.17		24.96		9.92		9.92
SDV	0.28		0.35		0.15		0.15
RSD	1.11		1.42		1.51		1.52

Dissolution Testing Requiring High Sensitivity Analysis


The trend in the pharmaceutical industry is toward more complex formulations, lower dosages, and more sophisticated drug delivery systems such as controlled release products. As a result, pharmaceutical researchers are moving away from the traditional analytical spectrophotometric techniques and are turning to HPLC for dissolution analyses and quality control. The excellent specificity of HPLC makes it a more sensitive and reliable analytical tool than spectrophotometry. In addition, the USP accepts HPLC as a method of analysis for dissolution samples.

Dissolution Testing of Oral Contraceptives

Sensitive and reproducible HPLC separations of 23 different commercial oral contraceptive tablets, containing norethindrone and ethinyl estradiol at very low levels, have been described by Papas, et al.¹ Their method involved dispersing the tablet in 10 ml of a mixture of methanol and water (80%/20%). This corresponded to working concentrations of 100 µg/ml and 5 µg/ml for a tablet containing 1 mg norethindrone and 0.05 mg ethinyl estradiol respectively.

Dissolution testing of oral contraceptives was modified by the USP² in 1988 to reflect the use of HPLC as a method of analyzing dissolution samples. Using this method the tablets are dissolved in 600 ml of dissolution media and the samples analyzed by HPLC. The formulation described above would provide working concentrations of 1.7 µg/ml for norethindrone and 0.08 µg/ml for ethinyl

Figure 2: Waters 860 Expert Ease System Suitability Report



System Suitability Report of Oral Contraceptive Standard Injection. The system suitability report for a single vial injection prints the chromatogram and reports the USP parameters such as resolution, plate counts, and tailing factors for all peaks identified in the calibration table. In addition, a summary report is also printed for multiple injections of standard vials, reporting means, standard deviations, and relative standard deviations of area, height, retention time, tailing factor, resolution and plate count.

estradiol. This represents a 60-fold dilution over the methodology originally described by Papas.

To demonstrate the sensitivity of Waters 484 Tunable Absorbance Detector, Waters Pharmaceutical laboratory performed the dissolution test outlined in the USP for norethindrone and ethinyl estradiol tablets.

In this experiment, the unique capabilities of the PowerLine Dissolution/HPLC interface were coupled with those of the 484 absorbance detector. The optics of the 484 detector are controlled by the PowerLine 600E system control module by means of an IEEE-488 communications bus. Wavelengths and sensitivity settings are

programmable, offering the chromatographer the ability to optimize the detector output.

A standard solution was prepared representing the concentration of a typical oral contraceptive dissolution sample. The HPLC analysis required detection of these compounds at concentrations below 1.7 µg/ml and 0.08 µg/ml for norethindrone and ethinyl estradiol respectively.

¹Papas, Andrew N., Marchese, Salvatore M., and Delaney, Michael T., *Rapid Determination of Norethindrone and Ethinyl Estradiol in Oral Contraceptive Tablets by Reversed-phase Liquid Chromatography*, LC, Vol 3 No. 4, (April, 1985)

²Norethindrone and Ethinyl Estradiol Tablets, USP XXI, Seventh Supplement, 2813 (1988)

The USP protocol suggests 100 μ l injection sizes which is equivalent to 170 ng of norethindrone and 8 ng of ethinyl estradiol injected onto the HPLC column. Given the poor spectral characteristics of these two compounds, analytical requirements dictate an HPLC method that can provide a stable UV detector response, at high sensitivity, when operating at 200 nm.

Dissolution data was reprocessed on the Waters 860 networking computer system. Accurate and reproducible integration was obtained on peaks representing as little as a two millivolt response.

Figure 2 shows a system suitability report generated by the Waters 860

Networking Computer System, illustrating a typical separation of a synthetically prepared oral contraceptive dissolution sample using a Waters 484 absorbance detector. Five replicate injections of a standard solution resulted in RSD's of 2.5% for both actives, below the 3.0% specification set in the USP monograph.

In addition, the system suitability report shows the theoretical plates for ethinyl estradiol at 10,551 (USP states 7,000 minimum) and the resolution between the two actives at 2.8 (USP states 1.5 minimum). The tailing factor for both actives was well below the maximum of 2.0 described in the USP monograph.

Conclusion

Analysis of dissolution samples by chromatographic techniques can put a burden on the HPLC system. Dilute solutions of active ingredients must be capable of being reproducibly analyzed. The Waters Automated Dissolution System is capable of providing the user with convenient methods setup, variable sample injection volumes and detector options to meet the demands for high sensitivity analyses.

The Waters 845 Chromatography Workstation

The Ideal Chromatography Workstation for the Pharmaceutical Laboratory.

The Waters 845 Chromatography Workstation is one of the most powerful laboratory data handling and control systems available today. It is a quantitative tool that enables the pharmaceutical chemist to measure area and/or peak height and to accurately calculate the true amount of material present in samples.

The Waters 845 combines the power of the Digital VAXstation 2000 or the VAXstation 3100 Desktop Computer with the latest enhancements of Waters Expert Ease Software. The 845's window environment enhances the system's multitasking capabilities.

Waters Expert Ease Software

Waters Expert Ease software is a direct access, menu driven software package optimized for gas, liquid, and ion chromatography applications. (There is an optional software package designed for gel permeation applications.) It automatically collects data, integrates peaks, stores methods, plots in real-time, runs statistics and generates graphics and custom reports.


The Waters Expert Ease software offers internal or external, single or multi-point calibrations up to 99 levels. It allows linear, quadratic, cubic, point-to-point and segmented curve fits. Calibration curves and regression statistics can be printed or plotted for presentation or documentation purposes.

The 845 offers default integration parameters to meet the needs of most pharmaceutical applications and integration timed events are provided for complex chromatographic separations.

Waters Easy-View software allows simple mouse driven manual integration and chromatogram zooming to see and accurately integrate even the smallest peaks.

Waters Custom Report Manager allows the user to tailor chromatographic report formats, report sequences and the report print location.

Waters Expert Ease Custom Convert capability allows transfer of valuable chromatographic results into RS-1, Lotus 1-2-3®, 20/20™ or ASCII



Waters 845 Chromatography Workstation featuring the VAXstation 3100 personal computer from Digital Equipment Corporation.

formats for further statistical analysis. Complete methods documentation is available for GMP/GLP purposes.

Waters LAC/E Technology

Waters Laboratory Acquisition and Communications Environment (LAC/E™) module provides advanced system control and data acquisition capabilities to increase your laboratory's productivity.

The LAC/E module controls all Waters pumping modules including isocratic, single-pump gradient and multi-pump gradients. It also controls Waters family of WISP sample processors as well as Waters detectors (410 Differential Refractometer, 484 and 490 Absorbance Detectors). LAC/E serial board options offer the ability to connect up to eight additional printers offering a more organized reporting structure for a multiple HPLC/GC laboratory.

Waters SAT/IN™ module provides two channels of high resolution A/D with separate time bases to meet the needs of capillary GC, as well as LC, IC and GPC data collection.

Introducing Waters Expert Ease System Suitability Software Option


Waters Expert Ease System Suitability software option provides the pharmaceutical chemist with critical HPLC separation data needed to validate chromatographic system performance. The software calculates plate count, resolution, and tailing factors, determines percent relative standard deviation for retention time, area, height, plate count, resolution and tailing factors for single vials or entire multi-methods.

Waters Custom Convert manager transfers summary reports into third party software formats (RS-1, Lotus 1-2-3, 20/20 and ASCII) for further statistical analysis and graphic presentations as shown in Figure 3.



Through advanced communication protocols, the Waters 845 Chromatography software provides unlimited methods development options and sophisticated multiple detector control.

Figure 3: Enhanced Graphic Display of System Suitability Results



Analytical results can be converted to a third party software format for enhanced visual displays. In this example, system suitability results from the analysis of a calibration standard containing hydrochlorothiazide and timolol maleate were converted to a format compatible with Access Technology's 20/20 software, and then graphed to demonstrate advanced reporting capabilities of the 845 system.

Waters Extends its Family of Sample Processors

Waters 715 and 1700 WISP Sample Processors for Advanced Sample Handling


Waters has introduced two new products to extend its line of auto injectors, the 715 Ultra WISP and the 1700 WISP Sample Processors. The 715 offers full, flexible programmability with easy-to-use menus and CRT display and an all-electric valve design. It is ideal for applications involving reagent additions and derivatization. It is compatible with automated workstations such as the Millilab™ 1A Sample Preparation Workstation.

The WISP 715 is compatible with Waters PowerLine HPLC system for single keyboard control and documentation of all operating system parameters. In addition, it features an

exclusive fluid path design that ensures accuracy and reproducibility for all injector volumes while eliminating cross-contamination and carry-over.

A choice of sample injection volumes, from 1 microliter to 2000 microliters and carousels with either 48 four-milliliter or 96 one-milliliter vials, are available.

The Waters 1700 WISP Sample Processor provides automated sample preparation with Waters Millilab Workstation, plus direct sample injection for on-line HPLC analysis. This innovative system performs filtration, solid phase extraction, liquid-liquid extraction, mixing, addition of internal standards, heating, cooling, evaporation, trace enrichment, plus many other repetitive sample cleanup routines completed unattended.




The CRT in the 715 WISP verifies system status and conveniently programs priority samples as needed without interfering with the run in progress. In addition, system messages and system diagnostic routines make troubleshooting and maintenance easy.

Pharmaceutical Notes Free Column Offer

We would like to receive editorial contributions from you. If you have an HPLC, IC or GPC application or an improved sample preparation technique that you would like to share with our readership, we would like to include it in a future issue of Pharmaceutical Notes.

In exchange, we will send you a free analytical reverse-phase C₁₈ column of your choice in a Nova-Pak™, μBondapak™ or Resolve™ chemistry if your contribution is accepted.

Please send your suggestions to:
Carole Wade-Clark MC
Waters Division of Millipore
34 Maple Street, Milford, MA 01757



How to receive future issues of Pharmaceutical Notes

Return the enclosed reply card today to be sure you don't miss a future issue of Pharmaceutical Notes. For further information concerning any of the articles in this issue, contact your nearest Waters technical representative or call 508-478-2000.

Waters Chromatography Division
Millipore Corporation
34 Maple Street
Milford, MA 01757
508-478-2000

U.S.A. Sales Offices

New England Waters Chromatography Division / Millipore Corporation / 34 Maple Street / Milford, MA 01757 / Tel. 508-478-2000

Chicago Waters Chromatography Division 1555 Mittel Boulevard / Wood Dale, IL 60191 / Tel. 312-766-6060

New Jersey Waters Chromatography Division 35 Airport Park / Suite 230 / Morristown, NJ 07960 / Tel. 201-285-1404

North Carolina Waters Chromatography Division / 117 Edinburg Drive / Suite 206 Cary, NC 27511 / Tel. 919-469-2501

San Francisco Waters Chromatography Division c/o Millipore Corporation / 448 Grandview Drive / South San Francisco, CA 94080 Tel. 415-952-9200

Southwest Waters Chromatography Division 2626 South Loop West / Suite 170 Houston, TX 77054 / Tel. 713-665 7310

Washington Waters Chromatography Division 3702 Pender Drive / Fairfax, VA 22030 Tel. 1-800-526-0771

Overseas Offices

Australia Millipore Pty. Ltd. / Waters Chromatography Division / Locked Bag 18 87-89 Mars Road / Lane Cove, N.S.W. 2066 Tel. Inwats (008) 222-111 / In Sydney Area (02) 427-0511 / Telex 25934

Austria Millipore Ges.m.b.h. / Hietzinger Hauptstrasse 145 / A-1130 Wien Tel. (0222) 82 89 26 / Telex 131464 Milliv A

Belgium and Luxembourg Millipore-N.V. Waters Chromatography Division / Raketstraat 60 / Rue de la Fusée 60 / 1130 Bruxelles-Brussels Tel. (02) 242.1740 / Telex 22798 milwat b

Brazil Millipore Industria e Comercio Ltda. Waters Chromatography Division / Rua Prof. Campos de Oliveira, 430 / CEP 04675-011 Postal 19065 / Sao Paulo / Tel. (011) 548-7011 Telex (11) 25078

Canada Waters Chromatography Division 3688 Nashua Drive / Mississauga / Ontario L4V-1M5 / Tel. 416-678-2151 Toll-free 1-800-268 4881 / Telex 06-968513

Denmark Millipore A/S / Waters Chromatography Division / Roskildevej 342 2630 Tåstrup / Tel. (42) 52 88 11 Telex 33544

Finland Millipore Oy / Waters Chromatography Division / Ruukinkuja 1 / 02320 Espoo Tel. 90-801 90 77 / Telex 12-2123 milpor si

France Millipore S.A. / Waters Chromatography Division / Boite Postale 307 / 78054 Saint Quentin / Yvelines Cédex / Tel. (1) 30.45.70.00 Telex 698 371 F

India Waters Instruments (India) Pvt. Ltd. 1st Floor Shankar House / 1 Rajmahal Vilas Extn. / Mehra Circle / Bangalore 560 080 Tel. (812) 34-9203 / Telex 8452028 SWPL IN

Italy Millipore S.p.A. / Waters Chromatography Division / Via Achille Grande, 23/20090 Vimodrone MI/Tel. (2)250781 Telex 312284 milspai/Telefax (2)2650324

Italy Millipore S.p.A. / Waters Chromatography Division / Viale del Caravaggio 107 / 00147 Roma RM / Tel. (06) 51.32.452 / (06) 51.31.593 (06) 51.31.512 / Telefax (06) 51.26.267

Japan Nihon Millipore Ltd. / Waters Chromatography Division / No. 5 Koike Bldg. 3-12, Kita-Shinagawa 1-chome / Shinagawa-Ku, Tokyo 140 / Tel. 81-3-471-7191 Telex 781-0232416i

Mexico Millipore S.A. de C.V. / Waters Chromatography Division / Avenida Ingenieros Militares 85 / Mexico 11230 D.F. Tel. 576-9688 / Telex 1777 442

The Netherlands Millipore B.V. / Waters Chromatography Division / Penningweg 33 Postbus 166 / 4870 AD Etten-Leur (N.B.) Tel. (0 16 08) 22000 / Telex 54 840 milwa nl

Norway Millipore AB / Waters Chromatography Division / Enebakkveien 133 / 0680 Oslo 6 Tel. (02) 67 82 53

Peoples Republic of China (Beijing) Nihon Waters Service Center Rm. 638, Xiyuan Hotel Erligou Xijiao, Beijing Tel. 8022018 / Telex 22685

Peoples Republic of China Millipore China Ltd. Waters Chromatography Division / Hope Sea Industrial Centre / 26 Lam Hing Street / Suite 105 / 1st Floor / Kowloon Bay / Hong Kong / Tel. 3-796 3736 Telex 3-7967374

Puerto Rico Millipore Corporation Puerto Rico Branch / Waters Chromatography Division KM 7.7 / State Road 172 / Cidra, Puerto Rico 00639 / Tel. 809-747-8444 Telex 385 72 76

Spain Millipore Ibérica S.A. / Waters Chromatography Division / Avda del Llano Castellano, 13 / 28034 Madrid Tel. (91) 729.03.00 / Telex 23545 milli e

Spain Millipore Ibérica S.A. / Waters Chromatography Division / Entenza 28 08015 Barcelona / Tel. (93) 325.96.16 Telex 50524 wtrs e

Sweden Millipore AB / Waters Chromatography Division / Box 233 / Långebergsgatan 30 421 23 Västra Frölunda / Tel. (031) 28 98 60 Telex 21064 milpor s

Sweden Millipore AB / Waters Chromatography Division / Grahmsvägen 15 172 46 Sundbyberg / Tel. (08) 98 89 60

Switzerland Millipore AG / Steinackerstrasse 11 8302 Kloten / Tel. 41-181-41-363 / Telex 56067

The United Kingdom and Ireland Millipore (U.K.) Limited / Waters Chromatography Division The Boulevard / Ascot Road Croxley Green, Watford Hertfordshire WD1 8YW Tel. 44-923 816375 Telex: 851-24191 MILLIPOR G

USSR, Middle East, Eastern Europe and Africa Millipore Ges.m.b.h. / Hietzinger Hauptstrasse 145 / A-1130 Wien / Tel. (0222) 82 89 26 Telex 131464 Milliv A

West Germany Millipore GmbH / Waters Chromatography Division / Hauptstrasse 81-85 6236 Eschborn / Tel. 06196-494-0 Telex 4072704 mil d

For all other countries: Waters Chromatography Division / Millipore Corporation / 34 Maple Street Milford, MA 01757 / USA / Tel. 508 478-2000 Telex 174166

Waters, Millipore, LACE/F, SAT/IN, Expert Ease, WISP, μ Bondapak, Resolve, Nova-Pak, Maxima and MilliLab are trademarks of Millipore Corporation.

Lotus 1-2-3 is a registered trademark of Lotus Corporation.

20/20 is a trademark of Access Technology, Inc.

Timolide is a registered trademark of Merck, Sharp & Dohme.

Pharmaceutical Notes is published by the Pharmaceutical Group of Waters.

Waters
Division of MILLIPORE