

**THE USE OF ON-LINE SAMPLE  
CONCENTRATION TO INCREASE  
THE SENSITIVITY OF CAPILLARY  
ELECTROPHORESIS**

**Presented by:**

**Michael Merion, Ruedi H. Aebersold, Martin Fuchs  
Waters Division of Millipore,  
Milford, MA USA  
HPCE'91**


# THE USE OF ON-LINE SAMPLE CONCENTRATION TO INCREASE THE SENSITIVITY OF CAPILLARY ELECTROPHORESIS.

Michael Merion. Ruedi H. Aebersold. Martin  
Fuchs.


Millipore Corporation. Waters  
Chromatography Division. Milford, MA  
01757.

Capillary Electrophoresis is a high resolution separation technique that is of growing importance in the bioanalytical laboratory. Using on-line UV detection, as little as 10pg of a peptide may be detected. However, because injection volumes are limited to 1-10nl, the lower limit for detection with respect to sample concentration is about 1ug/ml. In order to lower this limit we have constructed a fused silica capillary that contains a packed bed of C18 chromatographic material with a length of 1-2mm at the injection end. Large volumes of dilute peptide samples were injected using conditions where the sample binds to the C18. After sample loading, a small volume of organic solvent was injected and electrophoresis conducted in the normal fashion. Using this technique, peptides have been detected at original sample concentrations down to 1ng/ml.



PHOTOMICROGRAPH OF  
PRE-CONCENTRATION CAPILLARY



## Pre-Concentration Capillary



## EFFECT OF ORGANIC SOLVENT ON PLATE COUNT



### Conditions:

Capillary: 100  $\mu$ m x 60 cm, 1.5 mm Delta-Pak™  
Buffer: 25 mM Na Citrate, pH 3.5  
Eluent: 25 mM Na Citrate, pH 3.5/Indicated  
Organic 25/75  
Sample: 2  $\mu$ g/ml Neuromedin in Buffer  
Injection: 500 seconds, 8 kV  
Elution: 20 seconds, 7 kV  
Run: 8 kV



# Reproducibility Using A Preconcentration Capillary

## Experimental:

Capillary: 100  $\mu$ m x 60 cm , 1.5 mm Delta-Pak  
Buffer: 25 mM Sodium Citrate, pH 3.5  
Eluent: 10 mM Sodium Citrate, pH 3.5 /  
Acetonitrile (25/75)  
Sample: 2  $\mu$ g/ml Neuromedin in Buffer  
Injection: 500 seconds at 8 KV  
Elution: 20 seconds at 7KV  
Run: 30 minutes at 8 KV

<u>Migration Time</u>	<u>Area</u>
26.13	90392
25.92	89711
25.60	89531
25.41	92269
25.20	91398
mean = 25.65	90,660
S.D. = 0.336	1037
%R.S.D = 1.14	1.30

# EFFECT OF CAPILLARY INSIDE DIAMETER ON PLATE COUNT



In order to examine possible sample carry-over, a second elution with organic solvent was performed. No additional peaks were observed.




## Conditions:

Capillary: 75 µm x 60 cm, 0.6 mm Delta Pak™  
Buffer: 25 mM Na Citrate, pH 3.5  
Eluent: 25 mM Na Citrate, pH 3.5/ Acetonitrile 25/75  
Sample: 2 µg/ml Neuromedin in Buffer  
Injection: 100 seconds at 12 kV  
Elution: 10 seconds at 7 kV  
Run: 12 kV

## **Extending the pH Range for Pre-Concentration**

**Our first experiments employed a silica based packing material for pre-concentration. This limits the pH range to 3-7. In order to increase the utility of this technology, we next used a polymeric hydrophobic packing material. This extends the pH range to 3-10. The properties of this capillary are otherwise very similar to the silica based capillary.**






Conditions:

Capillary: 100  $\mu$ m x 60 cm, 1.5 mm Delta Pak™  
 Buffer: 25 mM Na Citrate, pH 4.0  
 Eluent: 25 mM Na Citrate/Acetonitrile  
           Ratio As Indicated  
 Sample: 10  $\mu$ g/ml Tryptic Digest Cytochrome c  
           in Buffer  
 Injection: 500 seconds, 8 kV  
 Elution: 20 seconds, 7 kV  
 Run: 8 kV

# ANALYSIS OF ANGIOTENSIN I AT pH 10.5




## Conditions:

Capillary: 75  $\mu\text{m}$  x 60 cm. 1 mm Polymeric  
Packing  
Buffer: 50 mM Na Borate, pH 10.5  
Eluent: 50 mM Na Borate, pH 10.5/Acetonitrile  
25/75  
Sample: 1  $\mu\text{g}/\text{ml}$  Angiotensin I in Buffer + 1M  
AccuPure™ Z1 Methyl  
Injection: 999 seconds at 15 kV  
Elution: 10 seconds at 7.5 kV  
Run: 15 kV

DETECTION OF TRACE IMPURITIES  
USING PRE-CONCENTRATION CAPILLARY

Angiotensin II 50  $\mu\text{g/ml}$

1% Impurity




Conditions:

Capillary: 75  $\mu\text{m}$  x 60 cm, 0.6 mm Delta-Pak™  
Buffer: 25 mM Na Citrate, pH 4.0  
Eluent: 25 mM Na Citrate, pH 4.0/Acetonitrile  
25/75  
Sample: 500 ng/ml Neuromedin + 50  $\mu\text{g/ml}$   
Angiotensin II in Buffer  
Injection: 500 seconds at 12 kV  
Elution: 10 seconds at 7 kV  
Run: 12 kV

# CONCLUSIONS

- \* Adding On-Line Sample Concentration to Capillary Electrophoresis Can Greatly Increase Sensitivity With Respect To Sample Concentration. 250x Demonstrated In This Paper.
- \* The Pre-Concentration Capillary Behaves in a Linear Fashion with Respect to Sample Injection.
- \* Separations Using the Pre-Concentration Capillary are Reproducible with Respect to Migration Time and Peak Area.
- \* Separations Using the Pre-Concentration Capillary Result in High Resolution
- \* Using the Pre-Concentration Capillary, Contaminants May Be Detected Well Below the 0.02% Level.
- \* The Use of Step Gradient Elutions With the Pre-Concentration Capillary Can Result In 2 Dimensional Separations.

# USE OF PRE-CONCENTRATION CAPILLARY FOR CE FRACTION ANALYSIS



## Conditions: Preparative Separation

Capillary: 75  $\mu$ m x 60 cm AccuSep™ Fused Silica Capillary  
 Buffer: 25 mM Na Citrate, pH 4.0  
 Sample: 200  $\mu$ g/ml Angiotensin I in Buffer  
 Injection: 20 seconds, Hydrostatic  
 Run: 15 kV


## Conditions: Fraction Analysis

Capillary: 75  $\mu$ m x 60 cm. 1 mm Polymeric Packing  
 Buffer: 25 mM Na Citrate, pH 4.0  
 Eluent: 25 mM Na Citrate, pH 4.0/Acetoni-  
 25/75  
 Sample: Collected Fraction + 1M AccuPur  
 Z1 Methyl  
 Injection: 999 seconds at 15 kV  
 Elution: 10 seconds at 7.5 kV  
 Run: 15 kV

## **USE OF THE PRE-CONCENTRATION CAPILLARY FOR ANALYSIS OF CE FRACTIONS**

**It is possible to use Capillary Electrophoresis as a preparative technique. However, because small sample masses (1 pg- 1 ng) are collected in several  $\mu$ ls of buffer, the concentration of the collected fraction is generally too low to be further analyzed using CE. The use of the pre-concentration capillary enables the user to analyze CE fractions directly.**

## INFLUENCE OF BACK PRESSURE ON ELECTROOSMOTIC FLOW



In order to determine the impact of the packing material on the rate of osmotic flow, 3 capillaries with dimensions of  $100 \mu\text{m} \times 60 \text{ cm}$  were tested. The pack pressure was determined using an HPLC pump with a flow rate of  $100 \mu\text{l}/\text{minute}$ . The rate of osmotic flow was determined using a neutral marker.

Capillary 1: No packing material


Capillary 2: 1 mm packed bed of Delta-Pak™

Capillary 3: 4 mm packed bed of Delta-Pak™

# Linearity Using A Preconcentration Capillary


## Experimental:

Capillary: 100  $\mu\text{m}$  x 50 cm , 0.66 mm Delta-Pak  
Buffer: 25 mM Sodium Citrate, pH 3.5  
Eluent: 12.5 mM Sodium Citrate, pH 3.5 /  
Methanol (50/50)  
Sample: 2  $\mu\text{g/ml}$  Neuromedin in Buffer  
Injection: 500 seconds at 5 KV  
Elution: 20 seconds at 7KV  
Run: 10 KV





# ANALYSIS OF ANGIOTENSIN I AT pH 10.5



## Conditions:

Capillary: 75  $\mu$ m x 60 cm. 1 mm Polymeric  
Packing

Buffer: 50 mM Na Borate, pH 10.5

Eluent: 50 mM Na Borate, pH 10.5/Acetonitrile  
25/75

Sample: 1  $\mu$ g/ml Angiotensin I in Buffer + 1M  
AccuPure™ Z1 Methyl

Injection: 999 seconds at 15 kV

Elution: 10 seconds at 7.5 kV

Run: 15 kV