

MILLENNIUM APPLICATION

BRIEF


TOPIC: Capillary Electrophoresis of N-methyl-2-pyrrolidinone (NMP) in Etchants

AUTHOR: Stuart Oehrle

Date: 01/3/95

Introduction:

Analysis of NMP (N-methyl-2-pyrrolidinone) in etchants was done. Samples of three etchants (one of which was not NMP based) were diluted and analyzed using the Q4000e CIA. Since NMP is a neutral compound MECC (micellar electrokinetic capillary chromatography) was used. The structure of NMP is shown below. All samples were diluted (10uL-100mL) with the running buffer. Plastic volumetric flask were used for all dilution's and sample prep. Analysis of these mixes showed NMP present. Figure 1 is an electropherogram of a methanol-water blank used for marking the void.. Figure 2 is an electropherogram of a NMP standard. Figure 3 is an electropherogram of the quartz etchant showing no NMP present. This was expected since its composition is a mix of HF and NH_4F . Figure 4 and 5 of the intermediate and final etchants, show the presence of NMP. This was expected since they are composed of HF and NMP. The main interest here was in monitoring the level of NMP present as well as other cations and anions present in the etchants. Changeover time from cation analysis to MECC was less than 10 minutes.



NMP

EXPERIMENTAL:

MECC Analysis

System:	Q4000E
Electrolyte:	Sodium borate/Boric acid/SDS
Capillary:	75um X 60cm
Injection:	Hydrostatic, 30 seconds
Run Voltage:	+20kV
Detection:	Indirect UV at 185nm
Data:	Millennium 2010 Chromatography Manager, Ver. 2.1 with CIA Option

FIGURES:

Figure 1: Electropherogram of methanol-water void marker

Figure 2: Electropherogram of NMP standard

Figure 3: Electropherogram of quartz etchant ($\text{NH}_4\text{:HF}$)

Figure 4: Electropherogram of intermediate etchant (NMP:HF)

Figure 5: Electropherogram of final etchant (HF:NMP)

For Sample: H2O/MeOH Blank Vial: 1 Inj: 1 Chan: SATIN

Date Processed 12/28/94 02:03:34 PM

Channel Descr: Direct UV

Northern Kentucky Laboratory
Sample Information

Project Name: CIA HF_01
SampleName: H2O/MeOH Blank

Vial: 1

Injection: 1

Channel: SATIN

Date Acquired: 12/28/94 12:23:57 PM

Scale Factor: 1.00

Acq Meth Set: CE_Control_MECC

Processing Method: HF_MECC_01a

Sample Type: Hydrostatic Unknown

Volume: 0.00

Run Time: 10.0 min


Date Processed: 12/28/94 02:03:34 PM

Dilution: 1.00000

SampleName: H2O/MeOH Blank Electrolyte: 12.5 Boric Acid/2.5 NaBorate/50mM SDS

RunVoltage: +20kV Detection: Direct UV 185nm, 0.3 TC

InjectionMode: Hydrostatic, 30sec Capillary: 75um X 60cm



Peak Results

#	Name	Migration Time (min)	Area (uV*sec)	Time Corr. Area	Height (uV)
1	MeOH	4.370	173730	39755.24	26335

For Sample: NMP- (0 Vial: 4 Inj: 1 Chan: SATIN

Date Processed 12/28/94 02:04:25 PM


Channel Descr: Direct UV

Northern Kentucky Laboratory
Sample Information

Project Name: CIA_HF_01
SampleName: NMP- (0.01-10mL)
Vial: 4
Injection: 1
Channel: SATIN
Date Acquired: 12/28/94 01:03:03 PM
Scale Factor: 1.00
Acq Meth Set: CE_Control_MECC
Processing Method: HF_MECC_01

Sample Type: Hydrostatic Unknown
Volume: 0.00
Run Time: 8.0 min
Date Processed: 12/28/94 02:04:25 PM
Dilution: 1.00000

SampleName: NMP- (0.01-10mL)
Electrolyte: 12.5 Boric Acid/2.5 NaBorate/50mM SDS RunVoltage: +20kV
Detection: Direct UV 185nm, 0.3 TC InjectionMode: Hydrostatic, 30sec
Capillary: 75um X 60cm



Peak Results

#	Name	Migration Time (min)	Area (uV*sec)	Time Corr. Area	Height (uV)
1	NMP	4.852	1335716	275310.82	251106

For Sample: Quartz Etch (0.Vial: 5 Inj: 1 Chan: SATIN

Date Processed 12/28/94 01:59:09 PM

Channel Descr: Direct UV

Northern Kentucky Laboratory
Sample Information

Project Name: CIA_HF_01
SampleName: Quartz Etch (0.02-10mL)
Vial: 5
Injection: 1
Channel: SATIN
Date Acquired: 12/28/94 01:14:41 PM
Scale Factor: 1.00
Acq Meth Set: CE_Control_MECC
Processing Method: HF_MECC_01

Sample Type: Hydrostatic Unknown
Volume: 0.00
Run Time: 8.0 min
Date Processed: 12/28/94 01:59:09 PM
Dilution: 500.00000

SampleName: Quartz Etch (0.02-10mL)
Electrolyte: 12.5 Boric Acid/2.5 NaBorate/50mM SDS RunVoltage: +20kV
Detection: Direct UV 185nm, 0.3 TC InjectionMode: Hydrostatic, 30sec
Capillary: 75um X 60cm




Table 'Peak Results' contains no data.

For Sample: Int. Etch (0.02 Vial: 6 Inj: 1 Chan: SATIN

Date Processed 12/28/94 02:01:13 PM

Channel Descr: Direct UV

Northern Kentucky Laboratory
Sample Information

Project Name: CIA_HF_01
SampleName: Int. Etch (0.02-10mL)
Vial: 6
Injection: 1
Channel: SATIN
Date Acquired: 12/28/94 01:26:20 PM
Scale Factor: 1.00
Acq Meth Set: CE_Control_MECC
Processing Method: HF_MECC_01


Sample Type: Hydrostatic Unknown
Volume: 0.00
Run Time: 8.0 min
Date Processed: 12/28/94 02:01:13 PM
Dilution: 500.00000

SampleName: Int. Etch (0.02-10mL)

Electrolyte: 12.5 Boric Acid/2.5 NaBorate/50mM SDS RunVoltage: +20kV

Detection: Direct UV 185nm, 0.3 TC InjectionMode: Hydrostatic, 30sec

Capillary: 75um X 60cm



Peak Results

#	Name	Migration Time (min)	Area (uV*sec)	Time Corr. Area	Height (uV)
1	NMP	5.078	1465806	288639.30	252715

For Sample: Final Etch (0.0Vial: 7 Inj: 1 Chan: SATIN

Date Processed 12/28/94 02:01:37 PM

Channel Descr: Direct UV

Northern Kentucky Laboratory
Sample Information

Project Name: CIA_HF_01
SampleName: Final Etch (0.02-10mL)
Vial: 7
Injection: 1
Channel: SATIN
Date Acquired: 12/28/94 01:37:58 PM
Scale Factor: 1.00
Acq Meth Set: CE_Control_MECC
Processing Method: HF_MECC_01


Sample Type: Hydrostatic Unknown
Volume: 0.00
Run Time: 8.0 min
Date Processed: 12/28/94 02:01:37 PM
Dilution: 500.00000

SampleName: Final Etch (0.02-10mL)

Electrolyte: 12.5 Boric Acid/2.5 NaBorate/50mM SDS RunVoltage: +20kV

Detection: Direct UV 185nm, 0.3 TC InjectionMode: Hydrostatic, 30sec

Capillary: 75um X 60cm



Peak Results

#	Name	Migration Time (min)	Area (uV*sec)	Time Corr. Area	Height (uV)
1	NMP	4.975	1113166	223751.95	208249