

Waters

Ion Chromatography Method


Alkali and Alkaline Earth Cations, Ammonium and Amines 2000

Required Instrumentation:

Part / Number

Alliance, 2690 Separations Module
(with Column Heater, Seal Wash, Degasser)
432 Conductivity Detector
Bus SAT/IN Module
Millennium 32

271013
043061
073645
Consult Waters



Analysis Conditions:

Column: IC-Pak Cation M/D
Eluent: 3 mM HNO₃ / 0.1 mM EDTA
Back Conductivity: 1250 ± 50 μS
Degas: Continuous
Flow Rate: 1 mL / min
Back Pressure: 2100 psi
Temperature: 30°C (Column Heater); 35°C (Detector)
Injection: 100 μL
Needle Wash: 12% AcCN in DI Water
Detection: Indirect Conductivity
Base Range: 2000
Attenuation: 100 μS / Volt Unattenuated
Polarity: Negative

Eluent Preparation:

- 1) Into a 1 liter plastic volumetric flask, add
-0.029 g of EDTA, free acid.
- 2) Dissolve in 500 mL of DI Water with the aid of a stirring bar.
- 3) Add 30 mL of 100 mM HNO₃.
- 4) Dilute to volume with DI Water.
- 5) Vacuum degas through a 0.45 μ m aqueous compatible membrane.
- 6) Store in a plastic container at ambient temperature. Discard after 1 month.

Standard Preparation:

It is recommended that certified 1000 ppm anion standards be used with this method. If unavailable, see Reagent Section for uncertified standard preparation.

Prepare at least 3 mixed analyte standards, using plastic volumetric flasks, within the expected range of the sample analyte concentration. This method is linear from 0.05 ppm to 20 ppm for lithium, sodium, and ammonium, and 0.05 ppm to 50 ppm for potassium, magnesium, and calcium. Above these concentrations the response is off scale. After the multi-point calibration curve has been validated, a single point calibration within the expected analyte concentration is appropriate for future calibrations.

This method can be used for the analysis of Rb, Cs, Sr, and Ba.

Sample Preparation:

Determine the expected range of analyte concentration and other anionic component in the sample matrix. Sodium should be less than 20 ppm for best results.

If necessary dilute the sample with DI Water.

If the sample contains high amounts of neutral organics, or is highly colored, then pass the diluted sample through a C₁₈ Sep-Pak Cartridge. Cations pass through unretained. There may be residual Na contamination from the cartridge.

Samples containing suspended solids should be filtered through a 0.45 μ m aqueous compatible disk prior to injection. Failure to filter solids results in the risk of increased column backpressure.

Sample pH should be within 2 to 7 for best results, especially the alkaline earth cations. Samples with pH less than 10 is appropriate for the alkali cations, ammonium, and amines.

For sample with pH less than 2, dilute the sample 1:10 with DI Water; or treat the sample with an Alltech Maxi-Clean TMIC-OH Cartridge to remove anions and neutralize pH.

Millennium Data Processing Method:

IC Processing Method using Peak Apex for Retention Time

Integration Peak Width = 30.0 Threshold = 25--40
 Min Area = 3000 Min Height = 500

 Inhibit Intg. = 0 to 2 min

Calibration Averaging = None RT Window = 5%

 Update RT = Never


 Peak Match = Closest

 Quant By = Peak Area

 Fit Type = Linear, for multi-point calibration
 or Linear Through Zero, for single point

Report Analyte Name
 Analyte Retention Time
 Peak Area
 Amounts

Method Detection Limits:






Based upon this representative chromatogram using a 100 μ L injection, the estimated detection limits, as ppb, at 3 times signal to noise (S/N) are:

Lithium = 1	Sodium = 5	Ammonium = 5
Potassium = 15	Magnesium = 10	Calcium = 15

Lower detection limits can be achieved by using a 250 μ L injection.

Examples of Use:



Stock Reagent Preparation:

It is difficult to prepare a stock eluent for this column; it is best to prepare fresh working eluent.

In a 1 liter plastic volumetric flask, add 0.029 g of EDTA (as the free acid, not its salts) in 800 mL of DI water, and a stirring bar. Place on a magnetic stir plate and stir for 10 minutes.

While stirring add 30 mL of 100 mM Nitric Acid and mix for 5 minutes.

Remove the stir bar and fill to the mark with DI water.

Filter through a 0.45 μm aqueous compatible membrane filter before use. There may be some remaining white crystals on the filter(EDTA); this does not effect the performance of the eluent. Discard the filter.

Waters

Capillary Ion Analysis Methods