
METHOD DEVELOPMENT ON QUANTITATIVE SPE-UPLC-MS/MS ANALYSIS OF A THERAPEUTIC MONOCLONAL ANTIBODY IN SERUM

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Solid phase extraction

Oasis[®] WCX and MCX uElution Plates (Waters) were used to cleanup the serum digest samples. Specifically, 400–750 μL of diluted serum digest was loaded onto the plate. The analytes were eluted twice with 25 μL of 40/1/55 ACN/TFA/water for WCX. 25 μL of 40/5/55ACN/NH₄OH/water for MCX.

LC conditions

System: ACQUITY UPLC[®] system (Waters)
Column: ACQUITY UPLC BEH300™ C18, 2.1 x 150 mm, 1.7 μm
Column temp.: 35°C
Flow rate: 0.3 mL/min.
Gradient: (A) 0.1% (v/v) formic acid in water
(B) 0.1% (v/v) formic acid in acetonitrile
0% B for 1 min, followed by a linear gradient up to 35% B
Cycle time: 20 min.

MS conditions

Waters Xevo[™] TQ MS and Xevo[™] TQ-S tandem quadrupole mass spectrometers
ESI potential: 3.5 kV
Source: 90 °C
Collision energy: range from 16 eV to 50 eV
Cone voltage: range from 12 V to 60 V
Dwell time: range from 1 ms to 100 ms

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Heavy Chain:

EVQLVESGGGLVQPGGSLR (monkey peptide)
LSCAASGFNIK (mouse peptide, has Cys)
DTYIHWR
QAPGK (human/lama peptide)
GLEVVVAR (human/mouse peptide)
IPTNGYTR (bacterial peptide)
YADSVK (human/lama peptide)
FTISADTSK
NTAYLQMNSLR (has Met)
AEDTAVYYCSR (human/lama peptide, has Cys)
WGDDGF SSASTK (has Met, too long: 26 AA)
EEMTK (too short: 5 AA, missed cleavage)

Light Chain:

ASQDVNTAVAWYQQK (mouse peptide)
LLIYSASFYSGVPSR
FSGSR

(Peptides in red indicate good options.)

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
Fig.2 Collision energy was optimized by varying CE from 16 eV to 50 eV (2 eV interval). The actual optimal CE was close to the predicted CE using the formula.


Fig.3 Cone voltage was optimized (from 12 V to 60 V at 4 V interval).

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
Fig.4 Calibration curve of a signature peptide without serum digest background. 5 μM tryptic digested trastuzumab was diluted in 0.1% FA to various concentrations.

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


Fig.5 Ion exchange and reversed phase mixed-mode SPE and μElution plate

Fig.6 A signature peptide detected on Xevo[™] TQ-S before and after SPE sample cleanup. SPE method not only reduced the background noise, but also increased signal substantially.

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FTISADTSK: 485.2 > 721.4

Fig.7 MRM quantification of trastuzumab in human serum on a Trizac Xevo TQ-S platform. MRM chromatograms for the spiked trastuzumab (1, 5, 10 nM) and the corresponding ¹³C¹⁵N-isotopically labeled peptide spiked at the same concentration in human serum.

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