


Characterization of Plant Extract Isoflavonoids Using LC-PDA-Mass Spec

Michael P. Balogh

LC/MS
Waters Corporation

Published in
LC•GC, Vol 15, No. 5, May 1997



© 2000 Waters Corp

Introduction slide to the analysis of isoflavonoids in red clover using the Waters Integrity LC/MS System.

Natural Product Analysis

Red Clover Isoflavonoid Characterization

Chromatographic Approach

Column: Waters Symmetry C₈, 3mm x 150mm

Mobile Phase: AcCN / Water Linear Gradient
15/85 to 36/64 over 40 minutes

Flow Rate: 400 µL/min

Detector 1: PDA Scan from 200 to 600 nm


Detector 2: Mass Spec, Electron Ionization
Scan from 50 to 500 m/z

© 2000 Waters Corp

Here are the optimized chromatographic conditions for the analysis. Both PDA and MS detection are employed.

Natural Product Analysis

Red Clover Isoflavonid Characterization




© 2000 Waters Corp


The chromatographic separation of red clover extract. The PDA plot (lower chromatogram) and mass spec TIC (upper chromatogram) is shown here. Note that all of the starred compounds were identified using the commercially available Wiley Library.

Some Common Flavonoids - their derivatives and secondary metabolites


Flavones



Acacetin




Diosmetin




Medicarpin

Isoflavones



Biochanin A




Genistein

Shown here are some of the structures of the flavones and isoflavones examined in this study

Natural Product Analysis

Red Clover Isoflavonoid Characterization




© 2000 Waters Corp

Identification of daidzein, diodmetin, and medicarpin in the red clover extract from the west coast.

Natural Product Analysis

Red Clover Isoflavonoid Characterization

Identification of Close-Eluting Peaks by MS




© 2000 Waters Corp.

A close up of the separation shows a very small peak (A) identified as the aglycone daidzein at the trailing edge of the glycoside ononin. The mass spectra for each are significantly different in this case, which allows direct identification by library match from a single experiment. The molecular ion of the large peak (ononin) at m/z 430 is followed by the loss of glucose (M⁺-162) forming the base peak of the spectrum. The smaller peak (daidzein) produces the molecular ion at m/z 254 as the base peak of the spectrum.

Natural Product Analysis

Red Clover Isoflavonoid Characterization

PDA Differentiation of Isomers



© 2000 Waters Corp

Comparison of the acquired peak with the UV spectrum of diosmetin standard, indicates that the peak is probably pratensein. Both diosmetin and pratensein have a formula weight of 300.27.

Natural Product and Nutraceutical Futures

- ★ Increased use of Mass Spectrometry to compliment Photo Diode Array detection for Product and Component characterization
- ★ Need for Validated Analytical Methods for Product Target Compounds
- ★ Need for Standardization among Suppliers and Manufacturers for product consistency
- ★ More FDA Regulation in the Future?

© 2000 Waters Corp

Waters foresees increased use of HPLC and LC/MS for the analysis of nutraceuticals.

More validated analytical methods are required.

Suppliers and manufacturers need to follow the same quality assurance procedures to assure product consistency.

Will the FDA regulate the nutraceutical industry more closely in the future?

Successful Analyses

Application	Sample Prep	Column	Technology
Sugars	C ₁₈ Sep-Pak	Carbohydrate	HPLC/RI
EED's	Oasis HLB	SymmetryShield RP ₁₈ & RP ₈	HPLC/PDA LC/MS
PAH's	Oasis HLB	LiChrosphere PAH	HPLC/PDA & Fluorescence
Carbamates	Oasis HLB	Carbamate	HPLC/PCFD LC/MS
Pesticides	Oasis HLB	Symmetry C ₁₈	LC/MS
PolymerAdditives	Org. Extract	Symmetry C ₁₈	LC/MS
Nutraceuticals	Org. Extract	Symmetry C ₈	LC/MS

© 2000 Waters Corp

In summary, this seminar has demonstrated successful analyses for a number of applications utilizing a variety of Waters products and technologies.

Agricultural Chemicals

Residues in Food & Water

Fine Chemicals

Pesticides
Polymer Additives



Functional Foods & Dietary Supplements

Nutrients
Nutraceuticals

Environmental Contaminants

Endocrine Disruptors

© 2010 Waters Corp

Conclusion: Waters offers solutions to a wide range of analytical challenges in the industrial market place.