

# NALIDIXIC ACID ANTIBIOTICS BY LC/MS

1.0 minute gradient, 2.1 x 20 mm *IS*<sup>TM</sup> Separation

## CONDITIONS

Column: Atlantis<sup>TM</sup> dC<sub>18</sub>, 2.1 x 20 mm *IS*<sup>TM</sup>, 3.0 µm, (P/N: 186002058)

Mobile Phase A: Water

Mobile Phase B: Methanol

Mobile Phase C: 1% HCOOH in Water

Flow Rate: 0.4 mL/min

Gradient:	Time (min)	Profile %A	Profile %B	Profile %C
	0.0	50	40	10
	1.0	30	60	10

Injection Volume: 2 µL

Sample concentration: 10 µg/mL

Temperature: 30°C

Instrument: Alliance<sup>®</sup> 2795 and Waters ZQ<sup>™</sup>

The top figure is the total ion current, followed by the extracted ion signals for each of the three analytes.

\*The "extra" peak in the cinoxacin panel is the isotope from oxolinic acid.

## MS CONDITIONS

Waters ZQ<sup>™</sup>

ES+

Capillary (kV) 3.5

Cone (V) 5.0

Extractor 3.0

RF Lens 0.1

Source Temp (°C) 150

Desolvation Temp (°C) 400

Cone Gas Flow (L/Hr) 50

Desolvation Gas Flow(L/Hr) 500

LM Resolution 15

HM Resolution 15

Ion Energy 0.5

Multiplier (V) 650

Compounds MW

1. Cinoxacin 262.2

2. Oxolinic Acid 261.2

3. Nalidixic Acid 232.2

