Poster #151 – ISC '96

21st International Symposium on Chromatography, Stuttgart, September 15th-20th, 1996

A Unique, New, Multipurpose Sorbent Engineered for **Reversed-Phase SPE**

Waters, Oasis, and Symmetry are trademarks of Waters Corporation. © Waters Corporation 1997



McDonald, P.D.*, Bouvier, E.S.P., Fisk, R.P., Glose, K.H., Cook, J.T., Iraneta, P.C., Capparella, M., Collins, P., Phillips, D.J., Cheng, Y-F., Brady, T.P., Doyle, T.V., Collamati, R.A., Walter, T.H., Bean, L.L., and Meirowitz, R.E.

Waters Corporation 34 Maple Street, Milford, MA 01757 USA



In the last thirty-five years, reversed-phase extraction on a solid support has gone through several stages of evolution. First, a low surface area, high pore volume, hydrophilic support matrix such as diatomaceous earth or kieselguhr, packed in large glass columns, was used to absorb an aqueous matrix, followed by elution of important analytes with an immiscible organic solvent [1]. In 1968, a new, hydrophobic, macroporous, large particle, lightly cross-linked, polystyrene polymer was used for the

extraction of steroids and their metabolites from 24-hour urine samples [2]. A major advance in miniaturization and efficacy occurred in 1977, with the invention of a small, pre-packed cartridge for SPE [3] containing an octadecylsilyl-bonded silica sorbent optimized for higher efficiency and recovery. Prototypes filled with a smallparticle, highly cross-linked, macroporous, divinylbenzene polymer [4] presaged introduction in the last five years of devices using similar polymer sorbents.

Abstract continued



While the use of more efficient, hydrophobic, silica- and polymerbased sorbents has been a major advance in the practice of reversed-phase SPE, there are still problems with the use of these materials. A major cause of difficulty has been the need to prewet these sorbents with a watermiscible organic solvent followed by conditioning with water or an aqueous buffer solution, before applying the aqueous sample matrix. If the sorbent dries out, even partially, in between the wetting/conditioning

and sample loading steps, retention of the analytes falls precipitously, capacity is severely compromised, and sample breakthrough may approach 100%.

A "next generation" polymeric sorbent (patent pending) has been designed and synthesized to solve these problems. It has a unique balance of lipophilicity, needed to effect reversed-phase, extractive adsorption, and hydrophilicity, required for water-wettability. This sorbent exhibits stronger retention and excellent recoveries

Abstract concluded



of a wider spectrum of analytes, even under conditions which cause phase-collapse of traditional C₁₈-bonded silicas or "drying" of hydrophobic polymers. The improved properties and dramatic performance advantages of this new replacement for traditional reversed-phase packings will be illustrated with examples, including the isolation of both highly polar and non-polar drugs from serum.

References:

- 1. P. K. Siiteri*, Steroids,* **2**, 687-712 (1963).
- H. L. Bradlow, "Extraction of Steroid Conjugates with a Neutral Resin," *Steroids*, **11(3)**, 265-272 (1968).
- 3. P.D. McDonald, R.V.
 Vivilecchia, and D.R. Lorenz, "Triaxially Compressed Beds," U.S. Patent #4,211,658.
- 4. D.J. Mackey and H.W. Higgins, J. Chromatography, 436, 243-257 (1988).

Problems in Reversed-Phase SPE



- Poor reproducibility
- Low recovery, often a result of breakthrough, esp., of polar compounds **D** Complicated methods development Tedious manipulations

A Cause of Poor Reproducibility/Recovery



A hydrophobic sorbent , once properly conditioned & equilibrated, must remain wet in order to retain analytes from an aqueous sample matrix.

"Wetting the sorbent. First, a solvent capable of wetting the alkyl chains is brought in contact with the sorbent. In the dry state, the bonded alkyl chains are twisted and collapsed on the surface; on contact with a suitable solvent, they are solvated so that they spread open to form a bristle. This ensures a good contact between the analyte and the bonded phase in the adsorption step. It is important that the sorbent remains wet in the following two steps [conditioning/equilibration & sample loading]. Failure to perform this stage will result in poor recoveries of the analyte."

p. 283 in: "Analyte Isolation by Solid Phase Extraction (SPE) on Silica-Bonded Phases, Classification and Recommended Practices, Technical Report, International Union of Pure and Applied Chemistry, Analytical Chemistry Division, Commission of General Aspects of Analytical Chemistry," M. Moors, D.L. Massart, and R.D. McDowall, Pure & Appl. Chem., 66(2), 277-304 (1994)



HLB = *Hydrophilic* – *Lipophilic Balance*

With a proper ratio of hydrophilic & lipophilic monomers, we have designed a porous polymer which not only retains a broad polarity spectrum of compounds, but also, once wetted with water, remains wet, even if the bulk fluid is removed by vacuum-assisted flow.

*Patent pending

Physicochemical Properties



This new polymer, now available packed in Oasis[™] HLB Cartridges [1cc, 30 mg, or 3 cc, 60 mg], has the following certified, typical properties:


Specific Surface Area: Average Pore Diameter: **Total Pore Volume:** Average Particle Diameter: Fines Content <10 μ m: Extractable Impurities:

 $833 \text{ m}^2/\text{g}$ 82 Å $1.34 \text{ cm}^{3/g}$ 31 µm 0.1% \leq ppm level

NOTE: A detailed Certificate of Analysis [COA] accompanies each package.

Effect of Drying - SPE Method






Drying Effect: HLB vs. C₁₈

		% Recovery			
		C ₁₈ 1cc, 100 mg		Oasis™ HLB Sorbent 1cc, 30 mg	
Compound	Concentration µg/mL	Wet	Dried 2 min	Wet	Dried 10 min
Acetaminophen	10	75.4	13.2	93.9	93.2
Ranitidine	10	94.2	64.9	95.6	95.7
Procainamide	10	91.8	43.9	96.6	96.8
Propranolol	20	89.6	44.7	92.3	92.8
Doxepin	5	84.7	69.6	91.8	92.3
	% RSD (n=3)	0.10-0.80	6.0-54	0.55-1.5	0.18-2.0

Drying: Oasis™ HLB Sorbent Unaffected



Procainamide	A Ranitidine	— Doxepin
Acetaminophen	→ Propranolol	



Ease of Use – Simple Procedure




With Waters Oasis™ HLB Cartridges: One sorbent/one protocol for many compound types – acids, bases, neutrals – makes methods development straightforward.

- □ Step 1: Wet w/ MeOH
- □ Step 2: Condition w/ water or buffer
- □ Step 3: Load sample (acidify, if needed)
- □ Step 4: Rinse w/ water or 5% MeOH
- □ Step 5: Elute w/ MeOH

Operation under continuous vacuum eliminates fiddling with stopcocks & worry about work interruptions!


Versatility + Uniformly High Recoveries

Wide variety of pharmaceutical compounds spiked in porcine serum at indicated number of μ g/mL were recovered by extraction with 1cc, 30 mg, OasisTM HLB Cartridges using procedure described above. [All % RSD's < 4, n = 6.]



14

Increased Retention for Polar Compounds



OasisTM HLB sorbent vs. trifunctional C₁₈-bonded silica Mobile Phase: 20 mM potassium phosphate, pH 7.0/ MeOH (95/5 v/v) Temp: 30° C. Flow: 1.0 mL/min Detection: UV@254 nm

			Retention	Retention Factor (k)		
		Compound	HLB	C ₁₈	Ratio	
						_
	1	Benzoic acid	7.36	2.05	3.6	
	2	Salicylic acid	34.0	4.97	6.8	
	3	Pyrogallol	32.8	2.60	12.6	
	4	Theobromine	35.9	15.5	2.4	
	5	Acetaminophen	69.0	9.93	6.9	
	6	Catechol	149	8.33	17.9	
ЛОС		соон он		CH ₃ N	NHCOCH3	
		OH OH				OH
1		2 ^{OH} 3	OF N CH ₃	4	Бн 5	6
			Ū			

Waters, Oasis, and Symmetry are trademarks of Waters Corporation. © Waters Corporation 1997

Tricyclic Antidepressants



A simple method, without amine modifiers, achieves high recoveries.

- Step 1: Condition Oasis[™] HLB Cartridge (1cc, 30 mg) with 1 mL MeOH
- Step 2: Equilibrate with 1 mL water
- Step 3: Load 1 mL acidified (20 μL H₃PO₄) & spiked porcine serum
- Step 4: Wash with 1 mL MeOH/water, 5/95 v/v
- Step 5: Elute cartridge with 1 mL MeOH
- Step 6: Evaporate; reconstitute in 200 μL 20/80 MeOH/buffer; inject 20 μL

Tricyclic Antidepressants




Recovery Data Comparison: HLB vs. C₁₈

		% Recovery (% RSD)	
		Oasis™ HLB	C ₁₈
Compound	Concentration	Cartridge	Cartridge
	(in porcine serum)	1cc 30 mg (n=5)	1cc 100 mg (n=3)
Nortriptyline	0.5 μg/mL	99.7 (2.3)	22.0 (2.3)
Doxepin	0.5 μg/mL	94.0 (1.3)	36.6 (3.9)
Amitriptyline	1.0 µg/mL	102.0 (2.5)	29.8 (8.4)

Tricyclic Antidepressants




HPLC Analysis: Peak 1: Nordoxepin [I.S.] Column: Waters Symmetry[®] C₁₈ Column, 3.9 x 150 mm CHCH₂CH₂NHCH₃ Mobile Phase: 20 mM phosphate, pH 7/MeOH (30/70 v/v) Temp.: 35° C.; UV Detection: 254 nm; Flow rate: 1.0 mL/min Upper Curve - Sample; Lower Curve - Serum Blank Peak 2: Nortriptyline CHCH₂CH₂NHCH₃ 0.001 AU 3 Peak 3: Doxepin $CHCH_2CH_2N(CH_3)_2$ Sample Peak 4: Amitriptyline Serum Blank $CHCH_2CH_2N(CH_3)_2$ 15 5 10 **Minutes** 20



column.	
Sample:	20 µL of reconstituted porcine serum extract
Mobile Phase:	1.2% glacial acetic acid/ACN (50/50 v/v)
Flow Rate:	0.6 mĽ/min
Temperature:	25° C.
UV Detection:	254 nm

SPE Spiked Serum Results					
Compound [µg/mL] %Recovery %RSD (n=6					
1	Naproxen	0.100 0.500 2.000	90.% 98.7 99.2	4.0% 1.9 0.85	
2	Ibuprofen	0.500 2.50	86.8 95.0	4.6 3.7	



HPLC Analysis Conditions:

Column:	SymmetryShield™ RP8, 3 x 150 mm
Sample:	20 μ L of reconstituted porcine serum extract
Mobile Phase:	20 mM ammonium acetate, pH 5/ACN, 95/5 (v/v)
Flow Rate:	0.6 mL/min
Temperature:	
UV Detection:	273 nm

SPE Spiked Serum Results					
Compound	[µ g/mL]	%Recovery	%RSD (n=6)		
1 Theobrom	nine <u>0.100</u> 0.500	103.%	3.4%		
	0.500	104.	2.4		
2 Paraxanth		<u> </u>	6.2		
	0.500	90.9	2.3		
3 Theophyll		100.	5.2		
	0.500	101.	2.1		
5 Caffeine	0.100	92.6	4.4		
	0.500	97.5	3.2		





A new polymer has been made with Hydrophilic-Lipophilic Balance to deliver unique benefits for SPE:

- It retains compounds with a wide range of polarities, even if the sorbent runs dry, for high, reproducible recoveries.
- A single, simple, rugged method for acidic, basic, and neutral analytes streamlines methods development dramatically.





- No longer are tedious, cumbersome, stopcock manipulations required to prevent cartridge drying. This makes SPE rapid & easy to perform, even on multiple vacuum manifolds operated in parallel for high sample throughput.
- Oasis[™] HLB sorbent is chemically & mechanically stable, extraordinarily clean, free of fines & silanol interactions for *predictable, consistent performance,* making *methods rugged and reliable*.

Want More Information?



For more information on Oasis™ HLB Extraction Cartridges, please contact your nearest Waters office.

Note added in proof: Oasis[™] HLB sorbent is now available in 96-well plate format (see next page).

Please visit our Internet site at (click on URL): <http://www.waters.com/oasis/welcome.htm> Try the Oasis™ Cartridge Challenge!

Sales Offices: Austria and Export (Central Europe, CIS, Middle East, India and Subcontinent) (43) 1 8771807, Australia (61) 2 9933 1777, Belgium and Luxembourg (32) 02 7261000, Brazil (55) 11 543 7788, Canada 800 252 4752, CIS 7 095 9 336 7000, Czech Republic (42) 2 472 8672, Denmark (45) 46 598080. Finland (358) 90 506 4140, France (33) 1 30487200, Germany (49) 06196 40 06 00, Hong Kong (852) 2964 1800, Hungary (36) 1 270 5086, India (91) 80 839 4799, Italy (39) 02 2500565, Japan (81) 3 3471 7191, Malaysia (603) 704 8600, Mexico (525) 524 7636, The Netherlands (31) 76 508 7200, Norway (47) 638 46 050, Peoples Republic of China - Beijing (86) 10 6506 0277, Poland (48) 22 33 4400, Puerto Rico (787) 747 8445, 790 2225, Singapore (65) 225 4855, Spain (34) 3 325 9616, Sweden (46) 8623 0090, Switzerland (41) 62 889 2030, Taiwan (886) 2 706 8766, UK and Ireland (44) 1 923 816700, All other countries: Waters Corporation U.S.A. 508 478 2000/800 252 4752

Oasis[™] HLB Sorbent Enables High Performance RP-SPE in 96-Well Plates²⁴

