## High Throughput LC/MS Quantitation for Drug Mixture

Kate Yu, Michael Balogh, Jeanne Li Waters Corporation, 34 Maple Street, Milford, MA 01757



## Waters Xterra<sup>®</sup> and Symmetry<sup>®</sup> Columns



# Introduction

In recent years, substantial effort have been made to increase the number of lead compounds produced in drug discovery. Because of the need of rapid turn around for discovery projects, the time allowed for method development is sharply decreased. And there is an ever increasing need to analyze more samples faster and cheaper.

Among the diverse analytical techniques available today, mass spectrometry seems to be the most versatile tool to cope with the demands of high throughput analysis. And LC/MS is the preferred method of choice in most cases.

In this paper, we described a a fast LC/MS method applied to a six drug compound model. The fast LC method was converted from a conventional LC method. By applying fundamental concepts of fast LC, we were able to reduce the cycle time from more than 20 minutes to 2.5 minutes. The analytes were detected by a single quadruple MS detector in electrospray mode by flow splitting. By applying"on-the-fly" positive and negative polarity switching, all six analytes in the mixture were analyzed within one injection. No carryover was observed. The quantitation limits of this "real life" assay model was studied so that the LOD, LOQ, linearity, precision and accuracy for each analyte were defined. The LC/MS quantitation procedure (including MS optimization, data acquisition, quantification and data reporting) can be automated by applying the QuanLynx<sup>TM</sup> (Part of MassLynx<sup>TM</sup> software).

## **Major Features**

#### Convert a conventional LC method to a fast LC method

- apply fast LC fundamentals
- utilize Waters Alliance<sup>®</sup> HT system an integrated solvent delivery manager designed for high throughput applications in LC and LC/MS.
- Reduce the cycle time for analysis to 2.5 minutes from more than 20 minutes
- Analysis carried on shorter Waters Xterra<sup>®</sup> C<sub>18</sub> columns (packed with 2.6 μm particles)

### • Electrospray MS detection on Waters ZMD<sup>TM</sup> 2000 MS Detector:

- "on-the-fly" positive/negative polarity switching allowed all compounds to be analyzed within one injection
- LC flow split to accommodate low flow requirement of electrospray to maintain or increase sensitivity

### • Quantitation performance evaluated:

- for each analyte, the LOD, LOQ, linearity, inter- and intra-assay precision and accuracy were defined
- No carryovers
- Automation by QuanLynx<sup>TM</sup>
  <sup>© Waters Corp. 2000</sup>

## Model Analyte

© Waters Corp. 2000



Bethamethasone MW 392 [378-44-9]



 $CH_2CH_2N(CH_3)_2$ 

Amitriptyline MW 277 [549-18-8]



Naproxen MW 230 [22204-53-1]



Prednisolone MW 360 [50-24-8]



Diphenhydramine MW 255 [147-24-0]



Ibuprofen MW 206 [15687-27-1]



# How To Reduce Analysis Cycle Time

# • In HPLC analysis, the time required for analyzing a sample is the cycle time

Cycle time is the time that required from injection to injection.

- sample withdraw, loading and injection time
- gradient time
- delay time
- equilibration time



- Shorten LC run time
  - Use shorter columns with smaller I.d. and smaller particle size
  - Increase flow rate
  - Apply higher column temperature
- Reduce Re-equilibration Time
  - Time needed for system equilibration: 3 x system volume + 5 x column volume
  - System equilibration can be accomplished in 0.2 minutes by switching the column off-line and apply high flow rate



- Apply pre-column volume concept
  - pre-column volume is a programmable software parameter in both Alliance<sup>®</sup> and Alliance<sup>®</sup> HT systems
  - it allows the injector to be loaded and gradient to be started, but sample is only injected onto to the column until the gradient reaches the column
  - the effect of pre-column volume is that the delay volume is at the end of the cycle, instead of beginning
  - the effect of pre-column volume is to compress the beginning of the chromatogram



### Make processes parallel

This can be accomplished easily with Waters Alliance HT solvent delivery system



# **Experimental Conditions**

## Conventional LC

### • HPLC: Waters Alliance<sup>®</sup> 2690

- Column: Waters Symmetry<sup>®</sup>  $C_{18}$  2.1 x 50 mm 3.5  $\mu$ m
- Column T: 30°C
- Mobile Phase: A: 10 mM NH4OAc in Water pH 5.0
  - B: 10 mM NH4OAc in AcN
- Flow Rate: 0.3 mL/min
- Injection Volume: 10 μL
- Gradient:

Time	<b>A%</b>	<b>B%</b>	Curve
0 min	75	25	
8 min	25	75	8
12 min	0	100	1
20 min	75	25	1

### MS: Waters ZMD<sup>®</sup> 2000 MS Detector

Ionization:	ESI+ or ESI-
Source T:	150°C
Desolvation T:	300°C
Capillary Voltage:	2.85 kv
Data Acquiring:	SIR (Single Ion Monitoring)

1	Scan Functions - S	IR.MDB								-		
Eile	Function Analog	<u>H</u> elp										
	olvent Delay ime (mins) 0.0 Pcl Probe Temp 20	0 Ar	halog Char	nnels— sed	2 🗖 Ur	nused	3 🗂 U	nused	4 🗖 U	nused		
Mir	\$ <mark>0 1 2</mark>	3	4	5	6	7	8	ģ	10	11	12	
1	SIR of 361.30 , Co	ne Voltage	13 V , (ES	6P+) ◀								-Prednisolone
2	SIR of 256.30 , Co	ne <mark>Voltag</mark> e	15 V , ( <mark>E</mark> S	6P+) 🗲								-Diphenhydramine
3	SIR of 393.30 , Co	ne Voltage	16 V <mark>, (E</mark> S	SP+)							_	-Betamethasone
4	SIR of 229.20 , Co	ne Voltage	16 V , (ES	6P-)			◀					-Naproxen
5	SIR of 278.30 , Co	ne Voltage	29 V , (ES	6P+)				● ←			-	-Amitriptyline
6	SIR of 205.20 , Co	ne Voltage	18 V , (ES	6P-)							•	—Ibuprofen

### **Extracted Mass Chromatograms**



## Fast LC

### • HPLC: Waters Alliance<sup>®</sup> HT 2790

- Column: Waters Xterra<sup>®</sup> 2.1 x 20 mm 2.6  $\mu$ m
- Column T: 40°C
- Mobile Phase: A: 10 mM NH4OAc in Water pH 5.0
  - B: 10 mM NH4OAc in AcN

- Gradient:

Time	A%	<b>B%</b>	Curve		
0 min	75	25			
0.7 min	55	45	6		
1.4 min	0	100	1		

- Pre-Column Vol: 100µL
- Flow Rate: 1.0 mL/min (Split flow, 150uL/min into MS)
- Injection Volume: 5 μL Full loop (4 times over fill)
- Fast Equilibration: System 5.0 ml/min 0.20 min
  Column 1.0 ml/min 0.35 min
  Wash Parameters: Injector Port: 10 s
- Wash Parameters: Injector Port:

Needle Exterior:

10 s

#### MS: ZMD<sup>®</sup> 2000 MS Detector:

ESI+ or ESI-
130°C
250°C
2.85 kv

🔪 Scan Functions	- DRUG	MIXSIR2.	MDB					_ 🗆 X	
File Function Anal	log <u>H</u> elp								
Full Scan SIR									
Solvent Delay	0.00	-Analog Cł	hannels —						
Desolvation Temp	20	1 🗖 U	nused	2 🗖 Unu	ised 3	🗖 Unuse	d 4 🗖	Unused	
L									
Mins Ó						i			N/ 000 001 050 000 070
1 SIR of 5 masse	s (ESP+)								M/Z 200, 361, 256, 393, 278
2 SIR of 3 masse	s (ESP-)								M/z 200, 229, 205
N//7	200	361	256	303	278	200	220	205	
	200	301	200	595	210	200	223	203	
Cone	+11	+16	+17	+17	+37	-11	-17	-16	
Voltage (V	()								
	'								

### **Extracted Mass Chromatograms**



## How Fast Should You Go

- How fast you should go will depend on the purpose of the analysis. Shorter run time does not necessarily results into shorter cycle time.
- In reality, there is a practical limit on how fast you can go. Faster is not necessarily better for your specific analysis. This is especially true for quantitative work.
- For quantitation, appropriate wash cycle is necessary to ensure that carryover is minimal from injector needles, injector valves and the LC columns. This will prolong the cycle time of the analysis.
- It also takes time for samples to be drawn and loaded. Water Alliance<sup>®</sup> HT allow this procedure to be hidden behind the gradient in parallel mode provided the gradient time is long enough.

# Quantitation with Fast LC

- The calibration curves shown below were constructed in the following way:
  - They are the results of 6 injections in 3 consecutive days.
  - Each day, two injections were made. Each injection was made from a separate vial.
  - Two sets of stock solutions were prepared and they were alternated between points.
- The limit of detection (LOD) and the limit of quantitation (LOQ) of each analyte was determined by measuring the signal to noise ratio (S/N).

- LOD: S/N = 3; LOQ: S/N = 10

- The QC samples were injected 6 times each day and repeated in 3 consecutive days.
  - Intra-assay CV (%):
  - Iner-assay CV (%)
  - Iner-assay Error (%)

within day precision, average of 3 days results

between day precision, take all 18 injections as one group of data

between day accuracy, take all 18 injections as one group of data

### Calibration Curves Naproxen, ESI-

#### Amitriptyline, ESI+

Correlation coefficient:  $r^2 = 0.998737$ Cal. curve: 22.9011 \* x + 1656.93 Cal. curve: 508.205 \* x + -104334 Response type: External Std, Area Response type: External Std, Area Curve type: Linear, Origin: Exclude, Curve type: Linear, Origin: Exclude, Weighting: Null, Axis trans: None Weighting: Null, Axis trans: None 200 400 200 0.0 0.0 400 N = 6N = 6pg/ul − pg/ul 10000.0 0.0 5000.0 5000.0 10000.0 0.0 © Waters Corp. 2000

Correlation coefficient:  $r^2 = 0.997367$ 

## **Quantitation Results**

	LOD (pg)	LOQ (pg)	R <sup>2</sup>	Intra- Assay CV (%) N=6	Inter- assay CV (%) N=18	Inter- assay Error (%) N=18
Pred.	55.2	184	0.992	3.77	4.40	1.62
Diphen	13.4	46.0	0.998	1.27	4.57	-8.00
Betha.	62.5	208	0.993	3.30	4.48	-0.00
Amitrip	8.42	28.1	0.999	2.06	2.63	-7.29
Napro.	46.1	154	0.997	2.72	4.60	14.5
lbu.	6.43	21.4	0.996	1.56	4.63	-0.66

All precision and accuracy data shown here were obtained at low QC level (3 x LOQ).

The LOD and LOQ was calculated based on total injection amount to the column. The LC flow was split about 1:7 ratio before the MS.

## The Effect of Carryover



No carryover was observed with the system and the method used for the project. Shown here are the example of one set of experiments.

Samples were running in the following sequence: Blank, Series of standards from low concentration to high concentration (1ng/ml to 10 ug/ml), and Blank.

Shown here are the comparison of the lowest signal detectable for each analyte with the blank injection (injected after the highest concentration of sample).

### Quantitation Automated by QuanLynx<sup>™</sup>

- The LC/MS quantitation procedure can be automated by utilizing QuanLynx<sup>™</sup> software from the MassLynx<sup>™</sup> package
- QuanLynx<sup>™</sup> allows MS condition optimization, data acquisition, data processing to be performed in one step

QuanLynx Method Editor Caffeine.QLM	$\times$
<u>F</u> ile	
Optimisation Acquisition Tune Files Adducts Losses Inlet Methods	_
Cone VoltageIonisation Modeminmaxstep55075507	
Collision Energy  MS Method Creation    min  max  step    10  50  6	
Optimisation Peak Detection Parameters — Fragment Size	
MS Quan Method: MS Quan Bro <u>w</u> se min: 50	
MSMS Quan Method: Browse	

# Conclusions

- Utilizing the Waters Alliance<sup>®</sup> HT 2790 solvent delivery system, by applying and coordinating the fast LC principles, the cycle time for the drug mixture was reduced to 2.5 minutes from more than 20 minutes.
- Quantitation results were obtained with no carryover.
- At low QC level, the intra-assay precision was less than 3.8%, the inter-assay precision was less than 4.7% and the inter-assay error was less than 14.5% for the whole methods.
- The LC/MS quantitation process (optimization, data collection, and calculation) can be automated with QuanLynx<sup>™</sup>.

## Acknowledgment

• Bonnie Alden from Waters Corporation for providing the Xterra columns