



# LC/MS Application Notes: EI/APCI Identification and Quantitation of Polymer Additives

Kate Yu, Waters Corporation, Milford, MA



[Instrumentation](#)

[Chromatogram](#)

[MS Spectrum](#)

[Calibration Curve](#)

[Quantitation](#)

## **Key Words:**

Polypropylene

Additives

EI/APCI

The analysis of polymer additives has always been important to some industries (polymer, food packaging, medical device, chemical, inks/coatings). However, the identification and quantitation of polymer additives has been rather elusive. Reports of analyzing polymer additives by single quadrupole LC/MS are scarce. This work demonstrated a how quantitative polymer additive analysis can be accomplished by the two complementing LC/MS techniques: EI and APCI.

# Instrumental Conditions

## HPLC

- Waters Alliance Solvent Delivery System
- Symmetry C8 (3.9 x 150 mm) with Acetonitrile/Water gradient
- 0.4 mL/min. flow rate with 10 uL injection

## EI

- Waters Integrity  
ThermaBeam™ System
  - Library search provides  
lead suspects


## APCI

Micromass PlatformLC System

Positive APCI

- Full scan provides [MW information](#)
- Selected Ion Recording (SIR)  
provides [quantitation](#)


# Figure 1: Positive Sample Identification




## Quantitation

- The quantitation was obtained by positive APCI SIR acquisition
- There were four standards commercially available
- Four calibration curves were obtained with a representative curve shown on [Figure 3](#)

|              | LOD<br>(ppb) | Linear Range<br>(ppb) | % (w/w) in Polymer<br>(ppb) |
|--------------|--------------|-----------------------|-----------------------------|
| NC-4         | 30.4         | 30.4 - 7600           | 0.07                        |
| Naugard-XL   | 14.4         | 14.4 - 3600           | 0.04                        |
| Octadecanol  | 27.2         | 27.2 - 6800           | 0.4                         |
| Irgonox 1076 | 12.8         | 12.8 - 3200           | 0.06                        |



**Figure 2: Identification of Peak #5**



**Figure 3: Calibration Curve of Naugard-XL by APCI+**