

Waters

Lab Highlights


PAH

81.600.030.006.005
.038

POLYNUCLEAR AROMATIC HYDROCARBON ANALYSIS (PAH) USING THE RCSS

The continuing interest in the polynuclear aromatic hydrocarbons (PAH's) because of their carcinogenic and mutagenic properties has initiated research into novel separations techniques. Their trace level analysis from air, water and marine sediments has necessitated utilization of several chromatographic techniques, as well as optimized combinations of column, detectors and extraction methods.

A tailored cartridge and Waters Radial Compression (RCSS) technology, combined with a highly sensitive fluorescence detector (M420), allows an isocratic separation of six major carcinogenic pollutants in about eight minutes. Chromatographic conditions are shown below:



Column: Radial-PAK PAH (8mm x 10cm)

Solvent: H_2O/CH_3CN (28/72)

Flow Rate: 4.0 ml/min.

Detector: Waters M420 (254nm excitation,
375nm emission) 128X

Sample Volume: 10 μ l

Bill Dark
Market Development Lab



Waters Associates Milford, Massachusetts 01757 / (617) 478-2000