


DETERMINATION OF VITAMIN D₂ AND D₃ USING RCSS

The determination of Vitamin D and the resolution of its various forms is of particular importance in nutritional studies. Vitamin D is added routinely to many foods and animal feeds and its determination has long been the object of analytical techniques.

A recent publication (1) by Thompson, Maxwell, et al describes an improved method for determination of both Vitamin D₂ and D₃ from milk and milk products and presents a baseline separation of Vitamin D₂ and D₃ in 15 minutes, using the RCSS. Previously, this separation was not accomplished using many other commercial C₁₈ columns.

The chromatogram and conditions are shown below



Column: Radial PAK C₁₈ (8mm X 10cm, 5μ)
Solvent: Methanol
Flow Rate: 1 ml/min.
Detector: Waters M440, 254nm
Sample Volume: 50μl

William Maxwell 9/22/81
Waters Scientific Ltd. (Canada)

- (1) J. N. Thompson, G. Hatina, W. Maxwell, S. Duval, Determination of Vitamin D in Fortified Milks, Margarine, and Instant Formulas by HPLC, Health and Welfare Canada, Tunneys Pasture, Ottawa, Ont. in press