

Waters

# Lab Highlights

LAH 0094 2/83  
Doc # M1031  
AN/FA,/OC,MD/FF/BT

## QUALITY CONTROL FOR CAFFEINE-FREE SOFT DRINKS




Growing public concern over caffeine intake has prompted the soft drink industry to introduce a number of caffeine-free products. A recent article in the Wall Street Journal (11/11/82) addressed this topic by claiming there will soon be a dozen caffeine-free brands available which will capture between 8% and 10% of the soft drink market. Since many major soft drink producers offer franchises to independent bottlers as a means of achieving nationwide distribution, the quality control and assurance of caffeine-free production are of major concern. The QA-1<sup>™</sup> Quality Analyzer and the proper chemistry can provide a rapid and reliable determination of caffeine levels in beverages.

Figure 1 shows the separation of common additives in a soft drink. Peak 1 represents 0.03% (w/v) sodium saccharin, Peak 2 is the response from 0.025% caffeine and Peak 3 is 0.03% benzoic acid. The remaining chromatograms were generated with samples of a brand name soft drink which were produced by an independent bottler with a franchise from the parent company and were purchased directly off the store shelf.

Figure 2 represents a diluted (7/10 v/v) sugar-sweetened, caffeine-containing soft drink. Caffeine is indeed present (Peak 2) at 0.01%. Figure 3 was generated from an equally diluted sample of supposedly sugar-free and caffeine-free soft drink. As expected, saccharin (Peak 1) is present at sweetening levels (0.048%). Surprisingly, however, caffeine in a caffeine-free product is found at the same level (0.01%) as a caffeine-containing product. This represents an apparent error in the formulation of this product.

One must realize that a QA-1<sup>™</sup> Quality Analyzer in the bottler's facility with routine batch-to-batch analysis could have detected this apparent error and prevented the distribution of caffeine-free products which actually contained caffeine.

Peak 1 = Saccharin  
Peak 2 = Caffeine  
Peak 3 = Benzoic Acid



See reverse side for conditions.

John Morawski

Waters

Division of MILLIPORE 34 Maple St. / Milford, MA 01757 / 617-478-2000

CONDITIONS

Instrument: Waters QA-1™ Quality Analyzer  
Column:  $\mu$ BONDAPAK™ C<sub>18</sub> Radial-PAK™  
Mobile Phase: 10% CH<sub>3</sub>CN/  
90% 2.0% Ammonium Acetate with  
2.5% Acetic Acid  
Elution Volume: 30 ml  
Flow Rate: 4.0 ml/min  
Detector: UV 254 nm @ 0.1 AUFS  
Injection Volume: 10  $\mu$ l  
Sample: Soft Drink diluted 7/10 (v/v) with eluent