

Waters

Lab Highlights

LAH 11
POLYMER-STYRENE
81.600.104.004.011
(.078)

ULTRASTYRAGEL (100Å) VS. SHODEX (A801)

ULTRASTYRAGEL WINS!

Polystyrene 450 is a mixture of low molecular weight oligomers that can be conveniently used to compare the resolving power of small pore size GPC columns such as 100Å Ultrastyrigel and Shodex A801. The resolution of PS450 with two 30 cm 100Å Ultrastyrigel columns (total length 60 cm) is far superior to that of one 50 cm Shodex A801. This improvement is much greater than that which can be ascribed solely to the extra length of Ultrastyrigel.

- All valleys between adjacent peaks (A,B,C,D) are much deeper with Ultrastyrigel
- Peak E is clearly resolved with Ultrastyrigel but is barely discernible as a shoulder with Shodex.
- Peak F is partially resolved with Ultrastyrigel but not at all with Shodex
- 100Å Ultrastyrigel is capable of resolving higher molecular weight species than Shodex A801 since there is relatively more blank baseline between the exclusion volume (V_e) and the start of elution (G) of the polystyrene sample in Figure #1.

Figure #1





Figure #2



Juris Ekmanis

9/23/81

