

**Waters**


# Lab Highlights

LAH 0117 6/83  
Doc # M0355, M0341  
AN/HO,ES/PM,QC/IC/AX

## THE CHROMATOGRAPHY OF UV-ABSORBING INORGANIC IONS

A previous Lab Highlight on nitrate analysis indicated the potential for separating and detecting UV-absorbing anions using Waters' equipment. This Highlight shows a separation of the most common UV-absorbing anions. Acetate has been included in the ion mixture. Other organic acids can also be analyzed under these conditions. Some matrices (e.g., meat pickling brines) contain weak organic acids as well as inorganic ions. It is possible that such samples could be analyzed using a single set of conditions.

### INORGANIC IONS ABSORBING AT 214 nm



Column:  $\mu$  BONDAPAK™ C<sub>18</sub>

Mobile Phase: 2L of 0.01M KH<sub>2</sub>PO<sub>4</sub>  
and 1 vial of low UV  
PIC® Reagent A, pH = 6.0

Flow Rate: 1.5 ml/min.


Detector: Waters Model 441 UV/VIS Detector,  
214 nm, 0.02 AUFS

1. Iodate
2. Acetate
3. Bromate
4. Nitrite
5. Bromide
6. Azide
7. Nitrate
8. Chromate
9. Thiosulfate
10. Iodide

400 ng each; except Acetate  
(1200ng)

The obvious limitation of UV detection is that many anions (e.g., chloride, sulfate) do not absorb at 214 nm. This selectivity of UV detection can be used to advantage. In the following example, the iodide content of table salt is measured in about four minutes. Using conductivity detection, the chloride peak would overwhelm the iodide. At 214 nm, only the iodide is detected. Another example of this beneficial selectivity is in the analysis of nitrate, bromide, or iodide in sea water.

### IODIDE IN TABLE SALT



Joan Newburger  
Roy Day