

A RAPID SCREENING PROCEDURE USING RADIAL-PAK™ RESOLVE C₁₈ FOR COMMON ANTICONVULSANTS IN SERUM

A method for the fast analysis of three common anticonvulsants in serum has been developed using a Radial-PAK™ RESOLVE™ C₁₈ column (5 mm I.D. X 100 mm) and C₁₈ Guard-PAK™ cartridge. The procedure, developed by Udo Rupprecht of Millipore, Australia, has the advantage of rapid sample preparation along with a short (under 6 minutes) analysis time. More than ten laboratories in Australia use this as a screening procedure for anticonvulsants.

The procedure monitors amounts of phenobarbital, phenytoin and carbamazepine along with an internal standard, 5-(4 methyl phenyl)-5-phenyl hydantoin (5-MPPH). The 5-MPPH is dissolved in acetonitrile at a concentration of 0.05 mg/ml. Sample preparation requires taking 200 µl of the internal standard (5-MPPH) and adding 200 µl of serum. The sample is then vortexed for 10 seconds, centrifuged for 10 minutes and 20 µl of the supernatant is injected. A typical chromatogram obtained from this procedure is shown in Figure 1.

The mobile phase consists of 300 ml of methanol, 180 ml of acetonitrile, and 0.88 g of KH₂PO₄ in 650 ml of water (pH = 6.0 with H₃PO₄ or KOH) at a flow rate of 3.0 ml/min. Detection of the peaks is monitored at 214 nm, 0.1 AUFS.




FIGURE 1: Chromatogram of three anticonvulsants in serum plus internal standard. Sample concentrations are (1) phenobarbital 100 µmoles/L, (2) phenytoin 65 µmoles/L, (3) carbamazepine 40 µmoles/L, (4) 5-MPPH, internal standard. Instrument: Model 320 QA-1™ Analyzer. Other chromatographic conditions are in the text.