

Waters


Lab Highlights

ANTIBIOTICS

81.600.014.001.025

GENTAMICIN FROM FERMENTATION BROTH

Gentamicin is one of several aminoglycoside antibiotics that have been clearly established to have major value in the management of aerobic Gram-negative bacillary infections. Current production procedures for gentamicin involve fermentation of *Micromonospora purpurea* or *M. echinospora* and variants thereof. HPLC methods for the separation and post-column o-phthalaldehyde detection of the three major components of gentamicin - C_1 , C_{1A} , and C_2 - have been developed and applied to the measurement of serum levels of the antibiotic by J. P. Anhalt. Interfering serum compounds were removed by ion exchange gel chromatography utilizing CM-Sephadex (C-25). The producers of the antibiotic must also monitor the three gentamicin components, but must do so from the fermentation broth sample matrix. The application of the CM-Sephadex cleanup procedure is demonstrated on fermentation broth samples submitted by ANSA Antibiotics Co.



Standard

Column: μ Bondapak C_{18}

Column

Eluent: 3% methanol
97% 0.2M Na_2SO_4 /
0.02M pentane sulfonic acid/
0.1% (V/V) acetic acid

Flow at 1 ml/min

Derivatizing Reagent:


- 1) 15g H_3BO_3 dissolved in 500 ml H_2O
pH adjust to 10.4 with KOH pellets
- 2) 300 mg o-phthalaldehyde dissolved
in 5 ml ethanol
+ 0.1 ml β -mercaptoethanol

Mix solutions 1 and 2


(1) J. P. Anhalt, Antimicrobial Agents and Chemotherapy 11 (4), 651-655 (1977)

Robert Burgoyne 12/18/81





Fermentation Broth



Extraction Product
of Fermentation Broth

POST COLUMN OPA

