

Lab Highlights

THE MONITORING OF AN ENZYMATIC REACTION BY LC AND RADIAL COMPRESSION

A simple, quick and efficient analytical method has been reported by Ladron DeGuevara *et al.* (1) for the monitoring of the enzymatic conversion of DL-phenylhydantoin to D-phenylglycine.

The D-phenylglycine produced by this patented enzymatic process is an important starting material in the production of β -lactam antibiotics such as semi-synthetic penicillins and cephalosporins.

Figure 1

Column: Radial-PAK™ μ BONDAPAK™ C₁₈ (8mm X 10cm)
Holder: Z-Module™ Radial Compression Module
Detector: 212nm with M490 Programmable Multiwavelength Detector
Mobile Phase: Gradient Elution (Curve No. 5, M680 Gradient Controller)
20 to 25% methanol in 20mM sodium phosphate buffer,
pH 6.8 over 9 minutes.
Flow Rate: 1.5 ml/min.




Figure 1 shows a 10 μ l injection of a standard mixture (2.5 μ g of D-phenylglycine plus 5 μ g of DL-phenylhydantoin). Figures 2 and 3 show the results obtained at the start and at the end of the enzymatic reaction, respectively.

Figure 2





Figure 3



This quick and efficient analytical method is now the standard method used at the Centro de Ingenieria Genetica y Biotecnologia, UNAM, for the monitoring of the enzymatic conversion of DL-phenylhydantoin into D-phenylglycine.

1. Ladrón De Guevara, O., Quintero, R., and Padilla, P., J. Chromatogr., 329 (1985) 428-433.