


ANALYSIS OF PICOMOLE QUANTITIES OF ADENOSINE


Adenosine, an important nucleoside in physiological systems, influences coronary and cerebral circulation, cyclic AMP regulation, neurotransmission, platelet aggregation, and immunodeficiency disease <sup>1</sup>. The analysis of physiological levels of adenosine compounds has been hampered by the necessarily low detection limits required. Detection is enhanced by first derivatizing adenosine with chloroacetaldehyde according to the reaction below to form a highly fluorescent compound which enables detection limits to be attained at the

Chloroacetaldehyde

Typical reaction:



2 picomole range. Although the excitation maximum for the adenosine derivative is approximately 275 nm, very sensitive detection is also obtained by utilizing the strong 254 nm excitation line of the mercury lamp and detection at 400 nm.

Column:  $\mu$ BONDAPAK<sup>TM</sup> C<sub>18</sub> (60 cm X 3.9 mm)Eluent: CH<sub>3</sub>OH/50 mM KH<sub>2</sub>PO<sub>4</sub> (25:75)  
pH adjusted to 7.5

Flow Rate: 1.0 ml/min.

Detector: Model 420  
Hg lamp  
Ex: 254 nm band-pass  
Em: 400 nm long-pass  
Gain: 128X

Robert Burgoyne

1. J. F. Kuttesch, et al., J. Liquid Chromatogr. 1(1), 97-109 (1978).