

**Waters**

# Lab Highlights

LAH 0243 5/85  
AN/LS/MD/PR/ST


## SEPARATION OF PROTEINS AND PEPTIDES ON $\mu$ BONDAPAK<sup>TM</sup> C<sub>18</sub>

### PART II

The  $\mu$ BONDAPAK<sup>TM</sup> C<sub>18</sub> column is widely recognized as having versatile selectivity and excellent batch-to-batch reproducibility. It is also used extensively in reversed-phase protein and peptide separations, most likely because it has a proportion of wide pores and has the correct "chemistry" of a bonded phase (1).  $\mu$ BONDAPAK<sup>TM</sup> C<sub>18</sub> column packing is still the most used packing material as evidenced by the large number of literature citations (2). In order to support these claims, a number of Lab Highlights have been assembled (see LAH 0242 and LAH 0244) to illustrate the use of  $\mu$ BONDAPAK<sup>TM</sup> C<sub>18</sub> columns for protein separations.


A few of the most impressive separations along with chromatographic conditions are shown below.

FIGURE 1: Oligopeptides on  $\mu$ BONDAPAK<sup>TM</sup> C<sub>18</sub> (Ref.3).




HPLC separation of mixture of seven oligopeptides with TEAF buffer: 500 ng each of bradykinin (B), angiotensin II (A), methionine-enkephalin (ME), eledoisin-related peptide (E) and leucine-enkephalin (LE); 1  $\mu$ g each of substance P (P) and somatostatin (SS). 26% CH<sub>3</sub>CN; 74% TEAF; 200 nm; 0.1 AUFS; 1.5 ml min<sup>-1</sup>; 1200 psi. TEAF is prepared by titrating 0.04 M formic acid with triethylamine to pH = 3.15

FIGURE 2: HPLC of a thymic peptide (used in the maintenance of immune balance and also growth, development and function of the lymphoid system.) (Ref.4).




HPLC analysis of synthetic thymosin  $\alpha_1$  (bovine) in a 30  $\times$  0.39-cm  $\mu$ Bondapak C<sub>18</sub> column (10  $\mu$ m, Waters Assoc.) at 35°C. The solvents used were: 0.05% TFA in Reservoir A and acetonitrile containing 0.05% TFA in Reservoir B. The peptides were eluted with 10% B for 10 min, followed by a linear gradient from 10 to 20% in 10 min and second gradient from 20 to 45% B in 45 min. The flow-rate was set at 1.5 ml/min, with a chart speed of 0.3 cm/min. The elution was monitored by UV absorbance at 210 nm (—). The solvent gradient is indicated by a dashed line (---).

FIGURE 3: Separation of Diastereoisomeric peptides bombesin from [D-Met<sub>14</sub>] bombesin (Ref. 5).



Separation of bombesin (1) from [D-Met<sub>14</sub>]-bombesin (2) by HPLC. Conditions:  $\mu$ -Bondapak C<sub>18</sub> column (120 x 0.7 cm); flow rate 6 ml/min; isocratic elution using 30% acetonitrile-10 mM ammonium acetate, pH 4.2.


FIGURE 4: Separation of Polypeptides (Ref. 6).



Gradient elution separation of the polypeptides angiotensin II 1, bovine insulin A chain 2, angiotensin I 3, bovine insulin 4, bovine insulin B chain 5, and porcine glucagon 6.


Chromatographic conditions: column,  $\mu$ Bondapak C<sub>18</sub> (10  $\mu$ m, 30 x 0.4 cm i.d.); flow rate: 2.5 ml/min; temperature: 18°; mobile phase 90-min linear gradient from 15% acetonitrile-water-15mM orthophosphoric acid to 50% acetonitrile-50% water-15mM orthophosphoric acid; sample loading 2  $\mu$ g each polypeptide.

FIGURE 5: Thyrotropin-releasing hormone analogs. Mobile phase: 0.2 N AcOH, 0.1% C<sub>7</sub>SO<sub>3</sub>Na. Flow Rate = 2.0 ml/min; Temperature = 60°C (Ref. 7).



DA-TRH = Pyr-His-Pro, Phe-TRH = Pyr-Phe-Pro-NH<sub>2</sub>, Glu-TRH = Glu-His-Pro-NH<sub>2</sub>

FIGURE 6: Endorphin Fragments of Rat Pituitary Glands. Mobile phase: 10 mM NH<sub>4</sub>Ac, pH = 4.15, 30% CH<sub>3</sub>OH. Gradient: 30% to 75% MeOH over 45 min at 2.0 ml/min. Temperature at ambient. (Ref. 8).



Extracts of rat pituitary glands were chromatographed. Samples isolated by evaporation of methanol at 60° using a Buchler vortex evaporator and remaining material lyophilized. Sample volumes up to 1 ml did not effect resolution. Samples were found to be active in subsequent radioimmunoassays.

1. F. V. Warren, Jr. and B. A. Bidlingmeyer, Anal. Chem. **56**, 6 (1984).
2. Waters Lab Highlight LAH 0242.
3. D. M. Desiderio and M. D. Cunningham, J. Liq. Chrom., **4**, 721-733 (1981).
4. T. L. K. Low et al., J. Chromatogr., **266**, 533-544 (1983).
5. J. E. Rivier and M. R. Brown, Biochemistry, **17**, 1766 (1978).
6. B. Grego and M. T. W. Hearn, Chromatographia, **14**, 589-592 (1981).
7. E. Spindel and R. J. Wurtman, J. Chromatogr., **175**, 198 (1979).
8. J. G. Loeber et al., Biochem. Biophys. Res. Commun., **86**, 1288-1295 (1979).