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LC/MS Conditions

1D LC-SRM


System :ACQUITY UPLC® I-Class (Waters)
Column :ACQUITY UPLC BEH300 C18 , 2.1 x 150 mm, 1.7 um
Column temp.: 35°C
Flow rate :0.3 mL/min.
Gradient :(A) 0.1% (v/v) formic acid in water
(B) 0.1% (v/v) formic acid in acetonitrile
linear gradient from 0 to 35% B in 10min

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Time Range	Valve Position	Analyte undergoes
0 - 6.6 min	1	high pH RP separation
6.6 - 9.9 min	2	transfer from 1 D to 2D
6.9 - 14 min	1	low pH RP separation

Figure 2: Heart cut configuration for 2-dimensional chromatography: (A) Sample loading and first dimension separation under basic conditions (pH 10.0); (B) Analyte transfer from the first dimension to the second dimension; (C) Separation in the second dimension under acidic conditions at pH 2.5.

First Dimension separation (pH=10):
Column (P/N 186003118): XBridge C18, 1.0 x 50 mm, 2.5 μ m particles, operated at 100 μ L/min, kept at 35 °C
Mobile phases: 20 mM ammonium formate in water (Solvent A)
20 mM ammonium formate in 90% ACN (Solvent B)
Gradient profile: from 0 to 40% B over 5 min

Second Dimension separation (pH=2.5):
Column (P/N 186003685): BEH300, 2.1 x 50 mm, 1.7 μ m particles, Mobile phases: 0.1% formic acid in water (Eluent A)
0.1% FA in ACN (Eluent B)
Gradient profile: from 0 to 30% B in 5 min (starts 7 min after inj.).

RESULTS

Trypsin Digestion Optimization


Figure 3: Protein to trypsin ratio optimization. Five digestion ratios (1:10, 1:20, 1:30, 1:50 and 1:100) were investigated for 25 nM trastuzumab and 25 nM ¹³C¹⁵N-isotopically labeled peptides digested in human serum. The samples were digested overnight (16 h) with Sigma trypsin (cat no T-6567) at 37 °C.


Figure 4: Effect of incubation time on digestion efficiency. Trastuzumab (25 nM) and ¹³C¹⁵N-isotopically labeled peptides (25 nM) were digested in serum for 15 min, 30 min, 1h, 3h, 6h and 16 h (overnight) with Sigma trypsin (cat no T-6567) at 37 °C.


Figure 5: Reproducibility of trypsin digestion. Serum was spiked with 5 nM trastuzumab and 5 nM ¹³C¹⁵N-isotopically labeled peptides and digested overnight (16 h) with Sigma trypsin (T-6567) at 37 °C. The RSD % for the ¹²C/¹³C peptide ratio was better than 6%.

2D LC-SRM Assay


Figure 6. Linearity of the 2D-SRM assay. Four digest concentrations covering a dynamic range of 500 fold (0.1 to 50 nM trastuzumab) were prepared in 20 mM ammonium formate (pH 10) and analyzed in replicate (n=4). The RSD for this assay was better than 2%.


Figure 7. RADAR monitoring of SPE-cleaned samples. (A) 1D LC-SRM separation under acidic conditions (pH 2.5); (B) 2D LC-SRM separation with heart cutting. The sample was prepared by spiking a trastuzumab digest (5 nM trastuzumab and 10 nM ¹³C¹⁵N peptides) into SPE-cleaned human serum digest.


Figure 8. Reproducibility of the 2D LC-SRM assay in the presence of the serum digest. A sample containing 5 nM trastuzumab and 10 nM ¹³C¹⁵N peptides was spiked into SPE-cleaned human serum digest. Peak area RSD was better than 5.0%.


Figure 9: Evaluation of signal suppression. Comparison of average peak areas recorded for the endogenous and isotopically labeled trastuzumab peptide FTISADTSK in the absence/presence of the serum digest matrix. A sample containing 5 nM trastuzumab digest and 10 nM ¹³C¹⁵N peptides was spiked into SPE-cleaned human serum digest.

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