イオンモビリティ-飛行時間型質量分析計を用いた Top/Middle-Downフラグメンテーション解析によるタンパク質配列解析および酸化部位の同定

○川瀬泰司¹、廣瀬賢治¹、Stephane Houel²、Catalin Doneanu²、Asish Chakraborty²、Andrew Tudor²、Nick Tomczyk²、Martin Palmer²、Weibin Chen² 1日本ウォーターズ株式会社、2Waters Corporation, Milford, MA, United States

OVERVIEW

METHODS

- Fragmentation of large therapeutic proteins (top /mid down analysis) tends to generate many different types (and charge states) of fragment ions, making the spectra interpretation and sequence deduction very challenging.
- The biopharmaceutical Platform Solution with UNIEL provides a workflow for protein top/middle down analysis including targeted MSMS acquisition, data processing and automatic report generating.
- In this work, we demonstrated the top/middle down analysis using the developed UNIFI workflow to characterize Trastuzumab (Tmab) sub - units on a Xevo G2-XS quadrupole time-of-flight mass spectrometer
- In addition, Waters SYNAPT G2-Si HDMS provides a unique instrument platform for top down protein analysis to simplify the terminal sequence confirmation using ion mobility separation (IMS) of product ions

UNIFI Top Down Analysis Workflow UNIFI Platform Solution UNIFI Process Only Synapt G2-Si HDMS Data Acquisition geted MSMS CID or ETD Data Acquisition BayesSpray algorit Data Pr Option to export to ProSightPTM ProSightPTM 63 UNIFI Figure 1: The Biopharmaceutical Platform Solution with UNIFI provides a workflow for protein top/middle

down analysis including targeted MSMS acquisition, data processing and automatic report generating. UNIFI process only option supports MassLynx data with complementary CID and ETD fragmentation

Sample Preparation:

FaBricator (IdeS) was used to cleave Transtuzumab at the hinge region, followed by limited reduction and alkylation using DTT and iodoacetamide (IAA) to reduce the inter-chain disulfide bonds and generate the light chain and two heavy chain fragments. Neat acetic acid was added to terminate all digestions. The digests were diluted to 0.5 µg/µL prior to LC MS analysis.



Figure 2: Schematic of Ides digestion followed by reduction to generate Fc/2, Light Chain (LC), and Fd subunits from monoclonal antibodies.

LC-MS and LC-MS/MS Analysis:

LC: ACQUITY UPLC H-Class Bio System Column: ACOUITY UPLC Protein BEH C4 Column (300Å, 1.7 μm, 2.1 mm X 50 mm) (P/N=186004495) Mobile Phases:

A: 0.1% TFA in Waters; D: 0.1% TFA in Acetonitrile Time(min) How Rela(mL/min) %A %6 %C %D Curve

Initial	0.400	95.0	n.o	ഹ	5.0	Initial
1.00	0.400	95.0	0.0	n.o	5.0	6
1.10	0.200	85.0	0.0	8.0	5.0	6
2.50	0.206	76.0	aa	0.0	26.6	6
12.50	0.200	58.0	0.0	0.0	42.0	6
13.00	6.200	5.0	6.0	0.0	95.0	6
13.10	6.400	5.0	8.0	n.o	95.0	6
14.90	0.400	5.0	aa	ດມ	95.0	6
15.00	0.490	95-0	0.0	0.0	5.0	6
20.00	0.400	95.0	0.0	0.0	5.0	6

MS Conditions C

S

Capillary: 3kV;
Sample Cone voltage: 150 v;
Source Temp: 500°C;
Desolvation Temp: 350 °C;
Desolvation Gas Flow: 800 L/h
Scan Rate 0.5 Sec; Mass Range: 500-4000 m/z
Targeted MS/MS Conditions
Scan Rate 1 Sec; Mass Range: 100-4000 m/z
Isolation Window: 2 m/z
CE: 32 eV

Software: UNIFI Scientific Information System V1.8 and MassLvnx 4.1

une care aus de la Precursor Ions Selected for Fc/2 Targeted MSMS Fd MaxEnt1 Deconvolution

Subunit MS Analysis with UNIFI (Xevo G2-XS)

and 22 b- and 9 v-ions were assigned for Ed respectively.

corresponding to approximately 11% of the theoretical

number of fragments.

Terrar Bernerate Bernerate Bernerate Berner Berner Berner Bernerate Be Fd+2 Alk (10 ppm)

Et al anna

RESULTS AND DISCUSSION

Na A. Stream Tra-

Report

Figure 3: Trab subunits MS analysis followed by top down analysis using the Biopharmaceutical platform solution with UNIFI. MaxEnt1 deconvolution confirms the molecular weight of the subunits and their modified forms. The most abundant precursor ions were selected for CID fragmentation followed by data processing and reporting. Fc/2 top down analysis data and UNIFI report is shown here.

sequence coverage

CONCLUSION **Option to Export to ProsightPTM Analysis** LC Fd СОСТВОЕНИЕ С СОСТ The Biopharmaceutical Platform Solution with UNIFI provides a торамираа соружают в соружают списати зама соружают соружают за из сори зама соружают соружают соружают в соружают соружают соружают соружают соружают обратоваторая транаа соружают обратоваторая транаа соружают обратоваторая транаа соружают обратоваторая транаа соружают сору соружают соружают соружают сору соружают соружают streamlined workflow from mAb intact/subunit MS analysis and subsequent protein top/middle down MSMS analysis. 84419448 10164710868851941 The developed UNIFI workflow enables automate targeted MSMS acquisition, data processing and report generating. Figure 4: LC and Fd top down analysis data (BayesSpray Top down data can be exported to ProsightPTM analysis for deconvoluted fragment ions with intensity) were exported additional confirmation as .PUF file followed by ProsightPTM analysis. Matched sequence coverage of Tmab LC and Fd subunits are The post-fragmentation IMS separation enabled the separation of shown here. 10 b- and 22 y-ions were assigned for LC different charge states of fragment ions, produced easily

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

Waters THE SCIENCE OF WHAT'S POSSIBLE.®



© 2015 Waters Corporation MKT15083