# ANALYSIS OF PLASMA 17-HYDROXYPROGESTERONE, ANDROSTENEDIONE, AND CORTISOL USING A NOVEL SOLID-PHASE EXTRACTION SORBENT FOR UPLC-MS/MS ANALYSIS

Shunya Sasaki<sup>1</sup>, Jonathan P. Danaceau<sup>2</sup>, Erin E. Chambers<sup>2</sup> <sup>1</sup>Nihon Waters K.K., Tokyo, Japan, <sup>2</sup>Waters Corporation, Milford MA, USA

# INTRODUCTION

Bioanalysis of drugs and endogenous molecules from plasma samples often involves using solid phase extraction (SPE) to clean up the matrix and concentrate the analytes of interest. Despite the advances in SPE sorbents and formats, residual phospholipids often remain in extracts and can interfere with analyses by causing ion suppression and by prematurely fouling analytical columns and MS sources. Using a novel SPE sorbent that is designed to remove phospholipids from samples, a panel of corticosteroids has been extracted from plasma. This sorbent is also water wettable, enabling extraction without the usual requisite preconditioning and equilibration steps. This has resulted in a method with excellent sensitivity that demonstrates consistent recovery, minimal matrix effects, and the elimination of >95% of phospholipids compared to simple protein precipitation.

## **METHODS**

Calibration standards were purchased from Sigma-Aldrich (St. Louis, MO, USA). All deuterated internal standards were purchased from Cerilliant (Round Rock, TX). Double charcoal-stripped plasma was used for calibrators and QC samples.

Sample Preparation: 150  $\mu$ L plasma samples were precipitated with 300  $\mu$ L of a solution of 4:1 MeOH and 89 g/L ZnSO<sub>4</sub> (V:V). After centrifugation, 300  $\mu$ L of supernatant was diluted with 900  $\mu$ L 4% H<sub>3</sub>PO<sub>4</sub> and directly loaded onto Waters' Oasis PRiME  $\mu$ Elution SPE plates. All wells were then washed with 25% MeOH and eluted with 2 x 25  $\mu$ L of 90:10 ACN:MeOH. The sample eluates were diluted with 25  $\mu$ L water and analyzed by UPLC/MS/MS.

Calibration curves were spiked from 1-500 ng/mL for cortisol and from 0.05-25 ng/mL for 17a-hydroxyprogesterone (17-OHP) and androstenedione (Adione) .



#### **Chromatographic Conditions**

| Table 1—UPLC Gradient |      |            |    |
|-----------------------|------|------------|----|
| Time                  | Flow | % <b>A</b> | %B |
| 0.0                   | 0.6  | 70         | 30 |
| 1.0                   | 0.6  | 50         | 50 |
| 2.0                   | 0.6  | 45         | 55 |
| 2.5                   | 0.6  | 5.0        | 95 |
| 3.5                   | 0.6  | 5.0        | 95 |
| 3.6                   | 0.6  | 70         | 30 |
| 4.5                   | 0.6  | 70         | 30 |
|                       |      |            |    |

#### Waters Xevo<sup>®</sup> TQ-S Conditions, ESI +

| Capillary Voltage:     |  |  |  |
|------------------------|--|--|--|
| Desolvation Temp:      |  |  |  |
| Cone Gas Flow:         |  |  |  |
| Desolvation Gas Flow:  |  |  |  |
| Source Temp:           |  |  |  |
| MRM transitions, ESI+: |  |  |  |



#### Table 2

| Compound               | RT   | MRM<br>Transitions        | Cone<br>(V) | Coll.<br>Energy |
|------------------------|------|---------------------------|-------------|-----------------|
| Cortisol               | 0.72 | 337.2>121.1<br>337.2>91.0 | 42<br>42    | 22<br>52        |
| Androstenedione        | 1.50 | 287.2>97.1<br>287.2>109.0 | 58<br>58    | 20<br>26        |
| 17a-OH<br>progesterone | 1.55 | 331.2>97.1<br>331.2>295.2 | 58<br>58    | 26<br>16        |



**Figure 1.** Chromatography of a calibration standard of plasma steroids. Cortisol—20 ng/mL. Other steroids—1 ng/mL.

# THE SCIENCE OF WHAT'S POSSIBLE.



**Figure 3.** Residual phospholipid traces for samples prepared with Oasis PRIME HLB (A) and by protein precipitation (B). > 97% of phospholipids were removed compared to PPT. A 20x magnification of the SPE extract is shown at upper right

#### QUANTITATIVE RESULTS



*Figure 4.* Calibration curves for cortisol, androstenedione and 17-OHP

|                     | R <sup>2</sup> | Mean %<br>Dev. | Calibration<br>Range<br>(ng/mL) |
|---------------------|----------------|----------------|---------------------------------|
| Cortisol            | 0.996          | 8.4%           | 1.0-500                         |
| Androstenedione     | 0.989          | 8.0%           | 0.05-25.0                       |
| 17a-OH progesterone | 0.993          | 9.7%           | 0.05-25.0                       |

Table 3. Calibration summary for all analytes.

|                   | Cortisol |      |  |
|-------------------|----------|------|--|
| QC Spike<br>Conc. | Accuracy | %CV  |  |
| 3.0               | 92.3%    | 5.4% |  |
| 30                | 94.8%    | 3.0% |  |
| 300               | 94.9%    | 6.0% |  |

|                   | Androstenedione |      | 17a-OH progesterone |      |
|-------------------|-----------------|------|---------------------|------|
| QC Spike<br>Conc. | Асс             | %CV  | Асс                 | %CV  |
| 0.15              | 94.4%           | 5.7% | 93.7%               | 6.5% |
| 1.50              | 95.0%           | 3.6% | 92.3%               | 5.6% |
| 15                | 95.4%           | 5.5% | 93.7%               | 6.5% |





**Figure 2.** Recovery and matrix effects for plasma steroids (N=4).

**Table 4.** Quality Control (QC) results for plasma steroids extracted from plasma using Oasis PRIME HLB (N=6).

#### CONCLUSIONS

- Direct loading and extraction of samples without conditioning and equilibration results in a faster extract with less solvent usage
- Near complete elimination of phospholipids compared to protein precipitation
- Consistent recovery with minimal matrix effects
- Excellent analytical sensitivity, linearity, accuracy and precision for all compounds

For research use only– Not for use with diagnostic procedures

## TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS