ASMS 2000 faxBACK

Please photocopy and fax back to: Micromass' International Marketing Response Centre

USA: (978) 524 8210 EUROPE: +31 (0) 294 419052 UK / International: +44 (0) 161 282 4400

	_
Please send me	
information on the	
following products	
from Waters and	
Micromass:	/

- 2 please tick relevant box(s)
- Waters' CapLC™ System
- Waters' Alliance[®] HT System
- Waters' ZMD 2000 / 4000 \square
- Micromass' LCT™

Name: Position:

Tel:

Organization: Address:

Application:

Immediately

- Micromass' MUX-technology™
- Micromass' Quattro LC™
- Micromass' Quattro Ultima™

- Micromass' Q-Tof™ 2
- Micromass' TofSpec 2E™

- Micromass' GCT™
- Micromass' AutoSpec[™] NT
- Micromass' OpenLynx™
- Micromass' QuanLynx™
- Micromass' MetaboLynx™
- Micromass' ProteinLynx™
 - Micromass' MassPrep Station™
- Micromass' M@LDI™
- Micromass' MicrobeLynx™
 - Micromass' ProteomeWorks™ System

MICROMASS UK Limited Floats Road Wythenshawe Manchester M23 9LZ Tel: +44 (0) 161 945 4170 Fax: +44 (0) 161 998 8915 http://www.micromass.co.uk

UK Sales Desk Tel: 0161 946 0565 USA Beverly MA, Tel: 978 524-8200 Canada Pte-Claire. Tel: 514 694-1200 EUROPE Weesp, Tel: +31 (0) 294-480484 Netherlands Weesp, Tel: 0294-480484 Austria Achau, Tel: 022 3672 668 Belgium Vilvoorde, Tel: 02-2534550 France Villeurbanne, Tel: 04 72 14 89 00 Germany Eschborn, Tel: 0800-1817249 Spain Barcelona, Tel: 34 932 466 696 Switzerland Sissach, Tel: 061-9730800 NORDIC Täby, Tel: +46 (0) 8-756 71 30 Sweden Täby, Tel: 08-756 71 30



For more information please contact your local Micromass office or visit our web site: www.micromass.co.uk





UNAMBIGUOUS CHARACTERISATION OF SYNTHETIC PATHWAYS BY EXACT MASS GC-MS ANALYSIS USING A GCT MASS **SPECTROMETER**

²V. Garner Manchester M23 9YJ, UK



MALDI-MS

ICP-MS

LC-ICP-MS

IR-MS





Poster No. MMP/83/LAG/V1.00 oduct information in this document was corr of printing, however specifications are subject to alteration without notice. Information contained herein should not be construed as a mitment either by Micromass Lim or Waters Corporation. All orders are accepted subject to our current Conditions of Sale For research use only Not for use in diagnostic procedur

Microsoft Windows NT is a registered trademark of the Microsoft Corporation. No attempt is made to supersede this or any other copyrights.

© Micromass Limited June 2000 Micromass is a division of Waters Corporation ¹A. Newton^{*}, ¹M. Green, ¹P. Hancock, ¹R.H. Bateman, ¹R. Carruthers, ¹A. Gilbert,

¹Micromass UK Ltd. Floats Road, Manchester M23 9LZ, UK. ²Hall Analytical,

48th ASMS, 2000

LONG BEACH CALIFORNIA

11th - 15th June, 2000

*To whom all correspondence should be addressed:

Tel: + 44 (0) 161 945 4170 Fax: + 44 (0) 161 998 8915 e-mail: anthony.newton@micromass.co.uk



Overview

PURPOSE

- Characterisation of synthetic pathways by exact mass and elemental composition determination with EI/CI/FI.
- Differentiation of compound sources/manufacture by exact mass GC-MS profiling.

METHOD

- Orthogonal Time of Flight mass spectrometer.
- Exact mass using single point internal reference correction.

RESULTS

- Unambiguous identification of the synthetic pathway for the effect compound Speedcure™.
- Characterisation of FAME extracts from differing batches.

Introduction

- The identification of minor components in synthetic products has important legal and financial implications in many areas of commercial and forensic chemistry.
- These include the profiling of manufactured products in order to protect unique synthetic pathways and the determination of trace impurities in illicit drugs of abuse in order to correlate seizures with the site of production.
- In this paper two particular applications of a GCT in industrial and forensic analysis are used to exemplify the benefits of full spectra exact mass measurement.

Instrumentation

• The orthogonal acceleration time-of-flight mass spectrometer allows fast acquisition of full spectra, with high sensitivity and elevated resolution (7000FWHM) **Figure 2**. This elevated resolution reduces mass interference's. Furthermore the precise, and stable relationship between ion arrival time and the square root of its mass allows good measurement accuracy with only a single rnal reference mass.

Figure 1 shows a schematic of the GCT system.

- Ions produced in the grounded ion source are accelerated to 40eV and focused into a parallel ion beam. As the ions traverse the pushout region, a sudden voltage publics is applied, ejecting a portion of the beam orthogonally. A single stage reflectron reflects the ions back to a dual microchannel plate detector (MCP). Ion arrivals are recorded using a 3.6 GHz time to digital converter (TDC).
- The ion beam is sampled orthogonally at up to 30,000 times / second. Individual time of flight spectra are summed before being transferred to a host PC. The instrument is fitted with a heated GC transfer line, heated reference inlet for the introduction of volatile calibration compounds and supports direct insertion probe operation.



GCT: Doublet from a mixture of PFTBA and C₄Cl₆



Chromatographic conditions ● Analysis of Speedcure[™]



- GC-FI: HP1 30M x 0.32mm ID x 0.25µM Injection: Splitless. Oven: 50/0.8 -250@10
- Profiling of fatty acid methy esters
- Rtx-5-MS 30M x 0.25mm ID x 0.25µM Injection: Splitless. Oven: 50/0.8 -275@25

Characterisation of Synthetic Pathway

Speedcure™ is one of a range of polymer curing agents. Chemically it is ethyl 4-N,Ndimethylaminobenzoate. Its industrial synthesis via a reductive formylation route is subject to patent protection. An alternative common route for the synthesis of this compound is direct alkylation

The reaction pathway for the reductive formulation route is shown in Figure 3.

Unique co-products, formed in minute but sufficient yield, allow the actual pathway to be identified



Speedcure[™] was analysed by exact mass GC-MS in EI+, CI+ and FI (Field ionisation) modes. Characteristic reaction intermediates were then identified.

Figure 4 shows the EI + GC-MS TIC produced.

Figures 5 and 6 show the EI+ spectra for the components at RT 9.15 and 11.9 mins act mass measurement, empirical formula and proposed structure

Table 1 shows the statistical results from 5 repeat exact mass measurements for the two components highlighted

Figures 7 and 8 show the complimentary exact mass molecular ion information produced by exact mass FI and Methane CI+



Discussion

Two characteristic trace components were identified and their elemental composition determined by exact mass measurement. Complimentary exact mass CI and FI data confirm the molecular ion assignment.

The two components identified at m/z 177.0794 and 207.0901 are unique to the reductive formylation route described in Figure 3, unambiguously identifying the synthetic route as reductive formylation.

The pathway for the direct alkylation route is shown in Figure 9 for comparison. No unique by-products associated with this reaction were dete In the field of forensic analysis of drugs of abuse materials are often seized as tablets

Exact mass GC MS profiling

which incorporate residues of processing reagents such as magnesium stearate. It is routine procedure to isolate the free fatty acids from a small portion of a tablet and convert to the corresponding methyl esters. As the stearate is obtained from natural sources in batch processes, the MS profile of this FAME extract can be used to relate seizures of tablets to site(s) of producti









RMS error ppm= 5.4 ppm



RMS error ppm = 4.0 ppm

EI+GC-MS Exact mass measurement results from five repeat injections of Speedcure™



Table 1

Two batches of Magnesium Stearate (X and Y) from differing sources were analysed using the GCT. The sample was first treated with 1M HCl. partitioned into hexane. and after drying esterified using BF3 in the presence of Methanol.





Results

Figures 10a and b show a comparison of the TIC's from both batches. Figure 11 shows the library search and elemental composition data for one of the

tty acid methy esters identified.

Figure 12 shows the the library search and elemental composition data for one of the es characteristic of batch X).

Figure 13 Shows accurate mass chromatograms (10mda window) of a characteristic nined for the dioates in X and Y.

Figure 14 shows the nominal mass (0.5da window) chromatograms of a characteristic ion determined for the dioates in X and Y.











Discussion As expected many of the fatty acid methy esters found were common to both batch ${\bf X}$ and Y. However only sample X was found to contain a series of dioates. Library search data was compared with the exact mass measurements obtained to confirm the assignment. The ion at m/z 124.0888 corresponding to the fragment $C_8H_{12}O$,

was found to be unique to these dioates. Exact mass chromatogram, 10mda window (Figure 13) offer a very high degree of selectivity compared to a nominal mass chromatogram, 500mda window (Figure 14). In this example the selectivity is equivalent to a magnetic sector mass spectrometer, in SIR mode, at 12,400 resolution 10% valley definition. The combination of the selectivity and specificity of exact mass and full library search using the GCT provides a clear and unambiguous differentiation of these two

batches of sample



• Exact mass measurement using GCT-oaTof MS allows unambiguous and cost effective characterisation of synthetic pathways

MMP/83/LAG/V1.00