## GENERAL INFORMATION FOR ELUTION SOLVENT, SORBENT AMOUNT AND SAMPLE SIZE

## SORBENT AMOUNT SELECTION FOR OASIS® HLB VS. C18 SILICA

C <sub>18</sub> -SILICA-BASED SORBENTS AMOUNT PER CARTRIDGE OR WELL	OASIS® HLB SORBENT AMOUNT PER CARTRIDGE OR WELL	
50 mg	10 mg	
100 mg	30 mg	
200 mg	60 mg	
500 mg	1 <i>5</i> 0 mg to 200 mg	
1000 mg	200 mg to 500 mg	

The table gives the amount of Oasis® sorbent you need based on the amount of silica-based sorbent currently being used.

## SELECTING ELUTION SOLVENTS BY RELATIVE STRENGTH

SOLVENT	SOLVENT TYPE	RELATIVE ELUTION STRENGTH*	COMMENTS
Methanol	proton donor	1.0	low toxicity
Acetonitrile	dipole-dipole	3.1	low viscosity
Tetrahydrofuran	dipole-dipole	3.7	check for peroxides before concentrating
Acetone	dipole-dipole	8.8	least toxic high UV-cut-off
Ethyl Acetate	dipole-dipole	high	least solubility in water
Methylene Chloride	dipole-dipole	high	very low solubility in water

\* High Purity Solvent Guide, Burdick & Jackson Laboratories, Inc. Solvent Properties of Common Liquids, L.R. Snyder, J. Chromatgr., 92, 223 (1974); J. Chromatogr. Sci., 16, 223 (1978)

## CHOICE OF OASIS® SORBENT WEIGHT BASED ON SAMPLE SIZE

SORBENT PER WELL	MAXIMUM MASS CAPACITY	TYPICAL SAMPLE VOLUME	TYPICAL ELUTION VOLUME
2mg	0.030 to 0.050 mg	5 to 100 µL	≤ 50 µL
5 mg	0.15 to 1 mg	10 to 100 µL	≤ 150 µL
10 mg	0.35 to 2 mg	50 to 400 µL	≤ 250 µL
30 mg	1 to 5 mg	100 to 1 mL	> 400 µL
60 mg	2 to 10 mg	200 µL to 2 mL	> 800 µL

Elution volume and capacity depend on the polarity of the analyte. Stronger solvents can reduce the elution volume. Low elution volume allows the evaporation step to be eliminated.