NOVEL POROUS RESINS WITH DUAL ION-EXCHANGE/REVERSED-PHASE RETENTION BEHAVIOR

<u>Jeng-Jong Lee</u>, J.E. O'Gara, D.P. Walsh, A.R. Pelissey, E.S.P.Bouvier, P.C. Iraneta, Henry Cheng, D.J. Phillips Waters Corporation, 34 Maple Street, Milford, MA 01757, USA

Introduction

In 1996, Waters designed and commercialized the Oasis® HLB sorbent, a water-wettable copolymer with hydrophilic-lipophilic balance (HLB) for SPE. Unique in its ability to retain both polar and non-polar compounds and their metabolites, with high reproducibility and recovery, Oasis® HLB sorbent has set the new standard for SPE. In order to meet the demand for extraction of basic or acidic drugs, Oasis® MCX, a new member of the Oasis® family has been developed. It has exhibited high selectivity and sensitivity for extracting basic drugs. Furthermore, Oasis® MAX for acidic drugs is also being developed.

Oasis[®] MCX and MAX are brand new porous resins having a hydrophobic fragment, a hydrophilic polar moiety, and an ion-exchange functional group. In this paper, we would like to discuss their synthesis, selectivity, reproducibility and recovery. In addition, an example of using Oasis[®] MCX as a universal sorbent for isolating acidic, basic, and neutral drugs will be reported.

Synthesis Strategy

Selecting a suitable matrix for surface derivation

Disadvantage of silica matrix:

- Undesirable silanol activity
- pH limitations
- Low capacity
- Hydrophobic collapse

Advantage of Oasis® HLB sorbent:

• Porous poly(divinylbenzene-co-N-vinylpyrrolidone) resin with excellent chemical stability.

• Rigid porous resin with high surface area which provides better mechanical stability and binding capacity than that of conventional poly(styrene-co-divinylbenzene) resins.

• One of the cleanest sorbents available in the market.

• Hydrophilic-lipophilic balanced sorbent, surface modification by sulfonation, chloromethylation, and amination can be conducted without using halogenated or organic solvents.

Stability of Oasis® HLB Sorbent



Oasis® HLB

zwitterion

No indication of hydrolysis of Oasis® HLB

Oasis[®] HLB has a five member ring amide moiety. Surprisingly, no zwitterion was detected by FTIR after Oasis[®] HLB was treated under the following severe conditions:

- 140 °C in 20% KOH solution overnight
- 95 °C in 10% HCl solution overnight
- 120 °C in conc. Sulfuric acid for 23 hours

Characteristics of Oasis® HLB Sorbent

$\langle \rangle$	Specific Surface Area:	810 m ² /g
n n	Average Pore Diameter:	80 Å
	Total Pore Volume:	1.3 cm ³ /g
	Mean Particle Diameter:	30 µm
	Percent Fines < 10 µm:	0.1%
\downarrow	Pyrrolidone:	20%
Oasis® HLB	Cross-linked:	> 30%

- Waters Oasis[®] HLB sorbent is patented (US 5,882,521).
- Having greater than 30% cross-linking, which is about 2 to 10 times higher than that of conventional poly(styrene-co-divinylbenzene) resins.
- A highly clean sorbent purified by Waters proprietary processes after suspension polymerization.

Synthesis of Oasis® MCX Sorbent



Oasis[®] MCX is prepared by treating Oasis[®] HLB in conc. H₂SO₄ at 25°C for 1 hour, which is a very environmentally friendly process. Whereas, conventional poly(styrene-co-divinylbenzene) ion-exchange resins are prepared at high temperature with halogenated or organic solvents.

• Since Oasis[®] MCX is prepared at low temperature without using solvents, no contaminants are generated during the sulfonation.

 Cleanliness test of Oasis[®] MCX shows the extractable components less than 6 μg/g, which is as clean as Oasis[®] HLB.

• Oasis[®] MCX is a unique porous resin with a hydrophobic fragment, a hydrophilic polar moiety, and an ion-exchange functional group.¹

Development of Synthesis Process Ruggedness and robustness

Factorial experiment design was used to study the following synthesis variables:

- Components of sulfonation reactants
- Concentration of reagents
- Temperature of reaction (-30 °C to 120 °C)
- Reaction time (1 minute to 24 hour)

Summarizing from the above studies, a synthesis process with only two variables is developed and their relations are as follows:

Degree of sulfonation (ion-exchange capacity, meq/g of $-SO_3H$)

= 0.53 meq/g + 0.018 meq/g °C x (temp. °C) + 0.017 meq/g h x (time, h); R² =0.996

 It will take more than 5 °C or 5 hours variation to have about 10% change in ion-exchange capacity (1.0 meq/g of -SO₃H) of Oasis[®] MCX sorbent.

• During the product development, twelve different batches were prepared with reactors ranging from 0.5 L to 80 L by several people in different locations. The $-SO_3H$ loading was 0.98 meq/g, $\sigma = 0.05$ meq/g.

Why Oasis[®] MCX Sorbent has 1.0 meq/g of Cation-Exchange Capacity?

Effect of ion-exchange capacity on retention behavior



 Prototypes ranging from 0 to 3.2 meq/g of -SO₃H were packed in HPLC columns to study the retention of neutral and basic analytes. Neutral compounds' retention decreases greatly when sulfonation is above 1 meq/g. Basic compounds' retention increases when sulfonation is increased to about 1 meq/g. However, retention decreases at higher levels of sulfonation.²

• At 1.0 meq/g sulfonation, Oasis[®] MCX has both maximum ionexchange and hydrophobic retention capabilities enabling the retention of acidic, neutral and basic drugs.

Drug- Sorbent Interactions



• The strong acid sulfonate ($-SO_3H$) groups on the polymer interact with basic and cationic solutes.

• The amide groups enable the resin to have polar interactions and hydrogen bonding capabilities with polar analytes.

• The hydrophobic fragments of divinylbenzene and vinylpyrrolidone enable the affinity towards nonpolar substances through hydrophobic interaction.

Extraction of Basic Drugs on Oasis® MCX Cartridges

Basic drugs ⁴	Conc. (µg/mL)	% Recovery	% RSD
Methadone	0.01	100	2.7
Methadone metabo	lite 0.04	98	1.0
Estazolam	0.2	97	1.5
Propranolol	0.08	100	2.7
Oxprenolol	0.08	92	3.2
Metoprolol	0.08	100	3.4
Verapamil	4.5	100	0.2
Nor-Verapamil	3.4	98	0.5
Verapamil-Methoxy	6.7	100	3.3
Codeine	0.12	99	1.2
Codeine (metabolite	e) 0.12	100	1.5
Ranitidine	0.24	93	2.6

Generic Oasis[®] MCX method for extraction of basic drugs from human urine: ³

- Load 3 ml spiked and acidified urine into a cartridge
- Wash with 2 mL of 0.1N HCl
- Wash 3 or 4 times with 2 mL of methanol
- Elute with 2 mL of 5% ammonium hydroxide in methanol
- Evaporate and reconstitute in 300 µL of 20% methanol

Without preconditioning Oasis[®] MCX cartridges, routine recoveries greater than 85% with RSDs less than 5% can be achieved across a broad range of basic drugs.³

Fractionation with Oasis® MCX Cartridges

Drugs ⁴	Conc. (ւ <mark>g/mL)</mark>	Fraction	% Recovery	% RSD
Acetaminop (neutral drug	hen 3)	2	Elute /	A 100	0.6
Barbital (weak acid o	drug)	2	Elute A	A 98	1
Amphetamir (basic drug)	ne	4	Elute E	3 96	1
Methamphe (basic drug)	tamine	4	Elute	B 90	4

Generic Oasis[®] MCX method for extractions of neutral, acidic, and basic drugs from plasma: ³

- Equilibrate cartridges with 1 mL methanol/1 mL water
- Load 1 ml spiked and acidified plasma into a cartridge
- Wash with 1 mL of 0.1N HCl
- Elute A: with 1 mL methanol
- Elute B: with 1 mL of 5% ammonium hydroxide in methanol
- Evaporate and reconstitute

Oasis[®] MCX fractionated samples into the combined acidic and neutral drug fraction (Elute A) and a basic drug fraction (Elute B) with high recovery.

Synthesis of Oasis® MAX Sorbent



Oasis® MAX (Patent pending)

 Designed to be a mixed-mode anion exchange/reversed-phase sorbent capable of retention of acidic, neutral and basic compounds and with high selectivity to cleanup acidic drugs.

- Prepared through an environmentally friendly process without halogenated or organic solvents.
- Currently in the ruggedness and robustness testing stage at Waters' ISO 9002-certified plant under cGMP standard and available soon.

References

 Jeng-Jong Lee, John O'Gara, "Novel Sulfonated Porous Resins for Solid Phase Extraction and Chromatography" U. S. Patent Application filed on June 12, 1998.

 James S. Fritz reported similar effect on sulfonated polystyrene-divinylbenzene resins; Chambers, T. K.; Fritz, J. S.; *Journal of Chromatography A*, 797 (**1998**) 139-147.

3. Oasis® Sample Extraction Products Applications Notebook, ©1999 Waters Corporation, Literature WN047.

4. Structure of analytes:

