

Michelle Wood<sup>1</sup>, Marleen Laloup<sup>2</sup>, Gert De Boeck<sup>2</sup>, Nele Samyn<sup>2</sup>, Vivienne Maes<sup>3</sup> and Maria del Mar Ramirez Fernandez<sup>2</sup>.

<sup>1</sup>Waters Corp., Micromass Technologies Centre, Manchester, UK; <sup>2</sup>National Institute of Criminalistics and Criminology (NICC), Brussels, Belgium; <sup>3</sup>Free University of Brussels, Belgium

### INTRODUCTION

Drug-facilitated crime (DFC) is a relatively new term for an old practice. According to legend 'Slipping a Mickey' relates to a practice made infamous by Chicago saloon owner Mickey Finn in the late 1890's. Finn secretly laced the drinks of his patrons, with drugs, in order to knock them unconscious. After which, he and his wife would strip them of their valuables. Victims would awaken later - remembering nothing!



Figure 1: The age-old practice of 'Slipping a Mickey'; the surreptitious administration of a drug for criminal gain.

Over the last few years DFC has been increasing. DFC includes robbery, and assault including sexual assault; drug-facilitated sexual assault (DFSA). A few years ago the Society of Forensic Toxicologists (SFT) formed a committee to address various aspects of DFSA and drew up a list of compounds which have/could be implicated in DFSA in the United States<sup>1</sup>. They include illegal drugs, prescribed medications and 'over-the-counter' preparations.

Owing to the diversity of the analytes involved, a variety of analytical techniques are usually required, including immunoassay, GC-MS, GC-FID and LC-UV. Our aim was to develop a simple, generic method to screen for these analytes using a single analytical technique based on LC-*oa*-ToF.

### METHODS

#### Reference library

A reference library was created for the compounds of interest which included those identified by the SFT DFSA committee, in addition to other analytes which are relevant in Europe. Where reference material was available (~70 compounds), standard exact mass spectra were acquired using a low aperture voltage to generate exact mass spectra and retention time (RT) information. Spectra were also acquired at a higher aperture voltage to generate fragments by collision-induced dissociation (CID) within the source. All spectra were quality checked using I-FIT software, which compares measured and theoretical isotope ratios<sup>2</sup>.

Where reference material was unavailable (~30 compounds) theoretical spectra were added to the library using the elemental formula.

#### Chromatography

ACQUITY UPLC™

Column: ACQUITY™ HSS C<sub>18</sub> (2.1 x 100mm, 1.8 μm)  
Column temp: 30 °C  
Injection vol: 10 μL  
Solvent A: 0.05 % formic acid  
Solvent B: Methanol

#### Mass Spectrometry

Mass spectrometer: Waters LCT Premier™  
Ionisation mode: Electrospray +ve  
Capillary voltage: 3kV  
Aperture 1: 10 and 50V  
Mass range: 50–600Da  
Resolution: 10,000 (W-mode)  
LockSpray™ reference: Leucine enkephalin  
[M+H]<sup>+</sup> = 556.2771

#### Sample Preparation

Urine samples were subjected to liquid:liquid extraction (LLE) under acidic and basic conditions as described below.

**Acidic:** 125 μL acetate buffer (pH 3.5) was added to 250 μL urine. After mixing, samples were extracted using 750 μL solvent mix (DCM:ether:hexane containing 0.5% isoamyl alcohol). Samples were again mixed for 2 min, followed by centrifugation at 3000 rpm for 5 min.

**Basic:** As above but using borate buffer (pH 9.5).

Supernatants from the acidic and basic extractions were pooled before drying, followed by reconstitution in mobile phase prior to LC-MS analysis.

#### Data processing

Data was processed automatically using Chromalynx™ software (Waters).

### RESULTS & DISCUSSION

Where reference material was available, pure standards were injected and acquired, in full scan mode, at low and high aperture voltages to create a library which included exact mass, fragmentation and RT. Figure 2 shows the library data for MDMA which elutes at 2.7 min.

The chromatogram obtained on injection of a urine sample, spiked with morphine, hydromorphone and 7-aminoclonazepam (7-AC) is given in Figure 3. This data clearly demonstrates the benefits of exact mass and the ability to differentiate from other compounds of similar mass. Figure 3A shows the extracted TIC for the nominal mass i.e.  $m/z$  286. All 3 compounds are clearly visible. Figure 3B shows exactly the same data but with the exact mass for 7-AC ( $C_{17}H_{12}N_2OCl$ ) i.e.  $m/z$  286.0747  $\pm$  5 mDa extracted. The exact mass for both morphine and hydromorphone is  $m/z$  286.1443 [ $C_{17}H_{19}NO_3$ ].




Figure 2: The library entries for MDMA (Ecstasy) at aperture settings of 10 (A) and 50V (B) respectively.

Data processing was achieved using Chromalynx™ which automatically deconvolutes the data. The most intense ions (up to a max. of 8) are extracted using a defined mass window (either in mDa or ppm). The resultant data was then searched against the newly created DFC library.

Figure 4 shows a typical result for an authentic urine sample. Dihydrocodeine (DHC), acetaminophen and EDDP were identified by the new LC-*oa*-ToF screening method. These findings were consistent with the data obtained using the conventional screening methods i.e. HPLC-UV and immunoassay.

Library hits were also scored according to the agreement between measured masses in the acquired spectrum and those which are theoretically possible based on the elemental composition of the compounds in the library. Additional confidence of identified analytes was achieved using I-FIT (Figure 4B and C).




Figure 3: Extracted mass chromatograms for the nominal mass (Figure 3A) and exact mass of 7-AC (Figure 3B).




Figure 4A: Results for an authentic urine sample following Chromalynx™ processing. Figure 4A shows the list of proposed candidates and the average match to the library. The spectrum view allows direct comparison of acquired spectra with library data.




Figure 4B shows the spectral data for acetaminophen which was also identified in this sample. The results for the proposed candidate 'acetaminophen' was submitted for elemental composition. Figure 4C shows the proposed elemental composition for a  $m/z$  152.0711, the result was  $C_8H_9NO_2$  which corresponds to that of acetaminophen.

### CONCLUSION

- A spectral library has been created for over 100 analytes which have been/may be implicated in DFC.
- Where reference material was available (~70 analytes) spectra were collected at low and high voltages to generate fragments.
- Where reference material was available identification was based on a combination of exact mass, fragment ions and RT.
- Where reference material was unavailable (~30 analytes) identification was based on exact mass.
- Data was automatically processed using Chromalynx™ software, which effectively detects chromatographic peaks, even in complex matrices.
- The exact masses of proposed candidates were compared and scored against theoretical mass.
- Additional confidence of identification for any proposed candidates was achieved using I-FIT software.

#### References

1. M. LaBrec, W. Andols, W. Hearn, R. Basak, E. Cone, B. Finkle, D. Fraser, A. Jenkins, J. Mayer, A. Nagraz, A. Patis, H.C. Wells, L. Raymon, M. Robertson and J. Steady. J. Forensic Sci. **44**, 227-230 (1999).
2. K. Hobbay and B. Bennett. Poster Presentation, WSC, Prague 2006.