



Experiments in Liquid Chromatography — Separations of FD&C Dyes

452020

This experiment is one in a series prepared by Waters Associates, Inc. to introduce liquid chromatography to undergraduates.

You are encouraged to copy the inside pages for your students use.

Identification of FD&C Dyes

- Isolation by a SEP-PAK C₁₈ Cartridge
- Separation by Thin Layer Chromatography (TLC)

Objectives

To provide an introduction to both "normal phase" and "reverse phase" chromatography and to demonstrate applications in the isolation and identification of Food, Drug, and Cosmetic (FD&C) dyes.

Introduction

Within the past few years, much interest has been generated over the safety of FD&C dyes. At present, seven of these dyes are allowed in foods: Reds 3 and 40, Yellow 5 and 6, Blues 1 and 2, and Green 3. (Two additional dyes, Reds 2 and 4, have been officially "delisted" although foods packaged before the latter part of 1976 may contain them). Information on these dyes may be found in references 1-7. Because of continuing concern over dye safety and the fact that present law does not require listing of specific dyes on labels, there is genuine interest in dye separation and identification. This experiment involves the isolation of these dyes by SEP-PAK C₁₈ Cartridge and subsequent analysis by TLC.

Choosing a Food

For ease of sample preparation, the following "foods" are ideal substrates for this experiment: a certified food color, soft drinks,

pickle juice, juice from maraschino cherries, colored cough medicines, mouthwashes, or gelatin type desserts. Solid foods require dissolution and filtering before analysis. The food should indicate "artificial color" on the label.

Preparing Standard Dye Solutions

Stock solutions of the dyes are prepared by dissolving trace amounts of the dye in distilled water.

Equipment and Supplies

- Standard dyes (Note 1)
- Colored foods
- SEP-PAK C₁₈ Cartridges (one/student)
- TLC developing solution: n-butanol/ethanol/water/concentrated aqueous ammonia (50:50:25:10)
- Silica TLC plates
- 2-propanol
- 1% acetic acid

Procedure: Dye Extraction

For reproducible results, it is necessary to separate the dye from interfering substances such as water, alcohol, flavorings, etc. before TLC analysis. Standard methods for this separation have traditionally involved either tedious ion exchange chromatography or time-consuming white wool extractions (8).

In this experiment, a SEP-PAK C₁₈ Cartridge is used to isolate and "trace enrich" the dye in seconds. The dyes are adsorbed onto the C₁₈ packing, allowing their concentration (trace enrichment) on the cartridge head while numerous interfering compounds (sugars, flavorings, preservatives, etc.) are eluted. The specific procedure is as follows:

Refer to Water Publication "SEP-PAK Cartridges for

Information for Instructor continued on page 4



Lab Work Sheet

Identification of FD&C Dyes

- Isolation by a SEP-PAK C₁₈ Cartridge
- Separation by Thin Layer Chromatography (TLC)

Objectives

To introduce the theory of liquid chromatography and to demonstrate applications in the isolation and identification of FD&C dyes.

Introduction

Within the past few years much interest has been generated over the safety of Food, Drug and Cosmetic (FD&C) dyes.

Presently, seven dyes are allowed in foods: Reds 3 and 40, Yellow 5 and 6, Blues 1 and 2, and Green 3. (Two other dyes, Reds 2 and 4, have been officially "delisted" although foods packaged before the latter part of 1976 may contain them.) In this experiment, you will isolate dyes from a food with a SEP-PAK C₁₈ Cartridge and then separate and identify them by TLC.

SEP-PAK C₁₈ Cartridges are small chromatographic columns for rapid sample preparation. Their packing material is made by chemically bonding a C₁₈ group to silica.

The procedure described here can be modified to analyze other food additives such as preservatives, artificial sweeteners, etc.

Preparation of Food Sample

For ease of sample preparation, liquids or easily liquified foods are recommended for analysis. Acidic foods, i.e. soft drinks, pickle brine, and juices may be used without special preparation. Gelatin products should be dissolved in water (0.1 g in 3 ml), acidified (1 drop of 1% acetic acid), and filtered before analysis. Alcohol containing materials (mouthwashes, cough medicines, etc.) should

be diluted 1:1 with water and acidified (1 drop of 1% acetic acid) before analysis to enhance dye retention on the SEP-PAK C₁₈ Cartridge.

Procedure: Dye Extraction

For reproducible results, the dye must be separated from interfering substances such as water, alcohol, flavorings, etc., before chromatographic analysis. A SEP-PAK C₁₈ Cartridge isolates and concentrates ("trace enriches") the dye in seconds. The dyes are adsorbed onto the C₁₈ packing while numerous, more polar compounds (sugars, flavorings, etc.) are eluted from this reverse phase cartridge. The specific procedure is as follows:

- 1) Use a syringe to prewet a SEP-PAK C₁₈ Cartridge with 3 ml of 2-propanol.
- 2) Flush 5 ml of 1% acetic acid through the cartridge. This improves dye adsorption by insuring that the sulfonic acid dyes remain in their free acid form.
- 3) Flush 1 ml of sample through the cartridge, discard the eluate.
- 4) Flush with 2-3 ml of 13% 2-propanol, collecting only the colored portion (which may contain Blue 2, Yellow 5 and 6,

and Red 40).

- 5) Flush with 2 ml of 50% 2-propanol, collecting only the colored portion (which may contain Red 3, Blue 1 and Green 3).

- 6) The eluates from Steps 4 and 5 will be used for chromatographic analysis.

Procedure: TLC Analysis

Spot a TLC silica plate with some of the dye extracted according to the directions of your instructor. Develop the plate with a solution of n-butanol/ethanol/water/concentrated aqueous ammonia (50:50:25:10). The R_f values of the seven standard dyes are found in Table I. Table II lists dyes found in common "foods."

Food for Thought

- Calculate R_f values for all colored components.
- Which dyes were in your food sample?
- Would a change in composition of the developing solution change the R_f values? Why or why not?
- Discuss how you might analyze a food for an artificial sweetener using a modification of this procedure.

Table 1. R_f Values for FD&C Dyes

Dye	R _f Value *	Comment
Red 3	0.77	Not used extensively
Red 40	0.65	Most common red dye
Yellow 5	0.47	Bright yellow
Yellow 6	0.57	Bright orange
Blue 1	0.57	Most common blue dye
Blue 2	0.61	Not commonly used — fades rapidly in light
Green 3	0.61	Not used extensively

*The R_f values may vary depending upon the activity of the silica, but their relative values should remain constant.


Lab Work Sheet

Table II. FD&C Dyes Present in Food Products


Food	FD&C Dye
Red cough medicine	Red 40, Yellow 6
Pickle brine	Yellow 5
Orange soda	Yellow 5 and 6
Grape soda	Red 40, Blue 1
Red mouthwash	Red 40
Green mouthwash	Yellow 5, Blue 2
Yellow food color	Yellow 5
Red maraschino cherry juice	Red 40

Structures of FD&C Dyes


BLUE 1




RED 40




BLUE 2




YELLOW 5




GREEN 3



YELLOW 6



RED 3



References

- (1) "Separation and Identification of Some FD&C Dyes by TLC", McKone, H.T., Nelson, G. J., *J. Chem. Educ.*, **53**, 722, (1976).
- (2) "Identification of FD&C Dyes by Visible Spectroscopy", McKone, H.T., *J. Chem. Educ.*, **54**, 376, (1977).
- (3) "Red Food Coloring — How Safe Is It?", *Consumer Reports*, **38**, 131 (Feb. 1973).
- (4) "Countdown on Color Additives", Hopkins, H., *FDA Consumer*, **10**, 5, (1976).
- (5) "Behavioral Disturbances Learning Disabilities and Food Additives", Feingold, B., *Chemtech.*, **5**, 264, (1975).
- (6) "Handbook of Food Additives", 2nd ed, Furia, T.E. (editor), The Chemical Rubber Co., Cleveland, OH, (1972).
- (7) "Additives for Eye Appeal", Damon, G. E. and Janssen, W.F., *FDA Consumer*, **7**, 15, (1973).
- (8) "High Pressure Liquid and Thin Layer Chromatography of Synthetic Acid Fast Dyes in Alcoholic Products", Martin G.E., et al., *J. Assoc. Off. Anal. Chem.*, **61**, 908, (1978).

SP 79020 LIB #: 852020
Waters, Division of • Waters Experiment, Waters Liquid Chrom
S85 (April) 1979

Scale SA FA
MktGrp MktSubGrp
Compd Grp ComSbGrp
Confidence: DM
ED :PM
PG
DY

Rapid Sample Preparation", pages 4 and 5 before using a cartridge.

- 1) Use a syringe to prewet a SEP-PAK C₁₈ Cartridge with 3 ml of 2-propanol.
- 2) Flush 5 ml of 1% acetic acid through the cartridge to improve adsorption of the dye on the C₁₈ packing.
- 3) Flush 1 ml of acidified food sample through the cartridge — discard the eluate (see Note 2).
- 4) Flush with 2-3 ml of 13% 2-propanol, collecting only the colored portion (which may contain Blue 2, Yellow 5 and 6, and Red 40).
- 5) Flush with 2 ml of 50% 2-propanol, collecting only the colored portion (which may contain Red 3, Blue 1, and Green 3).
- 6) The eluates from Steps 4 and 5 will be used for chromatographic analysis.

TLC Analysis

Spot a TLC silica plate with some of the extracted dye and develop with a solution of n-butanol/ethanol/water/concentrated aqueous ammonia (50:50:25:10).

Results

Table I lists the R_f values for standard dyes.

Summary

The identification of FD&C dyes by SEP-PAK Cartridge isolation and TLC analysis is an instructive means of introducing the theory of chromatography to students. The SEP-PAK C₁₈ Cartridge isolation shortens the preparative procedure considerably, uses fewer chemicals, and is less hazardous than standard methods. This experiment can be extended to include High Performance Liquid Chromatographic (LC) analysis. Information on the LC analysis of FD&C dyes may be obtained from Waters Associates, Inc., University Marketing Division, Maple Street, Milford, MA 01757.

Note 1: Dyes may be purchased from chemical supply houses or easily extracted by a SEP-PAK C₁₈ Cartridge from "foods" (particularly mouthwashes) that list specific colors on their labels.


Note 2: To insure retention of the dye on the column, the food sample should be acidified (one drop of 1% acetic acid) before analysis. In this way, the sodium sulfonate dyes are converted to the free acid forms which are adsorbed more strongly on the "reverse phase" C₁₈ packing.


- For certified food colors, dilute 3 drops of dye with 1 ml of water; for solids, such as gelatin products, dissolve approximately 0.1 g of sample in 3 ml warm water and filter if necessary.

- The dyes in alcohol-based matrices such as mouthwashes are retained better on the SEP-PAK Cartridge if the liquid is acidified and diluted 1:1 with water before isolation.

Note 3: The R_f values may vary depending upon the activity of the silica; however, the relative values should remain constant. Other developing solutions—iso-butyl alcohol/methanol/water/concentrated aqueous ammonia (75:25:25:10) result in a separation similar to that illustrated in Table 1.

Note 4: If you have a valuable teaching experiment that could incorporate liquid chromatography, and would like to share it with others, contact the University Marketing Div. at Waters Associates, Inc.

Waters Associates, Waters, , and SEP-PAK are trademarks of Waters Associates, Inc., Milford, MA 01757.

Waters Associates, Waters, , and THE LIQUID CHROMATOGRAPHY PEOPLE are servicemarks of Waters Associates, Inc., Milford, MA 01757.



Waters Associates, Inc. / Milford, Massachusetts 01757 / (617) 478-2000

AUSTRALIA Waters Associates Pty., Ltd. / 4th Floor, 82-96 Myrtle Street / Chippendale, NSW 2008 / Tel.: 02-699-1377

AUSTRIA Waters Ges.m.b.H. / Schönbachstrasse 13 / A-1130 Wien / Tel.: 0222-820176

BELGIUM Waters Associates NV / Pontbeekstraat 33 / 1720 Groot-Bijgaarden / Tel.: 02-4655088

CANADA Waters Scientific Ltd. / 6480 Viscount Road, Unit No. 4 / Mississauga, Ontario L4V 1H3 / Tel.: 416/677-0580

DENMARK Waters Associates AS / Byageren 11 / DK-4000 Roskilde / Tel.: 03361080

FEDERAL REPUBLIC OF GERMANY Waters GmbH / 6240 Königstein / Taunus / Herzog-Adolf-Strasse 4 / Tel.: 06174/4051

FRANCE Waters Associates SA / 18 Rue Goubet / 75940 Paris CEDEX 19 / Tel.: 200-67-01

JAPAN Nihon Waters Ltd. / Shuwa Kioicho Park Building / 3, Kioicho, Chiyoda-ku, Tokyo / Tel.: 264-8005

NEW ZEALAND Waters Associates (Pty.) Ltd. / P.O. Box 5565 / 83 Wakefield St. / Auckland 1 / Tel.: 770392

THE NETHERLANDS Waters Associates BV / Postbus 166 / Penningweg 33 / 4870 Ad Etten-Leur / Tel.: 01608-16350

SPAIN Waters Española SA / Calle de la Coruña 29 / Barcelona 26 / Tel.: 225-14-68

SWEDEN Waters Associates AB / Riddaregatan 15 / S-40261 Göteborg 17 / Tel.: 031-80-21-12

UNITED KINGDOM Waters Associates (Inst.), Ltd. / 324 Chester Road / Hartford, Northwich / Cheshire CW8 2AH / Tel.: 0606-883846