
**Same Chemistry – Lower Detection Limits**

Waters XPOSure Aldehyde Sampling Cartridges [Box of 20, P/N WAT047205] are convenient, reproducible sampling devices for quantifying aldehydes down to the ppb range in workplace and indoor air. XPOSure Cartridges contain DNPH-coated silica. When carbonyl-containing compounds are adsorbed from air samples onto this matrix, they are derivatized according to the following scheme:



This is the same derivatization reaction used in NIOSH Method 2532, EPA procedures TO-11 and IP-6A, and ASTM D5197.

**Low Background for High Sensitivity**

XPOSure Cartridges have the lowest guaranteed background levels for reliable quantitation.

*Background per Cartridge.*

- ☐ < 0.06 µg formaldehyde
- ☐ < 0.15 µg acetaldehyde
- ☐ < 0.38 µg acetone per cartridge.

*Quantitation Limits.*

- ☐ 15 minute Short Term Exposure Limit [STEL] measurements: 20 ppb in 22.5-liter air sample.
- ☐ 8 hour Permissible Exposure Limit [PEL] measurements: 10 ppb in 48-liter air sample.

*Recommended Maximum Capacity per Cartridge.*

- ☐ 70 µg formaldehyde

**Step 1. Collect Air Sample**

Connect the female tip of a cartridge to the intake port of a suitable air pump, and draw an air sample through the cartridge.

- ☐ Suggested sampling rate: 0.1 - 1.5 L/min
- ☐ Suggested air collection volume: 15 - 60 liters
- ☐ Suggested sampling time: 15 min for STEL, 8 hrs for PEL measurements

*Note:* The pressure drop across the XPOSure Sampling Cartridge permits flow rates of up to 1.5 L/min with personal sampling pumps.

**Step 2. Elute Sample from Cartridge**

- ☐ To prepare a sample for analysis, fill a 10-cc, Luer-tipped syringe with acetonitrile [CH<sub>3</sub>CN], and connect the female end of the cartridge to the syringe tip.
- ☐ Using the syringe plunger, push 10 mLs of CH<sub>3</sub>CN through the cartridge at a rate of ≤ 3 mLs/min. Collect eluate in a 10-mL volumetric flask.
- ☐ Fill the volumetric flask to the mark with CH<sub>3</sub>CN and mix.
- ☐ Analyze a portion of this solution by HPLC.

**Step 3. Analyze by HPLC**

Column: Symmetry® C<sub>18</sub> Column,  
3.9x150 mm [ P/N WAT046980]


Mobile Phase: 45/55 CH<sub>3</sub>CN/H<sub>2</sub>O (v/v)

Flow Rate: 1.3 mL/min (isocratic)

Detection: UV at 360 nm

Sample Injection Volume: 20µL

Lab Air Containing 20 ppb Formaldehyde – 22.5-L Air Sample



*Pamela Iraneta*