# MRM QUANTIFICATION OF HOST CELL PROTEINS IN PROTEIN BIOPHARMACEUTICALS

Catalin E. Doneanu<sup>1</sup>, <u>Weibin Chen<sup>1</sup></u>, Alex Xenopoulos<sup>2</sup>, Jim Murphy<sup>1</sup>, Keith Fadgen<sup>1</sup>, StJohn Skilton<sup>1</sup> and Jeff Mazzeo<sup>1</sup> <sup>1</sup> Waters Corporation, Milford, MA, USA; <sup>2</sup> Millipore Corporation, Bedford, MA, USA

#### **OVERVIEW**

A high-throughput LC-MRM assay for quantification of low-abundance HCPs in recombinant biopharmaceuticals in the range of 10 - 1,000 ppm was developed

#### INTRODUCTION

- Residual host cell proteins (HCPs) are commonly present in minute quantities in the final biopharmaceutical protein (typically expressed as ng/mg of recombinant protein or ppm concentrations)
- HCPs can potentially elicit an unpredictable immune response in patients
- Current analytical methods for measuring HCPs, (typically ELISA, gels, blots), are expensive, subjective, time-consuming to develop, and require identities of the prior knowledge about the contaminant proteins.

# WORKFLOW OVERVIEW

- Spike a mixture of 5 proteins in the concentration range of 10-1,000 ppm (ng/mg) in the protein biopharmaceutical
- Enzymatic digestion of the protein sample to peptides
- 2D-LC/MS<sup>E</sup> with IDENTITY<sup>E</sup> to discover contaminant proteins (HCPs)
- Label-free quantification using the best responding top 3 peptides from each protein for comparison against known concentrations of the spiked proteins
- VERIFY<sup>E</sup> for mining the MS<sup>E</sup> dataset for extracting appropriate MRMs for each HCP
- LC-MRM assay for targeted, high-throughput absolute quantification, after spiking the <sup>13</sup>C<sup>15</sup>N-isotopically labeled peptides in the protein biopharmaceutical

# **METHODS**

#### Sample Preparation

A chimeric anti-phosphotyrosine IgG1 monoclonal antibody (PTG1 mAb) was expressed in two different CHO cell lines and purified by Protein A chromatography using two different protocols following manufacturer recommendations (Protocol I and Protocol II). Four samples labeled A1, B1, A2 and B2 were expressed in DG-44 CHO cells, while two samples labeled C and D were expressed in CHO-S cells. Samples A1/A2, and B1/B2 were biological replicates, grown under identical conditions. Five protein standards (LA, PHO, ADH, BSA and ENL) were spiked in PTG1 preparations along with four <sup>13</sup>C<sup>15</sup>N isotopically labeled peptides and the resulting protein mixture was denatured with RapiGest (Waters), reduced/alkylated with DTT/IAM and enzymatically digested with trypsin (Promega). The amounts of spiked proteins loaded on-column were: 4,000 fmoles LA (bovine alpha-lactalbumin), 800 fmoles PHO (rabbit phosphorylase b), 320 fmoles ADH (yeast alcohol dehydrogenase), 80 fmoles BSA (bovine serum albumin) and 16 fmoles ENL (yeast enolase).

#### LC Conditions

A nanoACQUITY<sup>™</sup> UPLC<sup>®</sup> system (Waters) with 2D technology was used for peptide separations. A reversed-phase/reversed-phase (RP/ RP) method was developed that uses the pH of the mobile phases to change the selectivity of a peptide separation in two separate dimensions [1,2].

- First Dimension (1D) pH=10: 1.0 mm x 50 mm XBridge C<sub>18</sub> column (5 µm particles), 10 µL/min flow. Mobile phase: 20 mM ammonium formate in water (Solvent A) and ACN (Solvent B).
- **Online dilution (1:10)** of the eluent from 1D before analyte trapping onto the trap column.
- **Trap column**: A 5- $\mu$ m Symmetry C<sub>18</sub> trap (2.7  $\mu$ l volume) was used to trap peptides between the two LC dimensions.
- Second Dimension (2D) pH=2.4: 0.3 mm x 150 mm analytical column BEH C<sub>18</sub> 1.7  $\mu$ m, kept at 65 °C and operated at 12  $\mu$ L/min.

Fractions were eluted in ten steps (1: 10.8% Eluent B, 2: 12.4%, 3: 14.0%, 4: 15.4%, 5: 16.7%, 6: 18.6%, 7: 20.4%, 8: 25.0%,

9: 30.0%, 10: 50.0% B. Each step was mixed in a 1:10 ratio with 0.1% TFA in water (pH=2.1) before trapping. Low pH separations in the second chromatographic dimension used a 30 min gradient from 7 to 35% acetonitrile (0.1% FA-formic acid). See Figure 1 for the diagram of the 2D-LC system.

#### MS conditions

A) Discovery and Label-free Quantification: Data independent, alternate scanning LC/MS<sup>E</sup> experiments were performed with a SYNAPT G1 HDMS<sup>™</sup> mass spectrometer (Waters).

- Acquisition time was 0.5 sec MS and MS<sup>E</sup>, m/z range: 50-1990.
- Fixed CE at 5 eV for MS; CE ramp 15-35 eV for high-energy MS<sup>E</sup>.
- ProteinLynx Global Server (PLGS) 2.4 with Identity<sup>E</sup> Informatics[3].

**B)** Targeted Quantitation: MRM analyses were performed on two platforms: 1) Xevo<sup>™</sup> TQ MS instrument coupled with the ACQUITY UPLC system, 2) Xevo<sup>™</sup> TQ-S mass spectrometer coupled to a TRI-ZAIC nanoTile<sup>TM</sup> and a nanoACQUITY LC system.

- 2.1 x 150 mm BEH130 C<sub>18</sub> column (1.7 μm particles) operated at 35 °C with a flow rate of 300  $\mu$ L/min; 3-40% B over 15 min (B: acetonitrile with 0.1% FA; A: water with 0.1% FA).
- nanoTile<sup>™</sup> packed with BEH 1.7 µm C18 particles, operated at 45 °C using a flow rate of 3 µL/min; 3-40% B over 7 min, same mobile phases as above.
- Other parameters: ESI potential 3.5 kV, CV 37 V, source 90 °C, MS1/MS2 isolation window 0.75 Da (FWHM), adjustable collision energy depending of precursors/fragment (15-30 ev range) and 10 ms dwell time.

# TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS







Figure 1: Fluidic configuration for 2-dimensional chromatogra phy with on-line dilution: (A) Sample loading; (B) Peptide elution from 1D (in 10 steps) followed by trapping; (C) Peptide separation in the second dimension.



Figure 2. SYNAPT HDMS<sup>TM</sup> mass spectrometer coupled to the 2D nanoACQUITY UPLC<sup>®</sup> system.

# **RESULTS**

#### **HCP Identification Method**

- LC-MS<sup>E</sup> Data: a combined MS<sup>E</sup> dataset collected from 10-step 2D RP/RP separations performed in triplicate.
- Protein Identification Informatics: PLGS 2.4 featuring Identity<sup>E</sup> Bioinformatics.
- **Database:** protein database was compiled from 12,943 mouse Swiss Prot entries, 654 Golden hamster proteins, the sequences of LA, ADH, PHO, BSA, ENL, porcine trypsin, protein A (S aureus) and the heavy and light chain seguence of PTG1. The final randomized database contained 27,212 entries.



Table I. HCP concentrations (expressed in ppm or ng HCP/mg protein) measured across six PTG1 mAB preparations using the three best responding peptides in ESI-MS. Samples labeled A1, A2 and C were purified using **Protocol I**, while B1, B2 and D were processed using **Protocol II**. Samples A1/A2 and B1/B2 were biological replicates produced by DG-44 hamster cells, while samples C and D were expressed in CHO-S hamster cells grown under identical conditions. High abundance proteins are highlighted in RED (> 1,000 ppm), medium abundance are either in Y (500-1000 ppm) or GREEN (100-500 ppm) and low-abundance HCPs are highlighted in GREY (< 100 ppm). The mAb purity for each preparation is displayed at the bottom of each sample column.

Sample	e Sample	CHO Cell	Protein conc	No of HCPs	Common	HCP conc	PTG1
no	ID	Line	(mg/mL)	Identified	HCPs	(ppm)	purity (%)
1	A1	DG44	5.00	21	4	5945	94.05
2	A2	DG44	7.37	23	4	7190	92.81
3	B1	DG44	5.79	30	4	12897	87.10
4	B2	DG44	4.47	30	4	15302	84.30
5	C	CHO-S	3.88	5	4	1151	98.85
6	D	CHO-S	5.00	9	4	1020	98.98

Table II. Summary of HCP identifications provided by the MS<sup>E</sup> analysis.

# /oters THE SCIENCE OF WHAT'S POSSIBLE.™

Peptide	Protein	RT	MRM Transition		CV	CE
sequence	ID	(min)	Native	13C15N	(V)	(eV)
SLLNSLEEA <mark>K</mark>	Clusterin	12.2	552.3 -> 903.5	556.3 -> 911.5	35	21
STTTGHLIY <mark>K</mark>	Elongation factor	8.3	560.8 -> 831.5	564.8 -> 839.5	35	22
GAAQNIIPASTGAA <mark>K</mark>	pastgaa <mark>k</mark> GAD3-PDH		685.4 -> 702.4	689.4 -> 710.4	35	26
VLFDTGLVNP <mark>R</mark>	Nidogen-1	12.9	615.8 -> 1018.5	620.8 -> 1028.5	35	24

Table III. Experimental conditions for the MRM assay



Figure 3. MRM chromatograms recorded for four HCPs measured in sample B1. Each panel displays the best responding transitions for the native and the corresponding <sup>13</sup>C<sup>15</sup>N-isotopically labeled peptide belonging to (A,E) clusterin, (B,F) elongation factor 1-alpha, (C,G) glyceraldehyde 3-phosphate dehydrogenase and (D,H) nidogen-1. MRMs displayed in panels A-D were recorded using the analytical scale chromatography setup (2.1 mm column ID), while chromatograms presented in panels E-H were obtained on the TRIZAIC nanoTile.



Figure 4. Comparison between MS<sup>E</sup> and MRM quantification. TOF-based quantification  $(MS^{E})$  is comparing the precursor signals in ESI-MS generated by the three best ionizing peptides against the signals produced by a spiked protein with a known concentration. MRM quantification is based on using spiked <sup>13</sup>C<sup>15</sup>N-isotopically labeled peptides with known concentration. Protein concentrations (ppm) measured in six mAb preparations are shown for (I) clusterin, (II) elongation factor 1-alpha (III) glyceraldehyde 3-phosphate dehydrogenase and (IV) nidogen-1.

# CONCLUSIONS

- The LC/MS<sup>E</sup>-based assay allows the identification and quantification of low-abundance HCP contaminants in biopharmaceuticals over four orders of magnitude in concentration.
- TRIZAIC nanoTile in combination with the TQ-S tandem quadrupole shows great promise for high-throughput HCP quantification or monitoring (under 10 min runtime).
- Protein A purification of mAbs using different purification protocols produces different HCP patterns.

#### References

- 1. Gilar M, Olivova P, Daly AE, Gebler JC J. Sep. Sci, 2005, 1694.
- 2. Stapels MD, Fadgen K Current Trends in Mass Spectrometry, Spectroscopy supplement, March 2009.
- 3. Li GZ, Vissers JP, Silva JC, Golick D, Gorenstein MV, Geromanos SJ Proteomics, 2009, 1696.