MaxEnt: An Advanced Maximum Entropy Based Tool for Disentangling the Electrospray Mass Spectra from Biopolymer Mixtures. A Brief Description and Bibliography.

Authors: B.N. Green and R.W.A. Oliver

Introduction

The electrospray mass spectrum from a biopolymer (e.g. protein, glycoprotein, or oligonucleotide) is composed of a series of multiply-charged ions on the mass-to-charge ratio (m/z) scale. On this scale, the positions of the multiply-charged positive ions from a single protein or glycoprotein of molecular mass M_r are given by $m/z = (M_r+nH)/n$, where H is the mass of the proton and n is an integer in a series of consecutive integers. Typically, a 20 kDa protein generates a series containing some 10-20 multiply-charged molecular ion peaks. It follows that the electrospray mass spectrum from a mixture of proteins can be extremely complex, since each protein in the mixture will give rise to its own characteristic series of such peaks. Simplification of the raw data is thus mandatory before interpretation of the complex spectrum can be attempted. If each series from a given component of the mixture could be condensed into a single peak on a true molecular mass scale, the amount of data to be considered would be significantly reduced and the process of interpretation could begin.

Various methods have been described in the literature for carrying out this deconvolution procedure and their drawbacks discussed⁽¹¹⁾. In 1991, an advanced maximum entropy (MaxEnt) based procedure was developed by J. Skilling of MaxEnt Solutions Ltd, Cambridge, UK^(1, 3). Subsequently, this MemSys5 maximum entropy algorithm was incorporated into the Micromass MassLynx





software package and has been employed in this form to solve many problems in biochemistry.

Here, the benefits of MaxEnt over alternative methods of processing electrospray mass spectra, particularly from biopolymer mixtures, are first summarised. This summary is followed by a bibliography of over 50 papers in which MaxEnt has been employed to deconvolute the electrospray mass spectra from a wide range of biological samples.

The Benefits of MaxEnt

MaxEnt has the following distinct advantages over alternative methods of processing multiply-charged electrospray data:-

 MaxEnt automatically finds the molecular masses of the components present in a mixture.

With a completely unknown sample, a first MaxEnt analysis or zero-charge survey is made over a deliberately chosen wide molecular mass range in order to localise the mass range or ranges for a second definitive analysis. MaxEnt is computationally intensive and the processing time increases with the number of data points (channels) in the zerocharge (output) spectrum. Consequently, when making a wide mass range survey, it is usual to employ a relatively coarse mass scale (5-10 Da/channel) to keep the processing time to within 5 minutes or so. However, when this is done, the molecular weights are not calculated with full accuracy. Therefore, once the approximate masses of the components have been found from the survey, a second definitive run is made over a narrower output mass range (or ranges) using a finer mass scale (0.5-1 Da/channel) in order to obtain fully accurate molecular mass values.

An example of the latter is illustrated in Fig.1. Here, the original spectrum from the mixture of the proteins and glycoproteins obtained by denaturing the extracellular haemoglobin from the earthworm, *Lumbricus terrestris*^(26, 32) is shown in Fig. 1A. The native Hb has a molecular weight of 3.5 million and is composed of some 15 different proteins and glycoproteins ranging in mass from 16 to 53 kDa. The MaxEnt deconvoluted spectrum shown in Fig. 1B, where each component in the mixture is represented as a single peak on a true molecular mass scale, is much easier to interpret than the original data. These data were acquired on a Micromass Quattro II.



Figure 1. Electrospray mass spectra of the globins and subunits from the haemoglobin of the earthworm Lumbricus terrestris analysed in denaturing solvent. A, the original m/z spectrum and B, the MaxEnt processed spectrum on a true molecular mass scale. d are monomer globins and T glycosylated disulphide bonded trimers. Linker subunits (L), of which L1 is also glycosylated, are necessary for assembly of the various proteins into the native haemoglobin of molecular mass 3.5 million. In Fig. 1A, the figures after the commas indicate the number of charges on the ions, e.g. d1,13 means component d1 with 13 positive charges. (Copyright 1996 American Chemical Society (32)).

 MaxEnt enhances the resolution i.e. improves the ability to separate and accurately measure the molecular masses of otherwise unresolved components.

The arguments for needing resolution enhancement in the analysis of haemoglobins by electrospray mass spectrometry have been described in an excellent review by Shackleton *et al.*⁽⁸⁾. In essence, there is a basic limit to peak width, and hence resolution, which is determined by the isotopic distribution of the elements in the protein molecule. For the haemoglobin chains (15-16 kDa), this basic width at half peak height is 8 Da. With state-of-the-art quadrupole analysers, the instrumental contribution to peak width is small for the haemoglobin chains, and increases the overall width to a little less than 10 Da, thus giving rise to the practical limit for resolving two haemoglobin chains, without deconvolution, of 12 Da. The Shackleton review (pp 148-151) shows how MaxEnt extends the practical limit for resolving and accurately measuring the masses of two haemoglobin chains from 12 to 6 Da.

Fig. 2 shows an example of resolution enhancement by MaxEnt from the analysis of a heterozygote for the α -chain variant Le Lamentin, $[\alpha_{20} \text{ (His} \rightarrow \text{Gln})]$. Here, the variant α -chain (Mr=15117.4) is not resolved from the normal α -chain (Mr=15126.4) in the original data since the mass difference is only 9 Da (Fig. 2A and inset). After deconvolution by MaxEnt (Fig. 2B), the two α -chains are clearly resolved allowing their mass difference to be accurately determined.

3. MaxEnt improves the signal-to-noise ratio.

MaxEnt has the power to extract useful zerocharge spectra from noisy multiply-charged m/z data.

Although the majority of publications in the Bibliography involve the analysis of proteins, it should be stressed that MaxEnt is not restricted to processing biopolymer data but can be used to process any data containing two or more ions from the same molecule, provided they carry consecutive numbers of charges. Fig. 3 shows the result of employing MaxEnt to produce the isotopically resolved zero-charge spectrum of a poly-sulphonated compound (Mr=1815) from partly resolved multiply-charged data. The original spectrum (Fig. 3A and inset) shows predominantly an ion with 6 negative charges together with less intense ions having 5 and 4 charges. The MaxEnt processed spectrum (Fig. 3B and insets) shows the zero-charge spectrum with the isotope peaks almost fully resolved. Moreover, the isotope pattern closely resembles the expected isotope pattern for this molecule of elemental composition C₆₁H₄₀N₁₉O₂₆S₈Cl₃.







Figure 3. Electrospray data from a poly-sulphonated compound of molecular mass 1815. A, the original m/z spectrum and B, the zero charge spectrum produced by MaxEnt. Note that in the original data (Fig. 3A inset), the peak with six negative charges is only partly isotopically resolved, whereas in the MaxEnt spectrum, the isotope peaks are clearly resolved (Fig. 3B inset). Furthermore, the experimental isotope pattern closely resembles the isotope pattern expected for this molecule.

The following is a Bibliography of published papers in which MaxEnt has been used to deconvolute electrospray ionisation mass spectra.

- Maximum entropy deconvolution in electrospray mass spectrometry. [Human cord blood haemoglobin] Ferrige, A.G., Seddon, M.J. and Jarvis, S. Rapid Commun. Mass Spectrom. 5, 374-379 (1991).
- Quantified maximum entropy. [Animal haemoglobins]. Skilling, J. Am. Lab. 24, 32J-32M (1992).
- Disentangling electrospray spectra with maximum entropy. [Human haemoglobins]. Ferrige, A.G., Seddon, M.J. Green, B.N., Jarvis, S.A., Skilling, J. Rapid Commun. Mass Spectrom. 6, 707-711 (1992).
- The application of "MaxEnt" to high resolution mass spectrometry. [Low molecular mass organic compound]. Ferrige, A.G., Seddon, M.J., Skilling, J., Ordsmith, N. Rapid Commun. Mass Spectrom. 6, 765-770 (1992).
- Application of electrospray ionisation mass spectrometry with maximum entropy analysis to allellic fingerprinting of major urinary proteins [of mice]. Evershed, R.P., Robertson, D.H.L., Beynon, R.J., Green, B.N. Rapid Commun. Mass Spectrom. 7, 882-886 (1993).
- Maximum entropy deconvolution of heterogeneity in protein modification; protein [β-lactoglobulin B] adducts of 4-hydroxy-2-nonenal. Bruenner, B.A., Jones, A.D., German, J.B. Rapid Commun. Mass Spectrom. 8, 509-612 (1994).
- Glycopeptide [Horseradish peroxidase] mapping and structural analysis by electrospray mass spectrometry.Hutton, T., Green, B.N., Geisow, M.J. Peptides ; Chemical Structure and Biology. Proceedings of the 13th American Peptide Symposium,1993.

Edited by Hodges, R.S., Smith, J.A. ESCOM, Leiden pp 329-331, (1994)

- Mass spectrometry in the characterisation of variant hemoglobins. Shackleton, C.H.L., Witkowska, H.E. In Mass Spectrometry: Clinical and Biomedical Applications. Vol. 2, Edited by Desiderio, D.M., Plenum Press, New York, pp 135-199 (1994).
- Incorporation of norleucine at methionine positions in recombinant human macrophage colony stimulating factor (rM-CSF, 4- 153) expressed in *Esherichia coli*; Structural analysis. Randhawa, Z.I., Witkowska, H.E., Cone, J., Wilkins, J.A., Hughes, P., Yamanashi, K., Yasuda,S., Masui, Y., Arthur, P., Kletke, C., Bitsch, F., Shackleton, C.H.L. Biochemistry 33, 4352-4362 (1994).
- Modification of glutathione S-transferase 3-3 mutants with 2-(S-glutathionyl)-3,5,6trichloro-1,4-benzoquinone.
 Hong, J.L., Liu, L.F., Wang, L.Y., Tsai, S.P., Hsieh, C.H., Hsiao, C.D., Tam, M.F. Biochem. J. 304, 825-831 (1994).
- MaxEnt: An essential maximum entropy based tool for interpreting multiplycharged electrospray data. [human, animal and leech haemoglobins, Fc fragment from IgG; methane mono-oxygenase A]. Cottrell, J.C., Green, B.N. Fisons Instruments Application Note No 212, pp15 (1994).
- The cyclic structure of the enterococcal peptide antibiotic AS-48. Samyn, B., Martinez-Bueno, M., Devreese, B., Maqueda, M., Galvez, A., Valdivia, E., Coyette, J., Van Beeumen, J. FEBS Letts. 352, 87-90 (1994).
- Electrospray mass spectrometric characterization of the components of protein mixtures and its application to members of the chemokine family of interleukins. Kay, I., Hafizi, S., Barr, R.M., Mallet, A.I., Greaves, M.W. Biol. Mass Spectrom. 23, 267-271 (1994)

- Electrospray ionization mass spectrometric determination of the complete polypeptide composition of *Tylorrhynchus heterochaetus* hemoglobin. Green, B.N., Suzuki, T., Gotoh, T., Kuchumov, A.R., Vinogradov, S.N. J. Biol. Chem. 270, 18209-18211 (1995).
- Mass spectrometric composition, molecular mass and oxygen binding of *Macrobdella decora* hemoglobin and its tetramer and monomer subunits.
 Weber, R.E., Malte, H., Braswell, E.H., Oliver, R.W.A., Green, B.N., Sharma, P.K., Kuchumov, A., Vinogradov, S.N. J. Mol. Biol. 251, 703-720 (1995).
- Electrospray mass spectrometric analysis of resorcinol suicide substrate binding to cytochrome C peroxidase. Doerge, D.R., Divi, R.L., Hutton, T., Lowes, S. J. Mass Spectrom. (Special Issue) S136-S142 (1995).
- Electrospray ionisation mass spectrometry of lens crystallins ; Verification and detection of errors in protein sequences of bovine gamma crystallins. Kilby, G.W., Truscott, J.W., Aquilina, J.A., Sheil, M.M., Riley, M.L., Harding, J.J. Eur. Mass Spectrom. 1, 203-208 (1995).
- Purification and characterization of the imidazoleglycerol-phosphate dehydratase of *Saccharomyces cerevisiae* from recombinant *Escherichia coli*. Hawkes, T.R., Thomas, P.G., Edwards, L.S. Rayner, S.J., Wilkinson, K.W., Rice, D.W. Biochem. J. 306, 385-397 (1995).
- Some effects of post-translational N-terminal acetylation of the human embryonic ζ globin protein. Scheepens, A., Mould, R., Hofmann, O., Brittain, T. Biochem. J. 310, 597-600 (1995).
- Mass spectrometric analysis of rat liver cytosolic glutathione S-transferases : Modifications are limited to N-terminal processing. Yeh, H.I., Hsieh, C.H., Wang, L.Y., Tsai, S.P., Hsu, H.Y. Tam, M.F. Biochem. J. 308, 69-75 (1995).

- Purification, mass spectrometric characterisation, and covalent modification of murine glutathione S-transferases. Mitchell, A.E., Morin, D., Lame, M.W., Jones, A.D. Chem. Res. Toxicol. 8, 1054-1062 (1995).
- 22. An anticodon sequence mutant of *Escherichia coli* initiator tRNA : Possible importance of a newly acquired base modification next to the anticodon on its activity in initiation. Mangroo, D., Limbach, P.A., McCloskey, J.A., Rajbhandary, U.L. J.Bacteriol. 177, 2858-2862 (1995).
- The use of electrospray mass spectrometry to identify an essential arginine residue in type II dehydroquinases. Krell, T., Pitt, A.R., Coggins, J.R. FEBS Letts. 360, 93-96 (1995).
- 24. The use of mass spectrometry to examine the formation and hydrolysis of the phosphorylated form of phosphoglycerate mutase. Nairn, J., Krell, T., Coggins, J.R., Pitt, A.R., Fothergill-Gilmore, L.A., Walter, R., Price, N.C. FEBS Letts. 359, 192-194 (1995).
- Analysis of the ion binding sites of calmodulin by electrospray ionization mass spectrometry. Lafitte, D., Capony, J.P., Grassy, G., Haiech, J., Calas, B. Biochemistry 34, 13825-13832 (1995).
- Mass spectrometric composition and molecular mass of *Lumbricus terrestris* hemoglobin: A refined model of its quaternary structure. Martin, P.D., Kuchumov, A.R., Green, B.N., Oliver, R.W.A., Braswell, E.H., Emory, H., Wall, J.S., Vinogradov, S.N. J. Mol. Biol. 255, 154-169 (1996).
- 27. Posttranslational processing of recombinant human interferon-γ in animal expression systems.
 James, D.C., Goldman, M.H., Hoare, M., Jenkins, N., Oliver, R.W.A., Green, B.N., Freedman, R.H.
 - Protein Sci. 5, 331-340 (1996).

- 28. The assignment of charge states in complex electrospray mass spectra. [Monoclonal antibody components; *Tylorrhynchus heterochaetus* haemoglobin]. Almudaris, A., Ashton, D.S., Beddell, C.R., Cooper, D.J., Craig, S.J., Oliver, R.W.A. Eur. Mass Spectrom. 2, 57-67 (1996).
- 29. The multi-hemoglobin system of the hydrothermal vent tube worm *Riftia pachyptila*. II. The complete polypeptide chain composition investigated by maximum entropy analysis of mass spectra. Zal, F., Lallier, F.H., Green, B.N., Vinogradov, S.N., Toulmond, A. J. Biol. Chem. 271, 8875-8881 (1996).
- Electrospray ionization mass spectrometric composition of the 400 kDa hemoglobin from the pogonophoran *Oligobrachia mashikoi* and the primary structures of three major globin chains. Yuasa, H.J., Green, B.N., Takagi, T., Suzuki, N., Vinogradov, S.N., Suzuki, T. Biochim. Biophys. Acta. 1296, 235-244 (1996).
- 31. Analysis of complex protein and glycoprotein mixtures by electrospray ionization mass spectrometry with maximum entropy processing.
 [Components of *Tylorrhynchus heterochaetus* and *Lumbricus terrestris* haemoglobins]. Green, B.N., Hutton, T., Vinogradov, S.N. In: Methods in molecular biology, vol 61: Protein and peptide analysis by mass spectrometry. Edited by Chapman, J.R. Humana Press Inc., Totowa, NJ, pp 279-294 (1996).
- Giant hexagonal bilayer hemoglobins. Lamy, J.N., Green. B.N., Toulmond, A., Wall, J.S., Weber, R.E., Vinogradov, S.N. Fig. 1 reprinted with permission from Chem. Rev. 96, 3113-3124 (1996). Copyright 1996 American Chemical Society
- Electrospray mass spectrometric evidence for the occurrence of two major variants in native pig pepsin A. Green, B.N., Jones, A.T., Roberts, N.B.Biochem. J. 313, 241-244 (1996).

- 34. Overexpression, purification, and refolding of link module from human TSG-6 in *Escherichia coli:* effect of temperature, media, and mutagenesis on lysine misincorporation at arginine AGA codons. Day, A.J., Aplin, R.T., Willis, A.C. Protein Expression Purification. 8, 1-16 (1996).
- 35. The use of electrospray ionization MS to determine the structure of glycans in intact glycoproteins. [Recombinant γ-interferon; avidin from chicken egg white].
 Oliver, R.W.A., Green, B.N., Harvey, D.J. Biochem. Soc. Trans. 24, 917-927 (1996).
- Fluphenazine photoaffinity labelling of binding sites for phenothiazine inhibitors of trypanothione reductase. Yon, H., Chan, C., Garforth,J., Douglas, K.T., Bolgar, M.S. Gaskell, S.J., Fairlamb, A.H. Chem. Commun. 6, 973-974 (1996).
- 37. Expression of soluble human β-globin chains in bacteria and assembly in vitro with α-globin chains.
 Yamaguchi, T., Pang, J., Reddy, K.S., Witkowska, H.E., Surrey, S., Adachi, K. J. Biol. Chem. 271, 26677-26683 (1996).
- Inhibition of TEM-2-β-lactamase from *Escherichia coli* by clavulanic acid: Observation of intermediates by electrospray ionization mass spectrometry. Brown, R.P.A., Aplin, R.T., Schofield, C.J.
- Biochemistry 35, 12421-12432 (1996). 39. Globin and globin gene structure of the
- Dewilde, S., Blaxter, M., Van Hauwaert, M.L., Vanfleteren, J., Esmans, E.L., Marden, M., Griffon, N., Moens, L. J. Biol. Chem. 271, 19865-19870 (1996).
- Haloenol lactone is a new isozymeselective and active site-directed inactivator of glutathione S-transferase. Zheng, J., Mitchell, A.E., Jones, A.D., Hammock, B.D.
 - J. Biol. Chem. 271, 20421-20425 (1996).

- 41. Probing the nature of noncovalent interactions by mass spectrometry. A study of protein-CoA ligand binding and assembly. Robinson, C.V., Chung, E.W., Kragelund, B.B., Knudsen, J., Aplin, R.T., Poulsen, F.M., Dobson, C.M. J. Am. Chem. Soc. 118, 8646-8653 (1996).
- 42. Localization of the active site of type II dehydroquinases.
 Krell, T., Horsburgh, M.J., Cooper, A., Kelly, S,M., Coggins, J.R.
 J. Biol. Chem. 271, 24492-24497 (1996).
- 43. Contributions of asparagine at α97 to the cooperative oxygenation process of hemoglobin [Recombinant modified human hemoglobin].

Kim, H.W., Shen, T.J., Ho, N.T., Zou, M., Tam, M.F., Ho, C.

Biochemistry 35, 6620-6627 (1996).

44. Rat kidney glutathione S-transferase 1 subunits have C-terminal truncations. Yeh, H.I., Lee, J.Y., Tsai, S.P., Hsieh, C.H., Tam, M.F. Biochem. J. 314, 1017-1025 (1996).

 Characterizing abnormal hemoglobin by MS.
 Shackleton, C.H.L. and Witkowska, H.E.

Anal. Chem. 68, A29-A35 (1996).

- 46. A novel silent postranslational mechanism converts methionine to aspartate in hemoglobin Bristol (β67[E11] Val-Met→Asp). Rees, D.C., Rochette, J., Schofield, C., Green, B.N., Morris, M., Parker, N.E., Sasaki, H., Tanaka, A., Ohba, Y., Clegg, J.B. Blood 88, 341-348 (1996).
- Determination of the sites of 4-hydroxy-2-nonenal adduction to protein by electrospray tandem mass spectrometry. Bolgar, M.S., Gaskell, S.J. Anal. Chem. 68, 2325-2330 (1996).

- 48. Rapid characterisation of chemically modified proteins by electrospray mass spectrometry. [Hen egg white lysozyme complexes with biotin, iodine and Co(diAMsar); murine B72.3 monoclonal antibody, its reduced Fc fragment, its Fab and its biotinylated Fab fragment] Bennett, K.L., Smith, S.V., Lambrecht, R.M., Truscott, R.J.W., Sheil, M.M. Bioconjugate Chem. 7, 16-18 (1996).
- Identification of five hemoglobins in B6C3F1 mice by mass spectrometry and sequence analysis. Palm, L., Burka, L.T., Hojrup, P., Stevens, R.D., Tomer, K.B. Intl. J. Biochem. Cell Biol. 28, 1319-1326 (1996).
- Identification and purification of a distinct dihydrolipoamide dehydrogenase from pea chloroplasts. Conner, M., Krell, T., Lindsay, J.G. Planta 200, 195-202 (1996).
- Phosphoglycerate mutase from Schizosaccharomyces pombe: development of an expression system and characterisation of three histidine mutants of the enzyme. Nairn, J., Price, N.C., Kelly, S.M., Rigden, D., Fothergill-Gilmore, L.A., Krell, T. Biochim. Biophys. Acta 1296, 69-75 (1996).
- Significant hydrogen exchange protection in GroEL-bound DHFR is maintained during iterative rounds of substrate cycling. Grob, M., Robinson, C.V., Mayhew, M., Hartl, F.U., Radford, S.E. Protein Science 5, 2506-2513 (1996).
- 53. Hemoglobin Old Dominion/ Burton-upon- Trent, β143[H21]His→Tyr, Codon 143CAC→TAC. A variant with altered oxygen affinity that compromises measurement of Hb A1c in diabetes mellitus; structure, function, DNA sequence. Elder, G.E. *et al.* Mayo Proceedings (1997) (in press).

- 54. Quaternary structure of the extracellular haemoglobin of the lugworm *Arenicola marina*. A multi-angle-laser-lightscattering and electrospray-ionisationmass-spectrometry analysis. Zal, F., Green, B.N., Lallier, F.H., Vinogradov, S.N., Toulmond, A. Eur. J. Biochem. 243, 85-92 (1997).
- 55. Purification of the T-cell receptor ζ-chain. Covalent modification by
 4-(2-aminoethyl)-benzenesulfonyl fluoride. Sweeney, B., Proudfoot, K., Parton, A.H., King, L.M., Slocombe, P., Perry, M.J. Anal. Biochem. 245, 107-109 (1997).
- Monitoring papain digestion of a monoclonal antibody by electrospray ionization mass spectrometry.
 Bennett, K.L., Smith, S.V., Truscott, R.J.W., Sheil, M.M.
 Anal. Biochem. 245, 17-27 (1997).

Discussion of Bibliography

The wide range and exact nature of the biological mixtures analysed in which MaxEnt procedures have been employed in order to deconvolute their electrospray mass spectra can be readily ascertained by a study of the titles of the papers given in the bibliography. If this information was not given by the original authors of the paper then it has been added to the title in square brackets. At the present time, the most common mixtures studied by this method are those of various haemoglobins, possibly because, as Shackleton and Witkowska comment in paper 45, these proteins afford "dream molecules" for the mass spectrometrist. A study of the adjacent Index of Journals shows that MaxEnt applications have been published in a very wide range of journals, some 27 in total, with the two most popular journals being the two most important biochemical journals, namely the Journal of Biological Chemistry and the Biochemical Journal.

Finally, a study of the country of origin of the various papers forming the bibliography showed that at present MaxEnt is most widely used in the UK (39 papers) followed by American laboratories (18 papers). This study also showed that papers have been published from laboratories in Australia (3 papers), Belgium (2 papers), Denmark (3 papers), France (4 papers), Japan (3 papers), Taiwan (4 papers) and Spain (1 paper) and hence it may be safely concluded that the MaxEnt procedures have been adopted worldwide.

Acknowledgements.

The authors would like to thank Dr T. M. Reynolds, Clinical Chemistry Dept., Burton Hospital, Burton upon Trent, Staffordshire for permission to use the Hb Le Lamentin data and BASF plc., P.O. Box 4, Hexagon House, Blackley, Manchester for permission to use the poly-sulphonated compound data.

INDEA OF JOURN	ALS
Analytical Biochemistry	55, 50.
Analytical Chemistry	45, 47.
American Laboratory	2.
Biochemical Journal	10, 18, 19, 20, 33, 44.
Biochemical Society Transactions	35.
Biochemistry	9, 25, 38, 43.
Biochimica et Biophysica Acta	30, 51.
Bioconjugate Chemistry	48.
Biological Mass Spectrometry	13.
Blood	46.
Books and [application note]	7, 8, 31, [11].
Chemical Communications	36.
Chemical Research and Toxicology	21.
Chemical Reviews	32.
European Journal of Biochemistry	54.
European Mass Spectrometry	17, 28.
FEBS Letters	12, 23, 24.
International Journal Biochemistry	49.
and Cell Biology	
Journal of the American Chemical Society	41.
Journal of Bacteriology	22.
Journal of Biological Chemistry	14, 29, 37, 39, 40, 42.
Journal of Mass Spectrometry	16.
Journal of Molecular Biology	15, 26.
Mayo Proceedings	53.
Planta	50.
Protein Expression and Purification	34.
Protein Science	27, 52.
Rapid Communications in Mass Spectrometry	1, 3, 4, 5, 6.

MICROMASS UK Limited

Floats Road Wythenshawe Manchester M23 9LZ Tel: + 44 (0) 161 945 4170 Fax: + 44 (0) 161 998 8915 http://www.micromass.co.uk

UK Sales Desk, Tel: 0161 946 0565 USA Beverly MA, Tel: 978 524-8200 Canada Pte-Claire, Tel: 514 694-1200 EUROPE Weesp, Tel: 431 (0)294-480484 Netherlands Weesp, Tel: 0294-480484 Belgium Vilvoorde, Tel: 02-2534550 France Villeurbanne, Tel: 02-2534550 Germany Idstein, Tel: 0612-691561 Sweden Täby, Tel: 08-7567130 Switzerland Sissach, Tel: 61 973 0800



Certificate No 951387 TECHNICAL NOTE NO. TN202/ICA Version 1 The product information in this document was correct at the time of printing, however specifications are subject to alteration without notice due to the Micromass UK Limited philosophy of continual product innovation and development and any information contained herein should not be construed as a contractual commitment by Micromass UK Limited. All orders are accepted subject to our current Conditions of Sale, copies of which are available on request.

© Micromass Limited May 1998

