


MS Nose™ Atmospheric Pressure Chemical Ionisation Gas Phase Analyser (APcI-GPA)

Authors: R. Linforth, Department of Applied Biochemistry & Food Science, University of Nottingham, Leicestershire LE12 5RD UK
S. Preece, Micromass UK Ltd., Floats Road, Wythenshawe, Manchester M23 9LZ UK

Introduction

The Micromass MS Nose™ ApcI-GPA is a new interface* for a mass spectrometer which has been designed for the analysis of volatile compounds. The low dead volume, inert sample introduction system, delivers the volatile compounds into the source where they are ionised at atmospheric pressure by a high voltage corona discharge. The ions are extracted from the atmospheric pressure region through a sampling orifice (sample cone) into an intermediate vacuum region and, from there, they are transported into the analyser region of the mass spectrometer where they are separated on the basis of their mass to charge ratio. APcI-GPA is a powerful tool enabling the analyst to monitor both the profile of compounds in the gas phase, and changes in that profile in real time.

Suitable for a range of applications, the system can be used to analyse volatile compounds from a wide variety of sources including human breath, perfumes, headspace autosamplers etc.



*European Patent #19937




Figure 1: Schematic of the APcI-GPA interface

Key Features:

- **Temporal Resolution** - The low dead volume of the interface coupled with the rapid response rates of mass spectrometers allows rapid changes in volatile concentration to be monitored.
- **Sensitivity** - APcI ionises a wide range of compounds (acids, alcohols, esters, ketones, pyrazines etc) which can typically be detected at 10ppbv in the gas phase.
- **Simplicity** - APcI is a soft ionisation technique giving minimal fragmentation. Compounds typically form the protonated molecule, $[M+H]^+$, simplifying ion profiles and making data interpretation easier
- **Reduced sample preparation** - With APcI-GPA, gas phase samples can be easily introduced without oxygen or moisture removal, minimising sample work-up.
- **The Venturi** - Changing the venturi region adapter sets the range of gas sampling flow rate available. The venturi adjustment alters the position of the silica in the venturi region allowing fine adjustment of the sampling rate.
- **Flexibility** - The APcI-GPA can be used to monitor both real-time temporal changes in volatile concentration and the volatile profile at a specific point in time.
- **Transfer line** - The flexible heated transfer line links the gas phase sample to the interface and the mass spectrometer.
- **Quantitation** - The absolute concentration of compounds in the gas phase can be determined by comparing the signal intensities of compounds in the sample with those of standards introduced via the calibration inlet.
- **MS analysers** - The APcI-GPA interface can be fitted to a wide variety of Micromass® mass spectrometers providing a range of detection modes. The interface may be used with quadrupole instruments, or it may use the enhanced sensitivity and mass accuracy of a time-of-flight instrument or the added specificity of an MS-MS instrument.

Example Application: The Aroma of Food

The aroma compounds present in foodstuffs are key characteristic components of quality. Their release during eating and their transport to the nose affect our perception of food and the delivery of the correct aroma profile to the consumer is of major importance to the food industry. The compounds present in exhaled breath are present at low concentrations and can change rapidly, even over the course of a single exhalation. The APcI-GPA, with its high sensitivity and rapid response rates, has been used to follow the concentration of aroma compounds in the exhaled breath of people eating a range of foods.

The examples shown here are for the analysis of apples and bananas. APcI-GPA can be used for rapid sample aroma profiling - samples can simply be eaten and, whilst chewing, the consumer breathes in and out through an adapter through which the aroma compounds present in the breath can be monitored in real time. Figure 2 shows the intensity and rate of release of different compounds. Differences in flavour of the two fruits are due to the different compounds which are released during eating. For example, bananas contain isoamyl acetate (m/z 131) and ethyl acetate (m/z 89) which are not present in apples. Also, eating apples produces more (E)-2-hexenal (m/z 99) than bananas. The (E)-2-hexenal is produced by the action of enzymes released from the fruit when it is macerated, or chewed, which accounts for the profile for (E)-2-hexenal starting at a later retention time than the other compounds in the breath. The volatile compounds in the fruit can also be monitored by simply monitoring the headspace (Figure 3) for a few seconds to produce an aroma-ion profile.

For manufactured foodstuffs the effect of product formulation is often important (effect of fat content etc), yet this can be difficult to study in vitro. Measuring the aroma profile in the breath of individuals eating the food may show these differences in formulation and can be used in the development of processed foods.




Figure 2: The temporal profiles of ethyl hexanoate (m/z 145), isoamyl acetate (m/z 131), ethyl butyrate (m/z 117), (E)-2-hexenal (m/z 99) and ethyl acetate (m/z 89) in the breath of an individual eating apples and bananas. The ethyl hexanoate and ethyl butyrate were present in the apple before consumption and were released earlier than the (E)-2-hexenal which was generated by enzymes during eating, causing a lag. Isoamyl acetate and ethyl acetate are characteristic compounds which are present in the banana but not in the apple.




Figure 3: The aroma-ion profile of the headspace above a crushed apple and banana. The major ions (m/z 89, 117, 131, 145, 159) correspond to the masses of esters typical components of fruit aromas.

Relevant Publications

Time course profiling of volatile release from foods during the eating process

RST Linforth, KE Ingham & AJ Taylor in "Flavour Science: Recent Developments", Royal Society of Chemistry, pp361-368,

ISBN 0 85404 702 6 Eds AJ Taylor & DS Mottram 1996

Measuring flavour release from chewing gum

B Harvey

Kennedys Confection May 1997 pp18-20

Going Physical

L Tuley

Food Manufacture July 1997 pp20-22

Aromes mesures

Pour la Science No 241, p19 Novembre 1997,

H.This

Flavour release from foods: Recent developments

AJ Taylor & RST Linforth

Proceedings of the 5th Wartburg Aroma Symposium

Eds J-P Kruse & M Rothe pp131-142 1997

Relationship between sensory Time-Intensity

measurements and in-nose concentrations of volatiles

AJ Taylor & RST Linforth

in Flavour and Fragrances Ed KAD Swift

pp 171-182 1997 Royal Society of Chemistry

ISBN 0-85404-787-5

Dynamism in flavour science and sensory

methodology in "Interaction of food matrix with

small ligands influencing flavour and texture"

Proceedings of the COST96

meeting Valencia, 1998, Ed J Bakker,

European Commission pp 59-65

G Dijksterhuis, AJ Taylor, WE Brown, JR Piggott,

R Solheim, C Auloge, J Velisek, H Valentova,

S Skronakova, J Pokorny

Volatile release from food in the mouth during the

eating process in "Interaction of food matrix with

small ligands influencing flavour and texture"

Proceedings of the COST96 meeting Valencia, 1998,

Ed J Bakker, European Commission pp 78-81

RST Linforth, AJ Taylor, WE Brown

constituents in gases

European Patent Application EP 0819 937 A2 1998

RST Linforth & AJ Taylor

Volatile release from mint flavored sweets

Perfumer and Flavorist 23, 47-53 1998

RST Linforth & AJ Taylor

Effect of variety, time of eating and fruit to fruit

variation on volatile release during eating of tomato

fruits (*Lycopersicon esculentum*)

J Agric Food Chem 46, 2287-2292 1998

MS Brauss, RST Linforth, AJ Taylor

Measuring flavour release: an additional tool for

flavour development

Food Ingredients and Analysis May/June 1998

pp20-22

AJ Taylor & RST Linforth

Understanding flavour release: the key to better food flavour?

Nutrition & Food Science in press 1998

AJ Taylor & RST Linforth

Physical chemistry of flavour

International Journal Food Science and Technology

AJ Taylor in press Feb 98

Sensory perception is related to the rate of change of volatile concentration in-nose during eating of model gels

Chemical Senses submitted June 98

IN Baek, RST Linforth, A Blake & AJ Taylor

Altering the fat content in a model yogurt system;

Instrumental and sensory effects

Journal of Food Science submitted May 98

MS Brauss, RST Linforth, I Cayeux,

B Harvey & AJ Taylor

Simultaneous instrumental and sensory analysis of

volatile release from gelatine and pectin/gelatine gels

Food Chemistry submitted April 98

RST Linforth, I Baek & AJ Taylor

In-mouth measurement of pH and conductivity during eating

J Agric Food Chem submitted June 98

JM Davidson, RST Linforth & AJ Taylor

Orange volatiles: release from fruit and juice during consumption

J Agric Food Chem submitted July 98

A Hinterholzer, IN Baek, RST Linforth, AJ Taylor &

P Schieberle

MICROMASS UK Limited

Floats Road

Wythenshawe

Manchester M23 9LZ

Tel: +44 (0)161 945 4170

Fax: +44 (0)161 998 8915

<http://www.micromass.co.uk>

UK Sales Desk, Tel: 0161 946 0565

USA Beverly MA, Tel: 978 524-8200

Canada Pointe Claire, Tel: 514 694-1200

EUROPE Weesp, Tel: +31 (0)294-480484

Netherlands Weesp, Tel: 0294-480484

Austria Achau, Tel: 0223 672668

Belgium Vilvoorde, Tel: 02-2534550

France Villeurbanne, Tel: 04 72 14 89 00

Germany Eschborn, Tel: 0800 181 7249

Spain Barcelona, Tel: 34 932 466 696

Sweden Täby, Tel: 08-756 7130

Switzerland Sissach, Tel: 061 973 0800



Certificate No 951387

Micromass is a registered trademark of

Micromass Limited (Reg. U.S. Pat. & Tm. Off.).

APPLICATION NOTE No. AN 301/ICA VERSION 1

The product information in this document was correct at the time of printing, however specifications are subject to alteration without notice due to the Micromass UK Limited philosophy of continual product innovation and development and any information contained herein should not be construed as a contractual commitment by Micromass UK Limited. All orders are accepted subject to our current Conditions of Sale, copies of which are available on request.

© Micromass Limited, November 1998

Apparatus and methods for the analysis of trace



The Mass Spectrometry People