

Evaluation of the Stability of New Organic/Inorganic Hybrid Reversed-Phase HPLC Packings

B. Alden, J. Carmody, J. Grassi, C. Gendreau,
P. Iraneta, T. Walter


Waters Corporation, 34 Maple Street, Milford, MA 01757-3696

Abstract:

We have recently developed a new family of reversed-phase HPLC packings based on organic/inorganic hybrid particles. To evaluate the hydrolytic stability of these packings relative to state-of-the-art C₁₈-silicas, we have used several different test protocols. These protocols include chromatographic tests using mobile phases containing buffers that range from pH 1.2 to 11.5, as well as temperatures up to 60 °C. The results show that the packings based on hybrid particles exhibit exceptional stability in alkaline mobile phases and excellent stability in acidic mobile phases.


Hybrid Technology Particle

- Hybrid Organic/Inorganic materials contain both organic and inorganic components
- The hybrid particles described here were synthesized from $\text{Si}(\text{OEt})_4$ and $\text{CH}_3\text{Si}(\text{OEt})_3$:



- Hybrid particles combine:
 - efficiency and mechanical strength of silica
 - extended pH range and absence of base tailing of organic polymers

XTerra™ Bonded Phases



Bonded XTerra™ Particle

- XTerra™ bonded phases are based on the hybrid particles shown in the previous slide
- XTerra™ RP₁₈
 - Embedded carbamate C₁₈ phase with a trimethylsilyl endcap
 - 14.95% C
 - 2.32 $\mu\text{mol}/\text{m}^2$
- XTerra™ MS C₁₈
 - Trifunctional C₁₈ phase with a trimethylsilyl endcap
 - 15.45% C
 - 2.22 $\mu\text{mol}/\text{m}^2$

pH 10 Et₃N/ 50 °C Accelerated Test Protocol

Accelerated Test Procedure:


1. **Equilibrate:** 45 min, 1 mL/min, test mobile phase
2. **Test:** Mobile Phase: pH 7.0 K₂HPO₄ / MeOH 35:65 v/v, Flow Rate: 1 mL/min, Detection wavelength 254 nm
3. **Challenge:** 50 mM N(CH₂CH₃)₃·HCl buffer pH 10.0
Time: 60 min, Flow Rate: 2 mL/min
4. **Purge:** 10 min at 2 mL/min with water
10 min at 2 mL/min with methanol
5. **Cycle** back to step 1 and repeat until column fails

The column (4.6x150mm) is maintained at 50°C throughout the entire test.


- Under these aggressive conditions XTerra™ columns lasted significantly longer than state-of-the-art C₁₈-silica columns.

pH 10 Et₃N/ 50 °C Accelerated Test Chromatograms

Typical Chromatograms for C₁₈-Silica S





Typical Results - % of Initial Efficiency for Acenaphthene Reproducibility for 3 C₁₈-Silica S Columns



pH 10 Et₃N/ 50 °C Accelerated Test Efficiency vs. Time Curves

% Initial N (5 sigma) for Acenaphthene





pH 7 Na₂HPO₄/60 °C Accelerated Test Protocol

Accelerated Test Procedure:

1. **Equilibrate:** 60 min, 1.0 mL/min, test mobile phase 1
2. **Test Probe Mobile Phase1:** MeCN/H₂O (percentages vary with packing material), Flow rate: 1.0 mL/min, Detection wavelength: 254 nm, probes: uracil, acenaphthene
3. **Equilibrate:** 60 min, 1.5 mL/min, test mobile phase 2
4. **Test Probe Mobile Phase2:** 10 mM pH 7 Na₂HPO₄/MeCN (percentages vary with packing material), Flow rate: 1.5 mL/min, Detection wavelength: 254 nm, probes: tricyclic antidepressants
5. **Purge:** 20 min at 1 mL/min with water
10 min at 1 mL/min with acetonitrile
20 min at 1 mL/min with H₂O/MeCN, 90:10
6. **Column Aging Mobile Phase:** 50 mM pH 7 Na₂HPO₄/MeCN 80:20, Flow Rate: 1.0 mL/min
Time: 60, 120, 180, 240, 240... min
7. **Purge:** 20 min at 1 mL/min with water
10 min at 1 mL/min with acetonitrile
8. **Cycle** back to step 1 and repeat until column fails.


The column (4.6x150mm) is maintained at 60°C throughout the entire test.

pH 7 Na_2HPO_4 /60 °C Accelerated Test Chromatograms



pH 7 Na_2HPO_4 /60 °C Accelerated Test Efficiency vs. Time Curves

% of Initial Acenaphthene Efficiency (5 sigma)



Column Volumes of Phosphate Buffer


– Using aggressive phosphate buffers at high temperature, bonded hybrid packings resisted dissolution which leads to the premature failure of C_{18} -silica columns.

pH 11.5 Lifetime Test Protocol

pH 11.5 Lifetime HPLC Conditions:



- Column : XTerra™ RP₁₈, 5μm, 4.6x150mm with an XTerra™ RP₁₈, 5μm, 3.9x20mm guard column
- Mobile phase : 50 mM Pyrrolidine- HCl pH 11.5 / MeCN, 50:50
- Column Temperature : 30 °C
- Flow Rate : 1 mL/min.
- Detector : 215 nm
- Injection volume : 2 μL

pH 11.5 Lifetime Test Chromatograms




– Hybrid technology extended the high pH column lifetime of X Terra™ packing to 10 times longer than C18-Silica S.


pH 11.5 Lifetime Test Efficiency vs. Time Curves



1. Methamphetamine,
 pK_a 9.5




2. Propranolol,
 pK_a 9.5



3. Nortriptyline,
 pK_a 9.7

Low pH Stability of Hybrid Particle

Particle % C vs Time Exposed to 1 M HCl at 95°C



–The low pH stability of the methylsiloxane unit in the hybrid particle was tested using strongly acidic high temperature conditions. **No loss of the methylsiloxane units was detected.**


pH 1.5 TFA/50 °C Accelerated Test Protocol

Accelerated Test Procedure:

1. **Test:** MeOH/ 20 mM K₂HPO₄ pH 7.0, 65:35 v/v, 25 °C, Flow rate 1.0 mL/min
2. **Purge:** 1% trifluoroacetic acid (TFA) pH 1.2 , Flow rate: 2 mL/min
 - 12 minutes for 4.6x75mm columns
 - 24 minutes for 4.6x150mm columns
1. **Challenge:** 1% trifluoroacetic acid pH 1.2 at **50 °C (no flow)**
2. **Purge:**
 - water at 2.0 mL/min for 10 min
 - methanol at 2.0 mL/min for 10 min
3. **Test:** MeOH/ 20 mM K₂HPO₄ pH 7.0, 65:35 v/v, 25 °C, Flow rate 1.0 mL/min
4. Cycle back to step 1 and repeat until column fails.

pH 1.5 TFA/50 °C Accelerated Test Retention vs. Time Curves

% Initial Retention Factor (k) for Acenaphthene
vs Time Exposed to pH 1.2/50 °C



– Loss of retention is due to the hydrolytic loss of bonded phase.