# How to navigate this large-format poster:

- Choose the hand tool ; it becomes the article tool
  [ ] when your cursor rests anywhere on the poster.
- Click on the place where you want to start reading.
- To start at the beginning, press Ctrl [Windows] or Option [Mac OS] while you click anywhere.
- To scroll *forward* through the poster, continue **clicking**.
- To scroll backward, press Shift while clicking.
- To print, select Acrobat's *Fit to Page* option & choose largest available paper size.

# Click HERE to Read Poster on Next Page

# Waters

# Development and Evaluation of a High Peak Capacity Reversed-Phase HPLC Column for the Separation of Peptides under LC/MS Conditions

## Introduction

The separation and characterization of peptides derived from protein tryptic digest by reversed-phase HPLC coupled with mass spectrometry has become an important tool in the study of the proteome. This however, requires mass spectrometry compatible methodologies in order to eliminate ion suppression and a subsequent loss of sensitivity by the mass spectrometer. We have developed a high resolution, high peak capacity column that does not require TFA for the separation of peptides. A comparison of our Atlantis<sup>TM</sup> dC18 column to other commercially available columns under LC/MS compatible conditions was evaluated. The column was found to have increased peak capacity as compared to all others tested for peptide mapping.

## Experimental

Waters PDA 996

Waters Alliance<sup>™</sup> HPLC

#### Waters ZQ<sup>™</sup> Mass Detector

-Source = ESI (+), Capillary (kV) = 3.3, Cone (V) = 30 -Temperature (°C), Source = 150, Desolvation = 350-Gas flow (L/Hr), Cone = 50, Desolvation = 500-Scan Mode

#### Micromass Q-Tof II

-Source = ESI (+) -Capillary (kV) = 3.3-Cone(V) = 30-Collision Energy=30

Separation Method UV/MS:

-A: 0.1% Formic acid in water -B: 0.1% Formic acid in MeCN

#### or

-A: 0.020% TFA in water -B: 0.016% TFA in MeCN -Gradient from 0% B to 40% B in 45 min. -0.75 ml/min, split flow to ~0.2ml/min diverted to ZQ

#### Reversed-Phase Columns Tested

–Atlantis<sup>™</sup> dC18 4.6 x 50 mm 5µm -Wide Variety of Columns "Market" as Peptide Columns 4.6 x 50 mm 5µm

The Average USP Tailing Factor and Average Peak Capacity for a wide variety of Reversed Phased

#### Chemistries "Marketed" as Peptide Separation Columns by the Leading Manufaturers.

Columns tested using Formic Acid in a H<sub>2</sub>0 / ACN gradient (see methods) with a standard nine peptide mixture at 4.5µg mass

1044		
Column	Avg. USP TF	Avg. P (50%)
Atlantis™ dC18	1.23	309.2
Leading Competitor A (Chemistry 1)	1.22	300.2
Leading Competitor A (Chemistry 2)	1.27	273.7
Prototype	1.42	250.0
Leading Competitor E (Chemistry 1)	1.33	250.6
Leading Competitor C (Chemistry 1)	1.26	237.9
Leading Competitor A (Chemistry 3)	1.22	210.9
Leading Competitor D (Chemistry 1)	1.47	196.2
Leading Competitor B (Chemistry 1)	1.47	204.5
Leading Competitor B (Chemistry 2)	1.53	187.8
Leading competitor E (Chemistry 2)	1.79	152.7
Leading Competitor A (Chemistry 4)	1.81	127.9
Competitor F (Chemistry 1)	1.58	126.3
Leading Competitor A (Chemistry 5)	2.24	107.4

## Total Ion Chromatogram of Low pH Separation of Cytochrome C Tryptic Digest on Atlantis<sup>™</sup> dC<sub>18</sub>



Paul Rainville<sup>1</sup>, Jeff Mazzeo<sup>1</sup>, Ken Fountain<sup>2</sup>, Martin Gilar<sup>2</sup>, and Reb Russell II<sup>1\*</sup> <sup>1</sup>Waters Corporation, <sup>1</sup>Life Sciences Applications, <sup>2</sup>Life Sciences R&D Milford, MA, USA (\*send correspondence: reb\_russell@waters.com)



The Effects of Formic Acid compared with TFA on Peak Capacity between Reversed Phase Columns for the LC/MS analysis of peptides.

Mobile Phase A:H <sub>2</sub> 0 with 0.100% Formic Acid Mobile Phase B: ACN with 0.065% Formic Acid					
Peak #	Column A (300 Å)	Column A (100Å)	Atlantis™ dC18 (100Å)		
1	0.319	0.360	0.260		
2	0.436	0.296	0.300		
3	0.343	0.261	0.259		
4	0.304	0.244	0.234		
5	0.500	0.341	0.313		
6	0.732	0.580	0.492		
7	0.671	0.446	0.377		
8	0.499	0.319	0.269		
Peak Capacity*	95.9	126.4	143.8		

Mobile Phase A:H <sub>2</sub> 0 with 0.020% TFA Mobile Phase B: ACN with 0.016% TFA				
Peak #	Column A (300 Å)	Column A (100Å)	Polarity™ dC18 (100Å)	
1	0.323	0.28	0.326	
2	0.389	0.308	0.321	
3	0.283	0.91	0.267	
4	0.269	0.288	0.251	
5	0.400	0.409	0.362	
6	0.516	0.519	0.443	
7	0.414	0.439	0.377	
8	0.318	0.35	0.309	
Peak Capacity*	123.6	124.8	135.5	
%∆ in Peak Capacity* Formic/TFA	78%	101%	106%	

\*Peak Capacity =

17.50 22.50 25.00 27.50 30.00 32.50 15.00 20.00

#### Differences in Peak Shape, Peak Capacity, and Retention between Atlantis™ dC18 (100Å) and the Leading Peptide C18 Column (nominal 100Å) run in LC/MS conditions under different mass loads.

- A) Atlantis<sup>TM</sup> dC18 separating 4.5 μg of tryptic cytochrome C Peptides
- B) Atlantis<sup>™</sup> dC18 separating 45.0 µg of tryptic cytochrome C Peptides
- C) Leading Competitor B C<sub>18</sub> separating 4.5µg of tryptic cytochrome C Peptides
- D) Leading Competitor B  $C_{18}$  separating 45.0  $\mu$ g of tryptic cytochrome C Peptides



## Conclusions

- 100Å reversed-phase columns outperform 300Å columns for the vast majority of peptides obtained from enzymatic digestion of proteins under LC/MS separation conditions regardless of manufacturer.
- Atlantis<sup>TM</sup> dC18 and a very few other columns retain their superior performance characteristics under LC/MS conditions using Formic Acid.
- Atlantis<sup>TM</sup> dC18 retains it's superior performance characteristics under increasing mass loads while the majority of other columns quickly deteriorate
- Superior retention of Hydrophilic and all peptides