# Identification and Characterization of Oxidation and Deamidation of Proteins and Peptides by LC/MS(/MS)

# Waters

## Introduction

Characterization of proteins is an important analytical challenge that must be met in pharmaceutical and biotechnology companies today. Reversedphase chromatography coupled with mass spectrometry has rapidly become the technique of choice for the determination of post-translational modifications. We have developed a simple and rapid LC/MS method for identifying the presence and degree of oxidation and deamidation of intact proteins or peptides in a mixed sample using either single quadrupole or Q-TOF type mass spectral instrumentation.

## Experimental

### ntact Proteins

### Vaters Alliance™ LC and 996 PDA Detector

Vaters ZQ<sup>™</sup> Mass Detector or Micromass Q-TOF2<sup>™</sup>

- -Source = ESI (+): Capillary (kV) = 3.3, Cone (V) = 25
- -Temperature (°C): Source = 150, Desolvation = 350
- -Gas flow (L/Hr): Cone = 50, Desolvation = 500
- -Scan Mode

-Symmetry 300<sup>®</sup> C<sub>18</sub>, 4.6 X 50 mm, 5.0 μm –A: 0.1% Formic acid in water –B: 0.1% Formic acid in ACN -Gradient from 15% B to 80% B in 15 min -0.5 mL/min, Split Flow to ~0.2mL/min diverted to ZQ

### <u>Peptides</u>

### /aters ZQ™ Mass Detecto

- -Source = ESI (+): Capillary (kV) = 3.3, Cone (V) = 30-Temperature (°C): Source = 150, Desolvation = 350 -Gas flow (L/Hr): Cone = 50, Desolvation = 500
- -Scan Mode

### Aicromass Q-TOF2™

- -Source = ESI (+), Capillary (kV) = 3.3, Cone (V) = 30 Collision Energy=30 -Survey Scan MS 200-1600 m/z, MS/MS 50-1600 m/z
- -MaxEnt3 Processed Spectrum
- -BioLynx™ Peptide Sequencing

–Atlantis™ dC<sub>18.</sub> 4.6 X 50 mm, 3 μm -A: 0.1% Formic acid in water

- –B: 0.1% Formic acid in ACN
- -Gradient from 0% B to 40% B in 45 min
- -0.75 mL/min, split flow to ~0.2mL/min diverted to ZQ

### <u>Samples</u> Oxidatic

-Proteins and Peptides reacted for 4 hours using 3-30%  $H_2O_2$  (v/v)

-Heating sample at 65°C for 1 day to 1 month







Paul Rainville, Jeff Mazzeo, Uwe Neue, Reb Russell II\*

Waters Corporation, Life Sciences Applications, Milford, MA, USA \*Send Correspondence: reb\_russell@waters.com