

# **LC-MS** Analysis of Therapeutic and Diagnostic Oligonucleotides

Martin Gilar\*, Kenneth J. Fountain, John C. Gebler

Waters Corporation, Life Sciences R&D, Milford, MA, USA



#### Overview

- LC-MS method was developed for characterization of oligonucleotides.
- Method was applied for analysis of therapeutic and diagnostic oligonucleotides.
- Excellent LC separation was achieved
- Mobile phase is compatible with MS detection.

### Introduction

Synthetic oligonucleotides are utilized for diagnostic as well therapeutic purposes (antisense drugs). Synthesis of chemically modified oligonucleotides is challenging and often results in a product of limited purity [1,2]. An LC-MS method has been developed for the characterization of oligonucleotide-based drugs and diagnostic probes [3]. The method was utilized for analysis of phosphorothioate oligonucleotides, dually labeled probes and native oligonucleotides.

LC-MS is useful for nucleic acids analysis [4,5]. We used an optimized LC-MS system (see Methods) for analysis of synthetic oligonucleotides and their failure products.

Phosphorothioate as well as guanidine-rich antisense oligonucleotide drugs are known to be particularly difficult to analyze [6]. Due to the chaotropic properties of triethlylaminehexafluoroisopropanol (TEA-HFIP) buffers, we were able to successfully separate <60mer antisense oligonucleotides [3] and identify the failure oligonucleotides in a drug form. Guidelines for column selection and optimization of mobile phase composition are discussed. The ionpair reversed-phase HPLC separation benefits from the use of 2.5  $\mu$ m sorbent and elevated temperature. The separation performance was equivalent or better than ion-exchange HPLC and in some cases rivaled capillary gel electrophoresis separation. Using a 50 x 1 mm XTerra column, we were able to obtain molecular weight confirmation for ~ 1-10 pmole of oligonucleotides injected.

Figure 1: RP-HPLC offers a similar performance for the separation of oligonucleotides as capillary gel

Results



CGE: BioCAP<sup>™</sup> 75 µm, 27.5/34.5 cm, 30°C, polyethylene glycol sieving matrix, 15 kV run, 4s injection at 5kV, 15 kV run HPLC: XTerra<sup>®</sup> MS C18, 50 x 4.6 mm, 2.5 µm; 60°C, 0.5 ml/min, A: 10 % MeOH in 400 mM HFIP, 16.3 mM TEA, pH 7.92, B: 40 % MeOH in 400mM HFIP-16.3 TEA, pH 7.92, linear gradient from 16.7 % to 30 % B at 10 min, then from 30 % B to 47.15 % B at 30 min.

Figure 2: TEA-HFIP ion-pairing buffers allow for more sensitive oligonucleotide MS detection than

Figure 4: LC-UV-MS analysis of crude synthetic 25mer oligoG (oligodeoxyguanidine). Failure sequences are clearly separated and identified by MS.



XTerra<sup>®</sup> C18, 50 x 1 mm, 2.5 μm; A: 5% MeOH in 400mM HFIP: 16.3mM TEA, pH 7.9; B: 60% MeOH in 400mM HFIP: 16.3mM TEA, pH 7.9; gradient from 14.5% to 69.1% B in 120 minutes; 23.6 μl/min; 60 °C; UV 260 nm, 2.5 μl injection; 2.6 μg (318 pmole) total mass load.

**Figure 5**: LC-UV-MS analysis of 21mer dually-labeled diagnostic oligonucleotide 5' (FLSN)-CAC CTC CAG TGG AAA TCA AGT -(TAMRA) 3'. Oligonucleotide was synthesized in one step ("one pot" synthesis).



triethylammoniumn acetate mobile phases.



1 mg/ml solution of 25mer phosphorothioate oligonucleotide (7776.42 Da) was infused in suggested buffer in 50 % acetonitrile solution. The infusion flow rate was 10  $\mu$ l/min.

XTerra<sup>®</sup> MS C18, 50 x 1mm, 2.5  $\mu$ m; A: 8.6 mM TEA—100 mM HFIP, pH 8.25; B: methanol gradient starts from 15 % to 45 % MeOH in 30 minutes, 23.6  $\mu$ I/min, 50°C. Upper trace shows oligonucleotides labeled with TAMRA dye ( $\lambda$  max 556 nm), middle trace oligonucleotides labeled with fluorescein ( $\lambda$  max 500 nm), lower trace — total ion current. For more details see paper [3].

#### Methods

HPLC system:	CapLC <sup>®,</sup> Waters,
	equipped with a photodiode
	array detector
Column:	50 x 1 mm XTerra <sup>®</sup> MS
	C18, 2.5 μm
Column temperature:	50-60 °C
Mobile phase flow rate:	23.6 μl per minute.
HPLC conditions:	see figure captions

lon pairing buffers:

MS instrument:

MS conditions:

propanol (weak acid), TEA-HFIP buffer pH was 8.3. Alternatively, 16.3mM – 400 mM TEA-HFIP buffer, pH 7.9

> with a methanol gradient. ESI-TOF mass spectrometer, LCT, Micromass capillary 2000 V cone 18 V desolv. temperature 120 °C

**Figure 3:** LC-UV-MS analysis of 25mer phosphorothioate drug digestion. Metabolites are separated from parent drug and identified by MS.



XTerra<sup>®</sup> C18, 50 x 1 mm, 2.5 μm; A: 5% MeOH in 400mM HFIP: 16.3mM TEA, pH 7.9; B: 60% MeOH in 400mM HFIP: 16.3mM TEA, pH 7.9; gradient from 18.2% to 36.4% B in 40 minutes; 23.6 μl/min.; 60°C; UV 260 nm;

0.6  $\mu$ l injection represents 0.8  $\mu$ g of 25mer, which is 103 pmole; the 19mer was detected at the level of ~ 1.5 pmole.

## Conclusions

- Presented LC-MS method is simple, robust, and reproducible
- It allows for characterization of therapeutic and diagnostic oligonucleotides
- We successfully analyzed
  - phosphorothioate oligonucleotides
  - G-rich oligonucleotide sequences
  - singly and dually labeled diagnostic oligonucleotides
  - locked nucleic acids
- We achieved excellent LC separation with MS friendly mobile phases
- The LC-MS method is capable of detecting ~1-10 pmole of oligonucleotides

## References

- [1] W.J. Warren, G. Vella, Mol. Biotechnol. 4 (1995) 179.
- [2] M. Gilar, E.S.P. Bouvier, J. Chromatogr. A 890 (2000) 167.
- [3] M. Gilar, Anal. Biochem. 298 (2001) 196.
- [4] C.G. Huber, A. Krajete, Anal. Chem. 71 (1999) 3730.
- [5] A. Apffel, J.A. Chakel, S. Fischer, K. Lichtenwalter, W.S. Hancock, Anal.Chem. 69 (1997) 1320.
- [6] P.J. Oefner, J Chromatogr. B 739 (2000) 345.



spectra were deconvoluted

with MaxEnt1 software

Triethylamine ion-pairing

agent (8.6 mM) was buffered

with 100 mM hexafluoroiso-

was used. Oligonucleotides

were eluted from the column

