

Eric S. Grumbach, Diane M. Diehl, Bonnie Alden, Pamela Iraneta, Uwe Neue, Jeff Mazzeo
Waters Corporation, Chemistry Applied Technologies, Milford, MA, USA;
Presented at HPLC, Nice, France, 15th-19th June, 2003

OVERVIEW-

Hydrophilic interaction chromatography (HILIC) has been described as an alternative to traditional reversed-phase chromatography for the retention of very polar analytes. It is a *variation* of normal-phase chromatography without the disadvantages of using solvents that are not miscible in water.

In this mode, a high organic-low aqueous mobile phase is used with a polar stationary phase to achieve retention of very polar analytes. This highly organic (> 80%) mobile phase is ideal for enhancing sensitivity in mass spectrometry.

INTRODUCTION-

We have developed and optimized a stationary phase for HILIC that provides good retention of very polar basic analytes, including several active pharmaceutical ingredients and their metabolites and/or impurities. We have examined the effects of sample diluents on peak shape and sensitivity and will make recommendations based on this data. We have also compared the differences in LC/MS sensitivity between HILIC conditions (i.e. > 80% organic) on the HILIC column and 100% aqueous conditions on a reversed-phase column designed to retain polar analytes. Because of the high organic mobile phase starting conditions, plasma samples cleaned-up by protein precipitation or mixed-mode SPE can be directly analyzed by LC/MS without solvent evaporation and reconstitution.


Atlantis™
Columns



ATLANTIS™ HILIC SILICA –


- First HILIC column in a family of columns optimized for the retention of **highly** polar analytes.
- Silica column that has been optimized for HILIC.
- Atlantis™ HILIC Silica is ideal for the retention of very **polar bases**

RETENTION MECHANISM ON SILICA-




- Polar analyte **partitions** into and out of adsorbed water layer
- Charged polar analyte can undergo **cation exchange** with charged silanol groups
- **Combination** of these mechanisms results in **enhanced** polar retention

HILIC VS. REVERSED-PHASE RETENTION CHARACTERISTICS-




- Retention characteristics of cytosine were monitored under HILIC and reversed-phase conditions at pH 3, 4, and 5 with 10% incremental values of acetonitrile from 0 to 90%.
- The optimum operating range to achieve retention is above 70% organic for HILIC and below 20% organic for reversed-phase chromatography
- HILIC offers significantly more retention than reversed-phase chromatography for polar basic analytes in their respective optimum operating ranges for retention.

SOLVENT STRENGTH: HILIC-




In HILIC, utilizing a **less** polar solvent can **increase** the retention of polar analytes




Retention of later eluting peaks increases by substituting a portion of the mobile phase with a less polar solvent

COMPLEMENTARY SELECTIVITY TO REVERSED-PHASE CHROMATOGRAPHY-



Both techniques offer good retention and peak shape of the analytes. However, the **polar** metabolite elutes **after** the compound of interest on the Atlantis™ HILIC Silica column.

IMPORTANCE OF SAMPLE DILUENT-




Atlantis™ HILIC Silica
4.6 x 50 mm, 3 µm
10.0 µL Injection volume

Compounds Sample Conc.
1. 5-Fluorouracil 25 µg/mL
2. Uracil 25 µg/mL
3. 5-Fluorocytosine 25 µg/mL
4. Cytosine 25 µg/mL

Peak shape improves as % organic in the diluent increases.

Peak shape can be further improved by removing the aqueous portion of the diluent.



Atlantis™ HILIC Silica
4.6 x 50 mm, 3 µm
10.0 µL Injection volume


Compounds Sample Conc.
1. 5-Fluorouracil 25 µg/mL
2. Uracil 25 µg/mL
3. 5-Fluorocytosine 25 µg/mL
4. Cytosine 25 µg/mL

Peak shape improves as % acetonitrile in the diluent increases.

Peak shape is improved by increasing the non-polar portion of the diluent.

- In HILIC, it is important for the sample diluent to be 100% organic – remember, the organic mobile phase is the weak solvent.
- The sample diluent can significantly influence peak shape and peak area of your analytes of interest
- However, polar analytes often have low solubilities in organic solvents
- We ran an extensive series of experiments to determine a generic diluent
- **75:25 acetonitrile:methanol** is a useful generic diluent for most polar analytes
- This diluent is a compromise between solubility and peak shape

ENHANCED MS SENSITIVITY-



Atlantis™ dC₁₈
2.1 x 50 mm, 3 µm
1. Albuterol 100 pg/µL 21
2. Bamlanivir 50 pg/µL 14

Atlantis™ HILIC Silica Peak Area
2.1 x 50 mm, 3 µm
2. Bamlanivir 50 pg/µL 110085
1. Albuterol 100 pg/µL 19567


The retention of polar analytes by reversed-phase chromatography often requires non-volatile highly aqueous mobile phases which is not ideal for compound ionization by ESI-MS. Since HILIC requires a highly volatile high organic mobile phase, a dramatic increase in sensitivity can be observed.

The example above shows an increase in sensitivity by 3 to 4 orders of magnitude on the Atlantis™ HILIC Silica for these analytes. A lower limit of detection can be achievable because of the sensitivity that these highly volatile mobile phases promote.

LOW LIMITS OF DETECTION-

Atlantis™ HILIC Silica
2.1 x 50 mm, 3 µm
1. Acetylcholine
2. Choline

Concentration
100 fg/µL
100 fg/µL




Acetylcholine

Choline

Low limits of detection are possible by HILIC due to the higher sensitivity achieved from utilizing highly volatile mobile phases. A limit of detection of **100 fg/µL** of Acetylcholine and Choline was achieved using an Atlantis™ HILIC Silica column on a Waters® ZQ™ single quadrupole mass spectrometer

DIRECT INJECTION OF SPE ELUENT ONTO HILIC COLUMN-



- Evaporation and reconstitution are often the most lengthy steps of solid phase extraction. The elimination of these steps would prove useful by reducing sample preparation time and increasing throughput and productivity.
- The high organic SPE eluent is compatible with the high organic mobile phase conditions that HILIC requires, thus eliminating the need to evaporate and reconstitute the SPE eluent.
- If the SPE method is not optimized and there is residual matrix post SPE, the polar species from the matrix elute after the analytes of interest thus reducing potential suppression from matrix interferences.

CONCLUSIONS –

Atlantis™ HILIC Silica columns offer:

- Complementary selectivity to reversed-phase chromatography
- Retention of highly polar basic analytes not retained by reversed-phase chromatography
- Enhanced sensitivity in ESI-MS due to highly volatile mobile phases (> 80% organic)
- Lower limits of detection
- Shorter sample preparation procedures due to the elimination of the evaporation and reconstitution steps and directly injecting the eluent

WATERS CORPORATION
 34 Maple St.
 Milford, MA 01757 U.S.A.
 T: 508 478 2000
 F: 508 872 1990
www.waters.com

Waters, Micromass, Alliance, ZQ and Atlantis HILIC Silica are trademarks of Waters Corporation.
 All other trademarks are the property of their respective owners.
 ©2003 Waters Corporation

