

High-Throughput High-Sensitivity Oasis® μ Elution Technology

Ray Fisk, Jon Belanger, Pam Iraneta, Sue Serpa, Ziling Lu, Claude Mallet, Jeff Mazzeo, Uwe Neue and Yung-Fong Cheng, Waters Corporation, Milford, Massachusetts, USA.

Novel Oasis® μ Elution Design Enables Elution Volumes as Low as 25 μ L

The purpose of sample preparation is to selectively isolate analytes of interest from a complex sample matrix in a highly concentrated form prior to identification and quantification.

Protein precipitation (PPT) provides the advantage of simplicity and low cost. However, PPT results in a higher background and increased interference, which causes ion-suppression, reduced column lifetime and instrument downtime. The current trend in sample cleanup is to use the solid-phase extraction (SPE) technique. The SPE methods contain the steps of conditioning, loading, washing off of interferences, eluting analytes of interest, evaporating and reconstituting the eluate before the analysis (e.g., LC-MS/MS).

The evaporation and reconstitution steps not only take time and effort, but can also lead to loss of valuable sample. Therefore, the

Waters

ability to elute in very small volumes of solvent is desirable to minimize the amount of time required.

Currently, there are two technologies (Disc vs μ Elution) to minimize the elution volume for SPE (depicted in Figure 1). In the disc design (e.g., membrane/glass fibre disc), a small mass of sorbent is embedded into a thin structure with an aspect ratio (defined as the ratio of bed height, h , to bed diameter, D) of 0.13. The disc technology functions more like a filter and is designed to provide good flow properties and a low elution volume (in the range of 75–500 μ L). The hold-up volume is measured to be approximately 35–65 μ L.


In the μ Elution design, a much smaller sorbent mass (e.g., 2 mg) is packed into an internally tapered well with an aspect ratio of 1.15. The μ Elution technology functions more like a chromatography column and enhances the capture of target analytes. It also helps to prevent breakthrough during the load and wash steps. The hold-up volume for the μ Elution technology was determined to be as little as 15 μ L.

A small hold-up volume is important for achieving a low elution volume. Larger hold-up volumes result in larger elution volumes. This is illustrated schematically in Figure 1. After applying a vacuum, more solvent droplets are still held-up in the disc plate due to its large hold-up volume. However, this is not the case with the μ Elution plate since it has much smaller hold-up volume. Typical elution volumes are determined to be 75–300 μ L for the disc technology and only 25–100 μ L for the μ Elution plate.

Comparison of Protein Precipitation (PPT), Disc 96-Well SPE Plate, and Oasis® μ Elution 96-Well SPE Plate

In order to compare the different techniques, terfenadine-alcohol was added to raw uncentrifuged rat plasma at a concentration of 0.5 ng/mL. For the PPT method, sufficient acetonitrile (1 mL) was added to 250 μ L of sample solution to induce protein precipitation. This solution was centrifuged, and the supernatant was then evaporated to dryness. The sample was reconstituted with 25 μ L acetonitrile/isopropanol (40:60) with 2% ammonium hydroxide and finally further diluted with 50 μ L water. For SPE, 96-well

Figure 1: Comparison of disc and μ Elution Technology. Disc technology provides good flow properties with a lower capacity and a larger hold-up volume than the μ Elution Technology. Larger hold-up volumes require a larger elution volume because more liquid droplets remain in the hold-up volume.



The Oasis® μ Elution technology functions more like a chromatography column and enhances the capture of target analytes. It also helps to prevent breakthrough during the load and wash steps.

extraction plates were conditioned with 200 μ L MeOH, equilibrated with 200 μ L water. Then 250 μ L of sample solution was loaded into the wells, which were then washed with 5% methanol in water, and eluted with acetonitrile/isopropanol 40:60 with 2% formic acid (25 μ L for the Oasis® HLB μ Elution plate, 100 μ L for the disc plates). Finally, the eluate was diluted with water (50 μ L for μ Elution plate, 200 μ L for disc plates) through the plates, and 25 μ L was injected into the LC-MS/MS system.

The results are shown in Figure 2. Good peak shape and signal-to-noise (S/N) were observed from protein precipitation.

However, the disc technology provided an inferior S/N to the protein precipitation. Clearly, the best results were obtained using the Oasis® μ Elution technology.

Summary

The Oasis® μ Elution technology is designed to make an elution volume as low as 25 μ L possible, eliminating time-consuming and tedious evaporation steps involved in the SPE technique. Much better peak shape and signal-to-noise ratio could be easily achieved with the μ Elution plate compared with the disc technology or protein precipitation.

Waters S.A.,

BP 608, 78056 Saint-Quentin,
En Yvelines Cedex, France.

tel. +33 13048 7200, fax +33 13048 7201

website: www.waters.com.

Reader Service 259

Figure 2: Comparison of signal-to-noise ratio results from protein precipitation, Brand S and E disc 96-well extraction plates and Oasis® HLB μ Elution plate. Much better results were achieved with the Oasis® μ Elution technology than with protein precipitation or disc technology.

